

# Cabin Car Mk2



Pictured: Limited Edition Cabin Car Mk2<sup>3</sup>

<b>Class</b>	CLASS 3
<b>Weight Capacity</b>	178 kg / 28 stone
<b>Speed<sup>1</sup></b>	8 mph (UK & US Specification)*
<b>Ground Clearance<sup>1,2</sup></b>	75 mm (unloaded)
<b>Length<sup>2</sup></b>	1580 mm
<b>Width<sup>2</sup></b>	765 mm
<b>Height<sup>2</sup></b>	1690 mm
<b>Tyres</b>	13" / 330 mm
<b>Range<sup>1</sup></b>	75Ah: 22 miles** 100Ah : 30 miles**
<b>Motor</b>	Rear Wheel Drive Sealed transaxle 24 Volt DC Motor 1400w (4 Pole) motor
<b>Brakes</b>	Electromechanical
<b>Batteries</b>	Dual 12V Deep Cycle 75Ah & 100Ah Options
<b>Battery Charger</b>	8 amp off-board charger
<b>Max Slope<sup>2</sup></b>	12 degrees
<b>Seat Back Height<sup>2</sup></b>	737 mm / 29"
<b>Seat Width<sup>2</sup></b>	483 mm / 19"
<b>Turning Radius<sup>2</sup></b>	3250 mm

\* International Speeds:  
Available dependant  
on jurisdiction and local  
legislation

\*\*Range is measured on a  
smooth surface traveling at  
a constant speed. Real world  
mileage may be impacted by  
terrain and speed.

<sup>1</sup>Varies with user weight,  
terrain type, battery charge,  
battery condition, tyre  
condition.

<sup>2</sup>Due to manufacturing  
tolerances and continuous  
product improvement this  
specification can be subject to  
variance of +/-3%.

<sup>3</sup>Limited Edition Cabin Car  
pictured. Specification applies  
to all Mk2 variants unless  
specified.