







Reliable ultra-lightweight travel

The Ergo Lite is an ultra-lightweight wheelchair that weighs a mere 8.6 kg complete. The compact design is ideal for transportation and storage. It is equipped with Karma's unique S-Ergo seating system. The Ergo Lite is an ideal choice for those who need a simple, lightweight wheelchair.







S-Ergo Seating Design (Shaped to fit you)

The revolutionary S-Ergo system is the first s-shaped ergonomic seating in the world, designed according to the body's curve, which is capable of simultaneously relieving pressure, stabilizing seating position, and reduces the tendency to slide forward



P&P BRAKE (Dual Function Braking System) The dual function braking system allows carers to slow or stop the wheelchair.



Foldable Frame

Features

- * We believe that this is the lightest wheelchair of its type on the market.
- * S-Ergo system provides pressure relief and reduces sliding.
- * AEGIS, Microbe Shield is a permanently bonded antimicrobial barrier upholstery that protects you from odour, staining, and deterioration from bacteria, fungus, and other microorganisms.
- * Foldable backrest for easy transportation .
- * Fixed armrest and footrest hangers.
- * Folding frame.
- * Comfortable & durable upholstery.

Optional Accessories

- * Quick release wheels
- * Single wheel anti-tippers
- * Travel bag



Ergo Lite Specifications

Seat Width* Seat Depth (inch)	16/18* 16
Seat Width (mm)	400/450
Seat Depth (mm)	430
Seat Height (mm)	480
F/R wheels (inch)	F: 6 R: 14
Backrest Height (mm)	420
Armrest Height (mm)	210
Calf Length (mm)	395-410
Total Weight (kg)	8.6-8.8
Max. User Weight (ke	g) 100
Overall L x W x H (mm)	Min. 930*560*905 Max. 930*610*905
Overall Folded L x W x H (mm)	Min. 830*275*700 Max. 830*280*710

Remark : Karma reserves the right to modify information herein without further notice. " / " means optional : " - " means adjustable " means adjustable ±1cm ± 0.5kg 104-40615 Rev.02 Sep.2020



KARMA MOBILITY LIMITED TEL: 0845-630-3436 FAX: 0845-630-3736 Email : info@karmamobility.co.uk



Unit 6, Target Park, Shawbank Road, Redditch, B98 8YN U.K.