

Permobil

Progeo Joker R2

Individualized performance



Aluminum



Rigid



Adjustable

permobil

Individualized with room for adjustments

Progeo Joker R2 is a rigid wheelchair with a welded front frame that still offers seating and drive performance adjustability. The welded front end enhances the rigidity resulting in a more durable and agile ride. The option of a tapered frame provides optimal support for the legs with adequate clearance for the feet. This makes the Progeo Joker R2 an attractive solution for experienced wheelchair users.



Progeo Joker R2

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Front from 45 to 53 cm (adjustable) Rear from 38 to 46 cm (adjustable)
Footplate distance	From 35 to 49 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	Aluminium 90°
Total length	88 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 5,7 kg approx.
Wheelchair weight	Starting from 8,3 kg approx.
Maximum load	125 kg

OPAQUE
BLACK

PEARL
WHITE

OPAQUE
ORANGE

BRIGHT
RED

BRIGHT
SILVER

Personalise your ride*

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Actual product colors may vary.

permobil

permobil.com | customer_support_uk@permobil.com | T. +44 (0)1484 722 888



PROGEOJOKERR2DATASHEETUK211216