



eFOLDi[®]

eFOLDi Lite

The eFOLDi Lite is the worlds lightest, most compact, folding mobility scooter. Weighing only 15kg (without battery) and no larger than a small suitcase when folded. Amazingly easy to set-up and fold - fully airline compliant, the eFOLDi Lite is the perfect super lightweight travel companion.

Enjoy freedom, friendship and fun!



eFOLDi[®]

Features

Speed	4mph (6kmh)
Range	10miles / 16 Km
Carry capacity	120Kg / 265 lbs / 19 Stone
Climbing ability	1:10 Gradient (6 degrees)
Folded dimensions	65 x 55 x 37 cm / 26 x 22 x 15 inches
Brake	Electro magnetic
Wheels	Quick release rear
Ground clearance	7cm / 2.7'
Tyre type	Front solid / rear pneumatic
Battery	10AH Li-ion
Motor	180W 24V Hub Motor
Ignition	Electronic, push switch

Version: Lite 23062021

www.efoldi.com