Alumina scooter Range





Heaviest part just 10kg!

Alumina

- RANGE 15 miles
- BATTERY PACK 4.5kg

Alumina^{Pro}

RANGE 30 miles

BATTERY PACK 5kg

2 YEAR WARRANTY ON ALL MODELS



Come To Bits Easily

The Alumina range offers the innovation of both an Aluminium chassis and Lithium battery technology, making it the lightest car boot scooter with the longest range.

The Alumina offers up to 15 miles range with a lightweight 4kg battery pack. The Alumina Pro goes even further, with a 30 mile range (the longest of any car boot scooter) and only a 5kg battery pack. A smaller battery pack is also available which complies with the Civil Aviation Authority's regulations for air travel.

Both the Alumina and Alumina Pro is designed to be disassembled easily for transportation with a heaviest part weighing just 10kg (the lightest of any car boot scooter). Both are also fitted with all-round suspension, giving an impressive ground clearance. The long and flat floorpan and adjustable tiller also allow generous leg room and the soft grip tiller handlebar and a larger than standard seat make for an amazingly comfortable ride.



Full Suspension



Long Flat Pan For Leg Room



Soft Grip & Dual Controls



Large Comfortable Seat



High Ground Clearance & Large Wheels



Li-On Battery

CONTACT US www.motionhealthcare.co.uk 0844 257 8150

Alumina specifications



MODEL	Alumina	Alu mina ^{Pro}
LENGTH	90cm / 35"	90cm / 35"
WIDTH	50cm / 20"	50cm / 20"
HEIGHT	50cm / 20"	105cm / 41″
FOLDED LENGTH	N/A	N/A
FOLDED WIDTH	N/A	N/A
FOLDED HEIGHT	N/A	N/A
SEAT WIDTH	45cm / 18"	45cm / 18"
SEAT DEPTH	39cm / 15.5″	39cm / 15.5″
FRONT WHEEL	Solid 19cm / 7.5"	Solid 19cm / 7.5"
REAR WHEEL	Solid 19cm / 7.5"	Solid 19cm / 7.5"
TOTAL WEIGHT	73lbs / 33kg	75lbs / 34kg
HEAVIEST PART	22lbs / 9.98kg	22lbs / 9.98kg
BATTERY WEIGHT	9lbs / 4kg	11lbs / 5kg
MAX USER WEIGHT	18stone / 115kg	18stone / 115kg
RANGE	15miles / 24km	30miles / 48km
SPEED	4mph / 6kph	4mph / 6kph
BATTERY	24v 19.35ah Lithium	24v 37.2ah Lithium
MOTOR	24v 200w	24v 200w
GROUND CLEARANCE	5cm / 2″	5cm / 2″
MAX GRADIENT	6°	6°
TURNING RADIUS	125cm / 49″	125cm / 49″
BOX DEPTH	114cm / 45″	114cm / 45″
BOX WIDTH	53cm / 21″	53cm / 21″
BOX HEIGHT	53cm / 21″	53cm / 21″