

NEW
PRODUCT

THE COST EFFECTIVE
LITE SOLUTION.



Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST[®] —



Tel: 01869 324600

Email: sales@pride-mobility.co.uk

Web: www.pride-mobility.co.uk

(Version 1: March 2017)

Features and Specification



Introduction

Modern, lightweight and agile the Apex Lite from Pride® Scooters, is the latest addition to the popular Apex series.

The Apex Lite features 8" Alloy style wheels, detachable frame design and a comfortable folding seat aiding for a smoother more enjoyable journey. The Apex Lite is the perfect travel partner, helping you get to where you need to go.

Unique Features available with the Apex Lite



Easy to use handle bar controls



Chunky 8" alloy style wheels



Roomy backrest pocket



High-level tiller charger

Colours



Black














Blue



The Apex Lite can be dissembled into 5 pieces



Technical Specifications

Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
 21 Stone 6lbs (136 kg)	 4 Mph (6.4 km/h)	 3.4 cm ¹ (1.3 inch's)	 138.5 cm (54.4 inch's)	 102.2 cm (40.2 inch's)	 49 cm (19.2 inch's)	 Front: 20 cm Rear: 20 cm	 5.8 Stone 36.9 k.g.	 15 k.g. front section	 Folded Seat 10.33 k.g.	 x2 12 Ah deep cycle

1. Measured at Motor

Warranty

1. Life time warranty on frame
2. Two-year limited warranty on drive train, including: transaxle, motor and brake.
3. One-year limited warranty on all electrical parts
4. 12 month warranty for batteries

Notes

1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
3. Battery weight may vary based on manufacturer.
4. Information correct at time of printing
5. Images used for illustration purposes only.