

Invoice to:
 Name
 Street
 Town:
 E-Mail:
 Tel:

Deliver to:
 Name
 Street
 Town
 E-Mail:
 Tel:

Postal code:
 Fax:

Order date
 Purchase order:

Marked for:
 Sunrise order no:

Fax number: **+44 (0)845 605 66 89** Email: **orders@sunmed.co.uk**

Maximum user weight limit **125kg** Crash tested and passed

Measurements & adjustments
All measurements are in cm and <u>without cushion</u>
Seat depth 44cm as standard; adjustable from 41 to 46cm
Seat height 47.5cm (48cm with 12" wheels): adj. 42.5 to 50.5cm (1)
Back height 43.5cm as standard; adjustable to 41 and 46 cm

NEW STANDARD PACKAGE	
Colour	Brilliant Silver
Seat	Padded nylon
Backrest	13 degree half folding
Frame	Plus or minus 5cm centre of gravity adj
Back Uph.	7 Strap Tension adjustable
Armrest	Flip back - adjustable
Legrest	70 degree swing in/out
Footplates	Angle adjustable
Castor	8"x 1 1/4" solid
Rear Wheels	24" x 1" solid tyre / 12" transit
Brakes	Push to lock brake
Accessories	Positioning belt, tool kit, anti-tips (pair)

(1) Adjustment on Self propel version only

Base Price
£746



Please choose Self propel or Attendant Propel

BasiX 2 Self Propel - 24" Rear Wheels	
<input type="radio"/> 38cm seat width	074103-001
<input type="radio"/> 40.5cm seat width	074103-002
<input type="radio"/> 43cm seat width	074103-003
<input type="radio"/> 45.5cm seat width	074103-004
<input type="radio"/> 48cm seat width	074103-005
<input type="radio"/> 52cm seat width	074103-006

BasiX 2 Attendant Propel - 12" Rear Wheels	
<input type="radio"/> 38cm seat width	074104-001
<input type="radio"/> 40.5cm seat width	074104-002
<input type="radio"/> 43cm seat width	074104-003
<input type="radio"/> 45.5cm seat width	074104-004
<input type="radio"/> 48cm seat width	074104-005
<input type="radio"/> 52cm seat width	074104-006

SEAT UPHOLSTERY	£
<input checked="" type="checkbox"/> Standard padded nylon	
<input type="checkbox"/> Cushion 5cm	33
BACKREST	
<input checked="" type="checkbox"/> Tension adjustable back - 7 strap	
<input checked="" type="checkbox"/> Fixed push handles	
<input type="checkbox"/> Extension back (domed style headrest) <i>(Supplied as kit. Not factory fitted)</i>	132
ARMREST	
<input checked="" type="checkbox"/> Standard armrest with short & long arm pad <i>(tool height adjustable)</i>	
<input type="checkbox"/> Height adjustable armrest with short & long pad (sliding) <i>(Supplied as kit. Not factory fitted)</i>	83
LEGREST & FOOTPLATES	
<input checked="" type="checkbox"/> Angle adjustable footplate	
<input type="checkbox"/> Elevating legrest, pair <i>(Supplied as kit. Not factory fitted)</i>	277
<input type="checkbox"/> Amputee support pad, pair <i>(Supplied as kit. Not factory fitted)</i>	261
CASTOR	
<input checked="" type="checkbox"/> 8" x 1 1/4" solid	
WHEEL LOCKS	
<input checked="" type="checkbox"/> Push to lock brake	
<input type="checkbox"/> Lever extension <i>(Supplied as kit. Not factory fitted)</i>	59
ACCESSORIES	
<input type="checkbox"/> Crutch holder <i>(Supplied as kit. Not factory fitted)</i>	53
<input type="checkbox"/> Slide on tray <i>(Supplied as kit. Not factory fitted)</i>	132
<input type="checkbox"/> Step tube <i>(Supplied as kit. Not factory fitted)</i>	33

Note: Sunrise Medical strongly recommends the use of a step tube on all transit wheeled models, the 24" wheel model and also where attendant use is the predominant intended use.

Breezy BasiX 2 - Base Price	£746
Option price	£
Total price	£

**Customer Service Reference Only
Part Codes for options on BasiX Order Form**

BasiX²

IMPORTANT: CUSTOMER SERVICE MUST COMPLETE THIS INFORMATION:							
SEAT UPHOLSTERY - USE THE FOLLOWING IF REQUIRED							
Cushion 5cm	<input type="checkbox"/> 38cm seat width 064335-001	<input type="checkbox"/> 40.5cm seat width 064335-002	<input type="checkbox"/> 43cm seat width 064335-003	<input type="checkbox"/> 45.5cm seat width 064335-004	<input type="checkbox"/> 48cm seat width 064335-005	<input type="checkbox"/> 52cm seat width 064335-006	£33
BACKREST - USE THE FOLLOWING IF REQUIRED							
Extension back (domed style headrest)	<input type="checkbox"/> 40.5cm seat width 064920-001-ET	<input type="checkbox"/> 43cm seat width 064920-002-ET	<input type="checkbox"/> 45.5cm seat width 064920-003-ET	<input type="checkbox"/> 48cm seat width 064920-004-ET			£132
ARMREST - USE THE FOLLOWING IF REQUIRED							
Height adjustable armrest with short & long pad (sliding) (pair)			064475-100			£83	
LEGREST & FOOTPLATES - USE THE FOLLOWING IF REQUIRED							
Elevating legrest	<input type="checkbox"/> 38cm seat width	<input type="checkbox"/> 40.5cm seat width	<input type="checkbox"/> 43cm seat width	<input type="checkbox"/> 45.5cm seat width	<input type="checkbox"/> 48cm seat width	<input type="checkbox"/> 52cm seat width	£277
With footplate (Pair)	064897-101	064897-101	064897-102	064897-102	064897-103	064897-103	
Amputee support pad, pair	064280-000					£261	
WHEEL LOCKS - USE THE FOLLOWING IF REQUIRED							
Lever extension (Pair)			000861253			£59	
ACCESSORIES - USE THE FOLLOWING IF REQUIRED							
Slide on tray	<input type="checkbox"/> 38cm seat width 760848050.01	<input type="checkbox"/> 40.5cm seat width 760848050.01	<input type="checkbox"/> 43cm seat width 760848050.02	<input type="checkbox"/> 45.5cm seat width 760848050.03	<input type="checkbox"/> 48cm seat width 760848050.03	£132	
Crutch holder	067708					£53	
Step tube	064931					£33	