



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

JUNE 2023
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

Küschall K-Series 2.0

Standard Features

- Max User Weight 130kg (20.5 stone)
Crash Tested ISO 7176-Part 19
Active Rigid Wheelchair with Height Adjustable Rigid Footrest
Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
Full Hydroformed Aluminium Front Frame 75° Front Angle
24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Pneumatic Tyres
Choice of 10 Seat Widths 320mm to 500mm
Choice of 8 Seat Depths 350mm to 525mm
Range of Front Seat Heights from 450mm to 520mm
Range of Rear Seat Heights from 390mm to 490mm
Lockdown Folding, Angle and Height Adjustable Backrest
Black Nylon with Reflective Trim, Tension Adjustable Upholstery
Without Push Handles, Carbon Clothesguards
Black Aluminium Axle Tube with 3° Camber
Black Comfort Front Castor
Standard Brake Lever Bent
Tool Kit



Standard Base price: £2,540

Please make your selection with a cross in the box [] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with 2 columns: Code, Description. Includes BASE PRODUCT* and MYKÜSCHALL Special Request - Please include quotation reference in special description at the end of the form where applicable.

| FRAME* | | | | |
|--|--|---|---|------|
| DKA1040 | | V Style Standard Frame Angle 75° Aluminium | (NC 1475, 1490, 1495) | ◆ |
| DKA1050 | | Tapered Style Standard Frame Angle 75° Aluminium | (NC 1475, 1490, 1495, 1740, 1745, 1750, 3230) | ○ |
| DKA1060 | | V Style +5cm Frame Angle 75° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients | (NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760) | ○ |
| DKA1070 | | V Style Standard Frame Angle 90° Aluminium | (NC 1475, 1490, 1495) | ○ |
| DKA1080 | | Tapered Style Standard Frame Angle 90° Aluminium. Taper before the frame bend creating a parallel style front frame with a wider footrest space | (NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 5300, 5350, 5400, 5450, 5500) | ○ |
| DKA1090 | | V Style +5cm Frame Angle 90° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients | (NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760) | ○ |
| SUSPENSION | | | | |
| DKA1200 | | Seat Suspension, Weight Range 40-65 KG | (NC 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430) | £434 |
| DKA1210 | | Seat Suspension, Weight Range 65-85 KG | (NC 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430) | £434 |
| DKA1220 | | Seat Suspension, Weight Range 85-100 KG | (NC 1310, 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430) | £434 |
| SEAT WIDTH* | | | | |
| DKA1310 | | 320 mm | (NC 2050, 2610, 2630, 2710, 2720, 2748, 2750, 3420) | ○ |
| DKA1315 | | 340 mm | (NC 2710, 2720, 2748, 3420) | ○ |
| DKA1320 | | 360 mm | (NC 2710, 2748) | ○ |
| DKA1325 | | 380 mm | | ○ |
| DKA1330 | | 400 mm | | ◆ |
| DKA1335 | | 420 mm | | ○ |
| DKA1340 | | 440 mm | (NC 2630) | ○ |
| DKA1345 | | 460 mm | (NC 2630) | ○ |
| DKA1350 | | 480 mm | (NC 2630, 3230) | ○ |
| DKA1355 | | 500 mm | (NC 2630, 3210, 3230) | ○ |
| SEAT DEPTH* | | | | |
| DKA1410 | | 350 mm | (NC 1500, 1510, 2630, 2710, 2720, 2748, 7810, 7820, 7830) | ○ |
| DKA1420 | | 375 mm | (NC 2710, 2720, 2748, 7810, 7820, 7830) | ○ |
| DKA1435 | | 400 mm | | ◆ |
| DKA1450 | | 425 mm | | ○ |
| DKA1465 | | 450 mm | (NC 2630) | ○ |
| DKA1475 | | 475 mm | (NC 2630, 2748) | ○ |
| DKA1490 | | 500 mm | (NC 2630, 2710, 2748) | ○ |
| DKA1495 | | 525 mm | (NC 2630, 2610, 2710, 2748) | ○ |
| SEAT POSITION. Seat module can be moved forward and backward along the frame to achieve a more active or stable rear wheel position* | | | | |
| DKA1500 | | Position 1 = 140 mm MOST ACTIVE POSITION | | ○ |
| DKA1510 | | Position 2 = 125 mm | | ○ |
| DKA1520 | | Position 3 = 110 mm | | ◆ |
| DKA1530 | | Position 4 = 95 mm | | ○ |
| DKA1540 | | Position 5 = 80 mm LEAST ACTIVE POSITION | | ○ |

| BACK HEIGHT. Transportation kit is only possible with backrest heights of 405mm and above* | | | |
|---|--|---|---|
| DKA1600 | | Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if Matrx/Fixed backrest is required. | (NC 2100, 2220, 2230) <input type="radio"/> |
| DKA1605 | | 270 mm | (NC 2150, 2240, 2430, 7600) <input type="radio"/> |
| DKA1610 | | 285 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1615 | | 300 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1620 | | 315 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1625 | | 330 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1630 | | 345 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1635 | | 360 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1640 | | 375 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1645 | | 390 mm | (NC 2150, 2240, 7600) <input type="radio"/> |
| DKA1650 | | 405 mm | (NC 2150, 2240) <input type="radio"/> |
| DKA1655 | | 420 mm | (NC 2150, 2240) <input type="radio"/> |
| DKA1660 | | 435 mm | (NC 2150, 2240) <input type="radio"/> |
| DKA1665 | | 450 mm | (NC 2150, 2240) <input type="radio"/> |
| DKA1670 | | 465 mm | (NC 2150, 2240) <input type="radio"/> |
| DKA1675 | | 480 mm | (NC 2150, 2240) <input type="radio"/> |
| FRONT SEAT HEIGHT. For 75° frames, the maximum lower leg must be front seat height minus 30mm and minus 40mm for 90° frames* | | | |
| DKA1740 | | 450 mm | (NC 5300, 5350, 5400, 5450, 5500) <input type="radio"/> |
| DKA1745 | | 460 mm | (NC 5300, 5350, 5400, 5450, 5500) <input type="radio"/> |
| DKA1750 | | 470 mm | <input type="radio"/> |
| DKA1755 | | 480 mm | <input type="radio"/> |
| DKA1760 | | 490 mm | <input type="radio"/> |
| DKA1765 | | 500 mm | <input checked="" type="radio"/> |
| DKA1770 | | 510 mm | <input type="radio"/> |
| DKA1775 | | 520 mm | <input type="radio"/> |
| REAR SEAT HEIGHT* | | | |
| DKA1815 | | 390 mm | (NC 3490, 3495, 6030, 6040, 6430) <input type="radio"/> |
| DKA1820 | | 400 mm | (NC 3490, 3495, 6030, 6040, 6430) <input type="radio"/> |
| DKA1825 | | 410 mm | (NC 3490, 3495) <input type="radio"/> |
| DKA1830 | | 420 mm | <input type="radio"/> |
| DKA1835 | | 430 mm | <input type="radio"/> |
| DKA1840 | | 440 mm | <input checked="" type="radio"/> |
| DKA1845 | | 450 mm | <input type="radio"/> |
| DKA1850 | | 460 mm | <input type="radio"/> |
| DKA1855 | | 470 mm | <input type="radio"/> |
| DKA1860 | | 480 mm | <input type="radio"/> |
| DKA1865 | | 490 mm | (NC 6410) <input type="radio"/> |

| LOWER LEG LENGTH. Cushion height must be deducted* | | | |
|---|--|--|----------------------------|
| DKA1918 | | 290 mm | ○ |
| DKA1920 | | 300 mm | ○ |
| DKA1922 | | 310 mm | ○ |
| DKA1924 | | 320 mm | ○ |
| DKA1926 | | 330 mm | ○ |
| DKA1928 | | 340 mm | ○ |
| DKA1930 | | 350 mm | ○ |
| DKA1932 | | 360 mm | ○ |
| DKA1934 | | 370 mm | ○ |
| DKA1936 | | 380 mm | ○ |
| DKA1938 | | 390 mm | ○ |
| DKA1940 | | 400 mm | ○ |
| DKA1942 | | 410 mm | ○ |
| DKA1944 | | 420 mm | ○ |
| DKA1946 | | 430 mm | ◆ |
| DKA1948 | | 440 mm | ○ |
| DKA1950 | | 450 mm | ○ |
| DKA1952 | | 460 mm | ○ |
| DKA1954 | | 470 mm | ○ |
| DKA1956 | | 480 mm | ○ |
| BACKREST MATERIAL* | | | |
| DKA2000 | | Aluminium | ◆ |
| DKA2050 | | Titanium | (NC 2150, 2430, 7600) £206 |
| BACKREST TYPE* | | | |
| DKA2100 | | Standard Backrest with Backbrace, lockdown folding for easy transportation | (NC 2240, 2440) ◆ |
| DKA2150 | | Backrest with Lower Backbrace for use with Matrix fixed back. Backrest Tube Height is 285 mm, lockdown folding for easy transportation | (NC 2220, 2230, 2430) ○ |
| BACKREST UPHOLSTERY* | | | |
| DKA2220 | | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | (NC 2440) ◆ |
| DKA2230 | | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener | (NC 2440) £223 |
| DKA2240 | | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required | ○ |
| PUSH HANDLES* | | | |
| DKA2400 | | Without Push Handles | ◆ |
| DKA2410 | | Push Handles Fixed Short/Mini Hand Grip | £133 |
| DKA2420 | | Push Handles Fixed Long Hand Grip | £133 |
| DKA2430 | | Push Handles Height Adjustable, Rear Set. Maximum height of push handles is rear seat height + 600mm | (NC 3480) £254 |
| DKA2440 | | Push Handles Height Adjustable, Integrated in backrest | £219 |
| DKA2450 | | Push Handles Foldable | £198 |

| BACKREST ANGLE* | | | |
|--|--|---|--|
| DKA2520 | | 78° angle from the seat | ○ |
| DKA2540 | | 82° angle from the seat | ○ |
| DKA2560 | | 86° angle from the seat | ○ |
| DKA2580 | | 90° angle from the seat | ◆ |
| DKA2590 | | 94° angle from the seat | ○ |
| SEAT UPHOLSTERY* | | | |
| DKA2600 | | Standard Nylon Seat Upholstery with one 5cm central band (Loop) | ◆ |
| DKA2610 | | Tension Adjustable Seat Upholstery with one 5cm central band (Loop) | £45 |
| DKA2630 | | Light Seat Upholstery, mesh material to reduce weight | £69 |
| SEAT CUSHION | | | |
| DKA2700 | | Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk | (NC 2810, 2820, 2830, 2840, 2850) £66 |
| DKA2710 | | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk | (NC 2810, 2820, 2830, 2840, 2850) £121 |
| DKA2720 | | Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk | £410 |
| DKA2748 | | Contour Visco NG Cushion 100mm High - High Risk | (NC 2810, 2820, 2830, 2840, 2850) £169 |
| DKA2750 | | Vicair Active O2 lightweight cushion 60mm High - High risk | (NC 2810, 2820, 2830, 2840, 2850) £521 |
| SEAT BAG | | | |
| DKA2800 | | Seat Bag, Click in and removable | £74 |
| MATRX LIBRA CUSHION ACCESSORIES | | | |
| DKA2810 | | Fluid supplement sac (Pair) | £44 |
| DKA2820 | | Tapered Wedge Pair | £30 |
| DKA2830 | | Abductor | £23 |
| DKA2840 | | Lateral Side Wedge | £20 |
| DKA2850 | | Anterior Front Wedge | £27 |
| FOOTREST* | | | |
| DKA3200 | | Footrest Tube, Rigid - Aluminium Hydroformed Profile | ◆ |
| DKA3210 | | Footrest Tube, Rigid - Titanium (Rounded tube) | £185 |
| DKA3220 | | Carbon Footplate, Angle Adjustable | (NC 3290) £225 |
| DKA3230 | | Footplate Flip-to Side, Angle Adjustable | (NC 3290) £143 |
| FOOTREST COVER | | | |
| DKA3290 | | Footrest Cover | £37 |
| CALF STRAP* | | | |
| DKA3300 | | Calf strap mounted on the frame with hook and loop | ◆ |
| DKA3305 | | Calf strap mounted with bracket support | ○ |
| CLOTHESGUARDS/MUDGUARDS* | | | |
| DKA3420 | | Carbon Clothesguards, Flip up, Fold in | ◆ |
| DKA3440 | | Carbon Mudguards, Removable | (NC 3490, 3495) £307 |
| ARMREST | | | |
| DKA3480 | | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included | £241 |
| DKA3490 | | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included | £164 |
| DKA3495 | | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included | £164 |

| FRAME COLOUR* | | | |
|--------------------|--|---|------|
| DKA4000 | | Space Black (RAL 9005) | ○ |
| DKA4010 | | Anthracite Grey (RAL 7016) | ○ |
| DKA4020 | | Pure White (RAL 9016) (NC 4220) | ○ |
| DKA4030 | | Carmin Red (RAL 3002) | ○ |
| DKA4040 | | Azure Blue (RAL 5017) | ○ |
| DKA4060 | | Flashy Yellow (RAL 1016) (Additional 10 day lead time) | £31 |
| DKA4070 | | Atomic Violet (Additional 10 day lead time) | £31 |
| DKA4100 | | Black Glossy (RAL 9005) (Additional 10 day lead time) | £31 |
| DKA4130 | | Silver Grey (RAL 9006) (Additional 10 day lead time) | £31 |
| DKA4140 | | Royal Blue (Additional 10 day lead time) | £31 |
| DKA4150 | | Ruby Red (RAL 3003) (Additional 10 day lead time) | £31 |
| DKA4170 | | Dynamite Orange (Additional 10 day lead time) | £31 |
| DKA4180 | | Forest Green (RAL 6029) (Additional 10 day lead time) | £31 |
| DKA4050 | | Personal frame colour. Provide standard RAL colour code of your choice | £311 |
| FRAME DECAL TYPE* | | | |
| DKA4190 | | Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111) | £31 |
| DKA4195 | | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) | ◆ |
| DKA4197 | | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) | ○ |
| BOLD DECAL COLOUR* | | | |
| DKA4200 | | Black | ○ |
| DKA4210 | | Dark Grey | ○ |
| DKA4220 | | Light Grey | ○ |
| DKA4230 | | Red | ○ |
| DKA4240 | | Blue | ○ |
| DKA4250 | | Gold | ○ |
| DKA1111 | | No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag (NC 4190) | ○ |
| TAG COLOUR* | | | |
| DKA4300 | | Black | ○ |
| DKA4310 | | Red | ○ |
| DKA2222 | | No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag (NC 4195) | ○ |

| FRONT CASTORS* | | | |
|--|--|---|------|
| DKA5000 | | 75mm x 25mm (3" x 1") Supersport | ○ |
| DKA5050 | | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre | ○ |
| DKA5100 | | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre | ○ |
| DKA5150 | | 100mm x 25mm (4" x 1") Starec Black hub and tyre | £83 |
| DKA5200 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) | £93 |
| DKA5250 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £93 |
| DKA5300 | | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre | ○ |
| DKA5350 | | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre | ○ |
| DKA5400 | | 125mm x 25mm (5" x 1") Starec Black hub and tyre | £83 |
| DKA5450 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) | £93 |
| DKA5500 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | £93 |
| REAR WHEEL POSITION* | | | |
| DKA6020 | | Standard, Fixed Rear Wheel Position (NC 6100, 6110, 6120, 6130, 6140, 6230, 6240) | ◆ |
| DKA6030 | | Light Rear Wheel Extension (3 Positions) (NC 6070, 6130, 6140, 6230, 6240) | £113 |
| DKA6040 | | Rear Wheel Extension (5 Positions) (NC 6070) | £131 |
| CAMBER AXLE TYPE* | | | |
| DKA6050 | | Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°. Only active/scissor brake will lock 10° camber position (NC 6200, 6210, 6220, 6270, 7230) | £287 |
| DKA6060 | | Standard Aluminium Axle (NC 4444) | ◆ |
| DKA6070 | | Carbon Axle (NC 4444) | £212 |
| POSITION OF AXLE | | | |
| DKA6100 | | 1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position | ○ |
| DKA6110 | | 2nd Position from the front (NC 6230, 6240) | ○ |
| DKA6120 | | 3rd Position from the front (NC 6230, 6240) | ○ |
| DKA6130 | | 4th Position from the front (NC 6230, 6240) | ○ |
| DKA6140 | | 5th Position from the front (NC 6230, 6240) | ○ |
| CAMBER ANGLE* | | | |
| DKA6200 | | 0° Camber (NC 6270) | ○ |
| DKA6210 | | 3° Camber | ◆ |
| DKA6220 | | 6° Camber (NC 6240, 6250) | ○ |
| DKA4444 | | Without camber angle, select only with DKA6050 Vario-Axle | ○ |
| 2ND CAMBER AXLE TYPE. If a 2nd axle is selected, 1st axle must be placed in 1st position and the 2nd axle in the 5th position | | | |
| DKA6230 | | Standard Aluminium Axle | £232 |
| DKA6240 | | Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10° (NC 6250, 6260, 6270, 7230) | £287 |
| CAMBER ANGLE 2ND AXLE. Max of 3° difference in camber from 1st axle. If only 1 of the axles is a Vario-Axle, 2nd axle is only possible in 0° or 3° camber | | | |
| DKA6250 | | 0° Camber | ○ |
| DKA6260 | | 3° Camber | ○ |
| DKA6270 | | 6° Camber | ○ |
| DKA5555 | | Without camber angle, select only with DKA6240 Vario-Axle | ○ |
| QUAD QUICK RELEASE AXLE | | | |
| DKA6300 | | Quad Quick Release/Tetra Clip (NC 6550) | £85 |

| REAR WHEEL SIZE* | | | |
|------------------|--|---|---|
| DKA6410 | | 24" | ◆ |
| DKA6420 | | 25" | (NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900, 6910) ○ |
| DKA6430 | | 26" | (NC 6520, 6550, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770, 7240) ○ |
| REAR WHEEL TYPE* | | | |
| DKA6510 | | 24", 26" Kuschall Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim | ◆ |
| DKA6520 | | 24", 25" High Performance Wheel Black hub and rim | (NC 6660) £182 |
| DKA6550 | | 24", 25" Spinerger LXX, Black hub, rim and 18 spokes in Black | (NC 6660, 6910) £611 |
| SWO0907 | | 24", 25" Spinerger SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKuschall | £486 |
| SWO0909 | | 24", 25" Spinerger LX, Chrome hub and 12 spokes in Black. Coloured spokes are available via MyKuschall | £735 |
| REAR WHEEL TYRE* | | | |
| DKA6600 | | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection | ◆ |
| DKA6625 | | 24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant. | £100 |
| DKA6630 | | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) | £97 |
| DKA6660 | | 24" x 1 3/8" Solid Profile Black Tyre | (NC 7230) £66 |
| HANDRIM* | | | |
| DKA6700 | | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance | £223 |
| DKA6710 | | 24", 26" Steel, short tab with 1 hole, mounted at 3cm distance | £62 |
| DKA6720 | | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance | £192 |
| DKA6740 | | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip | £320 |
| DKA6750 | | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £320 |
| DKA6760 | | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion | £284 |
| DKA6770 | | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim | £356 |
| DKA6780 | | 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance | ◆ |
| DKA6790 | | 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting | £42 |
| SPOKEGUARD | | | |
| DKA6900 | | 24", 26" Transparent Spokeguard | £126 |
| DKA6910 | | 24", 26" Black Spokeguard | £126 |
| BRAKES* | | | |
| DKA7210 | | Standard Push to Lock Brakes, Lever Bent | ◆ |
| DKA7230 | | Performance Brakes, Push to Lock. Brake barrel flips away from tyre | (NC 7300) £78 |
| DKA7240 | | Active/Scissor Brakes | (NC 7300) £78 |
| DKA7250 | | Light/Compact Active, Aluminium Brakes | (NC 7300) £136 |

| BRAKE LEVER EXTENSION | | | |
|------------------------------|--|--|------|
| DKA7300 | | Brake Lever Extension | £49 |
| SAFETY OPTIONS | | | |
| DKA7500 | | Anti-Tipper Right (NC 7520) | £122 |
| DKA7510 | | Anti-Tipper Left (NC 7520) | £122 |
| DKA7520 | | Anti-Tipper Both sides | £244 |
| DKA7600 | | Transportation Kit Includes Standard Pelvic Belt and Tiedown Labels (NC 7800) | £74 |
| DKA7800 | | Standard Pelvic belt | £74 |
| BODYPOINT BELTS | | | |
| DKA7810 | | Bodypoint Posture Belt Evoflex Small | £81 |
| DKA7820 | | Bodypoint Posture Belt Evoflex Medium | £81 |
| DKA7830 | | Bodypoint Posture Belt Evoflex Large | £81 |
| OTHER OPTIONS | | | |
| DKA8010 | | Front Frame Protectors; 2 Piece. Removable with hook and loop | £37 |
| DKA8200 | | Cane Holder, Mounted Right Side | £75 |
| DKA8210 | | Cane Holder, Mounted Left Side | £75 |
| DKA8500 | | Air Pump | £13 |
| DKA8550 | | Tool Kit | ◆ |
| DKA8680 | | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | £320 |
| BODYPOINT ACCESSORIES | | | |
| DKA8800 | | Bodypoint Calf Strap | £57 |
| DKA8850 | | Bodypoint Mobility Bag | £55 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| DKA1310 | 320 mm |
| DKA1350 | 480 mm |
| DKA1355 | 500 mm |
| DKA1410 | 350 mm |
| DKA1420 | 375 mm |
| DKA1435 | 400 mm |
| DKA1475 | 475 mm |
| DKA1490 | 500 mm |
| DKA1495 | 525 mm |
| DKA1500 | Position 1 = 140 mm MOST ACTIVE POSITION |
| DKA1510 | Position 2 = 125 mm |
| DKA1520 | Position 3 = 110 mm |
| DKA1740 | 450 mm |
| DKA1745 | 460 mm |
| DKA1750 | 470 mm |
| DKA1755 | 480 mm |
| DKA1760 | 490 mm |
| DKA1815 | 390 mm |
| DKA1820 | 400 mm |
| DKA1825 | 410 mm |
| DKA1830 | 420 mm |
| DKA1835 | 430 mm |
| DKA2050 | Titanium |
| DKA2100 | Standard Backrest with Backbrace, lockdown folding for easy transportation |
| DKA2150 | Backrest with Lower Backbrace for use with Matrix fixed back. Backrest Tube Height is 285 mm, lockdown folding for easy transportation |
| DKA2220 | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener |
| DKA2230 | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener |
| DKA2240 | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required |
| DKA2430 | Push Handles Height Adjustable, Rear Set. Maximum height of push handles is rear seat height + 600mm |
| DKA2440 | Push Handles Height Adjustable, Integrated in backrest |
| DKA2610 | Tension Adjustable Seat Upholstery with one 5cm central band (Loop) |
| DKA2630 | Light Seat Upholstery, mesh material to reduce weight |
| DKA2710 | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk |
| DKA2720 | Matrix Libra Cushion, Black Smooth Cover 70mm High - Very High Risk |
| DKA2748 | Contour Visco NG Cushion 100mm High - High Risk |
| DKA2750 | Vicair Active O2 lightweight cushion 60mm High - High risk |
| DKA2810 | Fluid supplement sac (Pair) |
| DKA2820 | Tapered Wedge Pair |
| DKA2830 | Abductor |
| DKA2840 | Lateral Side Wedge |
| DKA2850 | Anterior Front Wedge |
| DKA3210 | Footrest Tube, Rigid - Titanium (Rounded tube) |

| | |
|---------|---|
| DKA3230 | Footplate Flip-to Side, Angle Adjustable |
| DKA3290 | Footrest Cover |
| DKA3420 | Carbon Clothesguards, Flip up, Fold in |
| DKA3480 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included |
| DKA3490 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included |
| DKA3495 | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included |
| DKA4190 | Bold Decal (Full coloured transfer on frame) |
| DKA4195 | Tag (Small KÜSCHALL brand transfer on front of frame) |
| DKA4200 | Black |
| DKA4210 | Dark Grey |
| DKA4220 | Light Grey |
| DKA4230 | Red |
| DKA4240 | Blue |
| DKA4250 | Gold |
| DKA1111 | No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag |
| DKA4300 | Black |
| DKA4310 | Red |
| DKA2222 | No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag |
| DKA5300 | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre |
| DKA5350 | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre |
| DKA5400 | 125mm x 25mm (5" x 1") Starec Black hub and tyre |
| DKA5450 | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) |
| DKA5500 | 125mm x 36mm (5" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) |
| DKA6030 | Light Rear Wheel Extension (3 Positions) |
| DKA6040 | Rear Wheel Extension (5 Positions) |
| DKA6070 | Carbon Axle |
| DKA6100 | 1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position |
| DKA6110 | 2nd Position from the front |
| DKA6120 | 3rd Position from the front |
| DKA6130 | 4th Position from the front |
| DKA6140 | 5th Position from the front |
| DKA6200 | 0° Camber |
| DKA6210 | 3° Camber |
| DKA6220 | 6° Camber |
| DKA4444 | Without camber angle, select only with DKA6050 Vario-Axle |
| DKA6230 | Standard Aluminium Axle |
| DKA6240 | Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10° |
| DKA6250 | 0° Camber |
| DKA6260 | 3° Camber |
| DKA6270 | 6° Camber |
| DKA6410 | 24" |
| DKA6420 | 25" |
| DKA6430 | 26" |
| DKA6510 | 24", 26" KÜSCHALL Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim |

| | |
|---------|---|
| DKA6520 | 24", 25" High Performance Wheel Black hub and rim |
| DKA6550 | 24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black |
| DKA6600 | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection |
| DKA6630 | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) |
| DKA6660 | 24" x 1 3/8" Solid Profile Black Tyre |
| DKA6700 | 24", 26" Supergrripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance |
| DKA6710 | 24", 26" Steel, short tab with 1 hole, mounted at 3cm distance |
| DKA6720 | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance |
| DKA6740 | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip |
| DKA6750 | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands |
| DKA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion |
| DKA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim |
| DKA6900 | 24", 26" Transparent Spokeguard |
| DKA6910 | 24", 26" Black Spokeguard |
| DKA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre |
| DKA7240 | Active/Scissor Brakes |
| DKA7300 | Brake Lever Extension |
| DKA7520 | Anti-Tipper Both sides |
| DKA7600 | Transportation Kit Includes Standard Pelvic Belt and Tiedown Labels |
| DKA7800 | Standard Pelvic belt |
| DKA7810 | Bodypoint Posture Belt Evoflex Small |
| DKA7820 | Bodypoint Posture Belt Evoflex Medium |
| DKA7830 | Bodypoint Posture Belt Evoflex Large |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):