RHealthcare

It's ability that matters



# Dash E-Fold Powerchair

The compact, power folding powerchair. Stylish with a powerful yet smooth driving experience.





### POWER FOLDING

Electric foldable frame by remote key with quick release and power-folding actuator.



#### **POWER DRIVE**

Powerful 2x250w motor and a removable lithium battery with a range of 9 miles. Add an extra 10 amp lithium battery to double your distance (optional).



#### SPACE SAVING

Can stand upright or flat when folded.

# Dimensions

| Total width     | 630mm    |
|-----------------|----------|
| Total length    | 1000mm   |
| Total height    | 900mm    |
| Seat width:     | 450mm    |
| Seat depth      | 480mm    |
| Backrest height | 520mm    |
| Backrest angle  | 8° - 30° |
| Armrest height  | 240mm    |
| Footrest length | 260mm    |
| Front wheels    | 8″       |
| Rear wheels     | 12.5″    |

## Other technical data

| Motor  | 2x250w                |
|--|-----------------------|
| Battery capacity                                 | 24v 10Ah Lithium      |
| Charger  | DC24v/2.0A            |
| Speed  | 4 mph                 |
| Travel distance<br>With extra battery (optional) | ≥9 miles<br>≥18 miles |
| Max. Climbable obstacle height                   | 100mm                 |
| Max. safe slope                                  | 10° - ±1°             |
| Turning diameter                                 | 900mm                 |



**Weight info** Chair weight 25.5kg Max user weight 125kg

# **Key features**

- ⊘ Foldable frame
- $\oslash$  Quick release for folding
- ⊘ Power folding actuators
- ⊘ Adjustable backrest in three positions
- ⊘ Detachable lithium battery



For enquiries:

**T:** 0845 146 0600 **E:** sales@rhealthcare.co.uk

**F:** 01246 453881

www.rhealthcare.co.uk



