




The Dual Wheel Auto Folding Scooter can be easily folded and unfolded by remote control in under 20 seconds, providing a quick and easy transport solution for users on the go.

-  **Copper**
HW009COP
-  **Purple**
HW009PUR
-  **Red**
HW009RD

Features

- ▶ Maximum range of 6 miles on full battery charge*
- ▶ Folds to a compact size within 20 seconds
- ▶ Compact scooter designed to fit into small car boots
- ▶ Lightweight folding scooter weighs in at only 21.5kg, making it ideal for cruises, flights and bus trips
- ▶ Scooter can be folded and unfolded using the remote control or buttons on the rear shroud of the scooter itself
- ▶ Comfortable and padded compact seat with auto-folding back
- ▶ Automatic electromagnetic brake system
- ▶ Height adjustable tiller to suit user requirements
- ▶ Solid puncture proof tyres
- ▶ Additional safety features include automatic speed reduction when cornering to aid stability and an auto-stop function to prevent the scooter from folding whilst in use
- ▶ On board charging**
- ▶ Anti-slip carpet in footwell for added safety
- ▶ Comes complete with padded carry bag to protect the scooter whilst in transit
- ▶ Excellent turning radius of 780mm
- ▶ Rear anti-tip wheels as standard
- ▶ Anti-roll back system safety device prevents scooter from rolling backwards on hills
- ▶ 24V 8.7ah lithium ion airline friendly batteries**
- ▶ Available in Copper, Red or Purple

Specifications

Width: 41cm (16")
Length: 89cm (35")
Height: 86 (34")
Max Speed: 4mph*
Max Climbing Angle: 10°
Ground Clearance: 1.7cm (0.67")
Batteries: 25V 8.7 Lithium Ion
Total Weight: 21.5kg (47lb)
Weight Capacity: 115kg (18st)

* Speed and range may vary depending upon user weight, type and incline of terrain, weather, battery charge and condition, operating speed and general driving situation.
**Batteries included as standard



Folds automatically to a compact size to easily fit into small car boots

Padded carry bag for extra protection during transportation



Dual Front Wheel for added stability



Leaflet Code: LL443

If you would like to view this product please contact your local stockist.

We have a continuous product improvement policy and consequently reserve the right to amend design and specification without prior notice. All sizes and weights stated are nominal. All errors and omissions excepted. Due to the limitations of photography and the printing process colours shown may not be 100% accurate.