

EDGE 3



Upgraded SRS (Smooth Ride Suspension) provides a smooth, comfortable ride and advanced stability. Optional iLevel® is available at 4.5 mph with up to 12" of power adjustable seat height. Standard features include front and rear LED fender lights, a USB charger and drive wheel colour accents.



GREEN MACHINE



Specifications

WEIGHT CAPACITY ⁴	21 stone 6 lbs. (136 kg)
MAXIMUM SPEED ¹	Up to 6 mph Up to 4.5 mph w/iLevel activated
GROUND CLEARANCE ²	69 mm (2.7")
TURNING RADIUS ²	521 mm (20.5")
OVERALL SIZE ^{2,6}	Length: 820 mm (32.23") without legrests Width: 622 mm (24.4")
DRIVE WHEELS	356 mm (14") pneumatic or solid available
CASTER WHEELS	Front: 152 mm (6") Rear: 152 mm (6")
SUSPENSION	SRS (Smooth Ride Suspension)
COMPONENT WEIGHTS ^{5,7}	Base: 138.6 lbs (63 kg) Batteries: 38 lbs (17.24 kg) each
DRIVETRAIN	Two motor, Mid-wheel 6@ drive
BRAKES	"Intelligent Braking" electronic regenerative, disc park brake
BATTERIES ³	Two 12-volt, deep-cycle 62 Ah (standard)
BATTERY CHARGER	8-amp, Off-board
MAXIMUM OBSTACLE CLIMBING ABILITY	75 mm (2.95")
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Tilt and iLevel®	Width: 12-22" Depth: 12-22"
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	16.5-18.5"
TRU-Balance® 3 Power Tilt ⁵	17.5-19.5"
TRU-Balance® 3 Tilt & iLevel® ⁵	17.5-19.5"

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Call for availability.


Features

- Upgraded SRS (Smooth Ride Suspension)
- iLevel® option available with up to 12" of seat elevation at walking speed (up to 4.5 mph)
- Fender lights and USB charger included
- Bold new colours

Options*

- iLevel® option allows 12" of seat elevation while at walking speed (4.5 mph)
- TRU-Balance® 3 Power Recline
- TRU-Balance® 3 Seating
- TRU-Balance® 3 Power Tilt

Colours*



* Please contact your representative for the full range of options and colours.