

Karma

mobility solutions

*Stylish, Lightweight,
Ergonomic Design.*



KM-9000

The KM-9000 is a lightweight modern, stylish wheelchair packed with innovative features and adjustments.

Features include ;

Karma's patented S-Ergo seating system, foldable handgrips and backrest angle adjustment.

These adjustments allow the user to set the wheelchair up to be as comfortable and easy to use as possible.

The folding frame makes it convenient to transport and store.



Pressure redistribution



Reduced sliding

FEATURES

1. Seat widths and depths:

14"x15", 16"x16", 18"x16"

14"x16", 16"x18", 18"x18"

2. Various configurations:


- Backrest height: 33cm~41cm or 41cm~48cm.
- Backrest angles: 85°, 90°, 95° (using tools).
- Seat to ground heights: 44cm~51cm (with different front and rear wheel combinations).
- Adjustable rear wheel position: Optimal driving performance can be achieved by adjusting the centre of gravity: 1. active (front) 2. standard (medium) 3. stable (rear)
- Angle adjustable footplate: One-piece footplate with length, height and angle adjustment

3. Folding, compact frame

4. Tension adjustable backrest



Adjustable rear wheel position



Backrest height & angle adjustment

OPTIONS

- Sport Rear wheels
- Magnesium alloy wheels
- Flip-back armrests(22cm height / 24cm height)
- Folding handgrip
- Anti-tippers
- Skirt guard



Folding handgrip

width	F/R Wheels												Max. User Weight
16"	5" - 24"	41	41	59	33	22 / 24	92	44~51	33~41 41~48	74~96	11.5 kg	115 kg	

• KARMA reserves the right to modify information herein without further notice •

± 1cm ±0.5kg

CE ISO 9001
Rev. 00 APR. 2012

KARMA Mobility Limited

www.karmamobility.co.uk e-mail:info@karmamobility.co.uk
Unit 6, Target Park, Shawbank Road, Redditch, B98 8YN U.K.
TEL: 0845-630-3436 FAX: 0845-630-3736