

<b>LENGTH</b>	<b>100cm / 39"</b>
<b>WIDTH</b>	<b>51cm / 20"</b>
<b>HEIGHT</b>	<b>89cm / 35"</b>
<b>FOLDED LENGTH</b>	<b>N/A</b>
<b>FOLDED WIDTH</b>	<b>N/A</b>
<b>FOLDED HEIGHT</b>	<b>N/A</b>
<b>SEAT WIDTH</b>	<b>43.5cm / 17"</b>
<b>SEAT DEPTH</b>	<b>39cm / 15.5"</b>
<b>FRONT WHEEL</b>	<b>Solid 20cm / 8"</b>
<b>REAR WHEEL</b>	<b>Solid 20cm / 8"</b>
<b>TOTAL WEIGHT</b>	<b>77lbs / 35kg</b>
<b>HEAVIEST PART</b>	<b>23.5lbs / 10kg</b>
<b>BATTERY WEIGHT</b>	<b>6.6lbs / 3kg</b>
<b>MAX USER WEIGHT</b>	<b>18stone / 115kg</b>
<b>RANGE</b>	<b>9miles / 15km</b>
<b>SPEED</b>	<b>4mph / 6kph</b>
<b>BATTERY</b>	<b>25.9v 11.6ah Lithium</b>
<b>MOTOR</b>	<b>24v 250w</b>
<b>GROUND CLEARANCE</b>	<b>3.5cm / 1.4"</b>
<b>MAX GRADIENT</b>	<b>6°</b>
<b>TURNING RADIUS</b>	<b>131cm / 52"</b>
<b>BOX DEPTH</b>	<b>108cm / 42.5"</b>
<b>BOX WIDTH</b>	<b>56cm / 22"</b>
<b>BOX HEIGHT</b>	<b>55cm / 26.6"</b>
<b>PACKAGED WEIGHT</b>	<b>83lbs / 38kg</b>

