

EXCEL ROADSTER DX8 Specification Sheet

Height adjustable headrest

Captain swivel seat with folding backrest

Width adjustable armrests

Angle adjustable Tiller

Removable basket with carry handle

Full LED lighting system

Spacious footwell for optimum seating posture



Simple dashboard with battery indicator

Rear LED lighting for improved visibility and extra safety

Unique two-tone shroud design with alloy wheels

Available Colours



Full Specifications

Seat Width	46 cm	Maximum Speed	8 mph
Seat Depth	40 cm	Range	35 km (22 miles)*
Seat Height (From Ground)	53 cm	Turning Radius	153 cm
Total Length	120 cm	Ground Clearance	9 cm
Total Width	61 cm	Battery	2 x 12V / 35Ah
Total Height	99 cm	Charger	24V / 4Ah
		Total Weight (Including Battery)	86 kg
		Maximum User Weight	136 kg (21 Stone)

KEY: * = Varies depending on weather and terrain conditions



www.vanosmedical.co.uk