

User Information	Dealer Information
Order No.	Order No.
Order Date:	Order Date:
Client Name:	Client Name:
Phone:	Phone:



Standard Features

- 130 kg Maximum Weight Capacity
- Self Propulsion with 20", 22" or 24" Q-Release Rear Wheels
- Transit Propulsion with 14" Fixed Rear Wheel
- 5", 6", 7" or 8" Caster Wheels
- Seat Widths: 380/ 405/ 430/ 455/ 480/ 505mm
- Seat Depths: 405/ 430/ 455mm
- Angle Adjustable Backrest: -4°/ 0°/ 4°/ 8°/ 12°/ 16°
- Backrest Heights: 430mm/ 510mm
- 80° Swing in/Swing out Legrest with Flip-up Footplates
- Pull to Lock Brakes

◆ Standard feature ○ Free of Charge

Seat

Seat Width and Seat Depth

Depth \ Width	Width				
	405mm(16")	430mm(17")	455mm(18")	480mm(19")	505mm(20")
405mm(16")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
455mm(18")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
430mm(17") (1" Extended Tube + Strip)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
455mm(18") (2" Extended Tube + Strip)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Seat-to-Floor Height

Front Caster \ Rear Wheel	5"	6"	7"	8"	Caster Position
	14"	<input type="checkbox"/> 18" / L2	<input type="checkbox"/> 18" / L1 <input type="checkbox"/> 19" / L3	<input type="checkbox"/> 19" / L2 <input type="checkbox"/> 20" / L3	
20"	<input type="checkbox"/> 15.5" / S1 <input type="checkbox"/> 16" / S2 <input type="checkbox"/> 17" / S2 <input type="checkbox"/> 18" / L2	<input type="checkbox"/> 18" / L1 <input type="checkbox"/> 19" / L3	<input type="checkbox"/> 19" / L2 <input type="checkbox"/> 20" / L3	<input type="checkbox"/> 20" / L2	
22"	<input type="checkbox"/> 15.5" / S1 <input type="checkbox"/> 16" / S2 <input type="checkbox"/> 17" / S2 <input type="checkbox"/> 18" / L2	<input type="checkbox"/> 18" / L1 <input type="checkbox"/> 19" / L3	<input type="checkbox"/> 19" / L2 <input type="checkbox"/> 20" / L3	<input type="checkbox"/> 20" / L2	
24"		<input type="checkbox"/> 19" / L3	<input type="checkbox"/> 19" / L2 <input type="checkbox"/> 20" / L3	<input type="checkbox"/> 20" / L2 <input type="checkbox"/> 21" / H2	

Seat Upholstery	
<input type="checkbox"/> Padded Seat Upholstery (Only available for seat depth with 1" or 2" extended tube)	○
Backrest	
Backrest	
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) 430mm Back Height Folding Backrest	○
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) 510mm Back Height Folding Backrest	◆
<input type="checkbox"/> Angle Adj. Backrest (-4° to +16°, by Step of 4°) 380mm -530mm Back Height Adj Fixed Backrest(Non holdable)	£60.00
Backrest Upholstery	
<input type="checkbox"/> Standard Backrest Upholstery	◆
<input type="checkbox"/> Backrest Upholstery for Height Adjustable Backrest	£80.00
<input type="checkbox"/> Strap Tension Adjustable Backrest (Including Detachable Padded Cover)	£50.00
Legrests/Footplates	
Legrest RH	
<input type="checkbox"/> 80° RH Swing in/Swing out Legrest	◆
<input type="checkbox"/> Elevating Legrest (0° to 80°) RH Swing in/Swing out	£80.00
<input type="checkbox"/> Amputee Legrest RH	£65.00
<input type="checkbox"/> 90° Straight Legrest RH w/ Flip-up RH Angle&Depth Adjustable Footplate	TBC
Legrest LH	
<input type="checkbox"/> 80° LH Swing in/Swing out Legrest	◆
<input type="checkbox"/> Elevating Legrest (0° to 80°) LH Swing in/Swing out	£80.00
<input type="checkbox"/> Amputee Legrest LH	£65.00
<input type="checkbox"/> 90° Straight Legrest LH w/ Flip-up LH Angle&Depth Adjustable Footplate	TBC
Footplate RH	
<input type="checkbox"/> Flip-up RH Tube-in-center Footplate	◆
<input type="checkbox"/> Flip-up RH Angle Adjustable Footplate (Not available for 15" seat width)	£20.00
Footplate LH	
<input type="checkbox"/> Flip-up LH Tube-in-center Footplate	◆
<input type="checkbox"/> Flip-up LH Angle Adjustable Footplate (Not available for 15" seat width)	£20.00
Armrests	
Armrest RH	
<input type="checkbox"/> Flip-back PU Armpad RH 220~240mm Height	◆
Armrest LH	
<input type="checkbox"/> Flip-back PU Armpad LH 220~240mm Height	◆
Front Castors/Forks	
Castors	
<input type="checkbox"/> 125mm x 30mm Solid Tyre (5")	○
<input type="checkbox"/> 150mm x 30mm Solid Tyre (6")	○
<input type="checkbox"/> 175mm x 30mm Solid Tyre (7")	○
<input type="checkbox"/> 200mm x 30mm Solid Tyre (8")	◆
<input type="checkbox"/> 200mm x 50mm (8" x 2") Solid Tyre (each)	£10.00
Forks	
<input type="checkbox"/> Standard Fork	◆
<input type="checkbox"/> 2" Wide Long Fork (While 200mm x 50mm Tire is ordered) for 8" castor	£10.00

Rear Wheels	
Rear Wheel	
<input type="checkbox"/> 14" Solid Tyre Plastic Rim	◆
<input type="checkbox"/> 14" Pneumatic Tyre Plastic Rim	○
<input type="checkbox"/> 20" Solid Tyre Spoked Rim with Ergo Handrim	○
<input type="checkbox"/> 20" Pneumatic Tyre Spoked Rim with Ergo Handrim	○
<input type="checkbox"/> 22" Solid Tyre Spoked Rim with Plastic Handrim	○
<input type="checkbox"/> 22" Pneumatic Tyre Spoked Rim with Plastic Handrim	○
<input type="checkbox"/> 24" Solid Tyre Spoked Rim with Ergo Handrim	◆
<input type="checkbox"/> 24" Solid Tyre Spoked Rim with Plastic Handrim	○
<input type="checkbox"/> 24" Pneumatic Tyre Spoked Rim with Ergo Handrim	○
<input type="checkbox"/> 24" Pneumatic Tyre Plastic Rim with Aluminum Handrim	○
Rear Wheels Axle	
<input type="checkbox"/> Fixed Axle (Only available for 14" rear wheel)	○
<input type="checkbox"/> Quick Release Axle (Standard version)	◆
Brakes	
<input type="checkbox"/> Pull to Lock Brakes (Only available for 20", 22" and 24" rear wheel)	◆
<input type="checkbox"/> Attendant Brakes	◆
Safety Options	
<input type="checkbox"/> Swing-away Height Adjustable Anti-tippers Pair	£30.00
Accessories	
<input type="checkbox"/> Lap Belt (Safety Belt)	£15.00
<input type="checkbox"/> Heel loops	£6.60
<input type="checkbox"/> Rigidifying Push Bar	£30.00
<input type="checkbox"/> Lap Tray	£45.00
<input type="checkbox"/> Cane Holder	£37.50
<input type="checkbox"/> Foldable Headrest	£65.00
Price	
Total price	