



reddot award 2014
winner

NEW



The colourful route to an effortless, safe drive

The 4 wheel *Invacare Colibri* is a stylish, simple and colourful micro lightweight scooter, designed for those who enjoy a totally independent lifestyle. It is now available in a 3 wheel option to offer additional leg room.

Thanks to the unique patented *Invacare LiteLock™* system the scooter can easily be taken apart without the need for tools and fits neatly into any car boot. Plus, with its small footprint the *Invacare Colibri* offers greater manoeuvrability in tight, or limited spaces.

Trusted Solutions, Passionate People



Yes, you can.®



Features and Options



LiteLock system - easy to take apart and reassemble

The ingenious and unique patented *Invacare LiteLock™* system automatically releases once the battery box has been removed. It takes less than a minute to take the scooter apart - without tools. Its front section can be simply lifted and placed in the boot of a car, along with the other sections. It's just as easy to reassemble and the scooter is secure.



Effortless Dismantling

It takes less than a minute to take the *Invacare Colibri* apart ready for transportation.



Option of a bigger wheel set and battery pack

Perfect for travelling longer distances, on bumpier terrain, such as cobbled streets. Choose where to charge your battery: directly on the *Colibri* or detached from it.



Increased comfort

The soft, padded seat provides comfort during driving. For increased comfort there is also the optional 'easy fit' seat suspension kit.

Technical data

<i>Invacare Colibri</i>	465 mm	400 mm	370 - 420 mm	385 mm	505 mm
4 wheels	1010 mm	44 kg (12 Ah) 49 kg (18 Ah)	136 kg	2 x 12 Ah 2 x 18 Ah	2200 mm
3 wheels	1010 mm	42 kg (12 Ah) 46 kg (18 Ah)	136 kg	2 x 12 Ah 2 x 18 Ah	1940 mm
<i>Invacare Colibri</i>	45 mm	6° / 10.5%	200 W	Up to 16 km	8 km/h



LiteLock™ system

Invacare Limited
Pencoed Technology Park -
Pencoed - Bridgend CF35 5AQ -
United Kingdom
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
Sales Order E-mail:
ordersuk@invacare.com
www.invacare.co.uk

Invacare Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS -
County Dublin - Ireland
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2016 Invacare International Sàrl
All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

Colibri - UK - 07/2016 - 1576204

Shroud colours




Sapphire Blue



Jasper Red



Jade Green



Marble White



Diamond Silver

1/ Compliant with the directive on the Restriction of the use of certain Hazardous Substances in electrical and electronic equipment - RoHS (2011/65/EC)



Follow us on Twitter:
[@Think_Mobility](https://twitter.com/Think_Mobility)



Like us on Facebook:
[Invacare UK](https://www.facebook.com/Invacare_UK)



Follow us on YouTube:
[Invacare Ltd - UK](https://www.youtube.com/Invacare_Ltd_UK)



Yes, you can.®