






permobil Order Form - Progeo JOKER

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	UK (English)	
User	Company stamp	
End user wheight (kg)		Rigid frame wheelchair CRASH TEST APPROVED
Language user manual	Responsible	£ 2592,00 Retail price - VAT not included

Wheelchair measures						
Seat width (LS)	☐ 33	☐ 36	☐ 39	☐ 42	☐ 45	☐ 48
Seat depth (PS)	☐ 35	☐ 37,5	☐ 40	☐ 42,5	☐ 45	☐ 47,5 ☐ 50
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.						
Footplate distance (DP)	 (cm) ±0.5cm - adjustment steps 1.5cm					
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)						
Setting	☒ Standard ☐ Prudential ☐ Active ☐ Extreme					
Backrest Height (TS)	Std. Value 33,5 (min 24.5-max 47) adjustment steps 1.5cm					
camber	☒ 0° ☐ 2° ☐ 4°					
Front height STD 49				Rear height STD 42		



	☆ Aluminum rear adjustment plate		
	☆ Removable mudguard ABS		
	☐ Removable ABS fibre straight sideguard		
	☐ Removable carbon fibre mudguard 03946 315,00 ☐ Removable Carbon fibre straight sideguard 08918 315,00		
	○ FULL CARBON rear adjustment plate with integrated mudguard	07646	939,00
	☆ Version mudguard (Available only with rear wheels 24")		
	☐ Version straight sideguard		

JOKER		☆ Aluminium footplate, rake adjustable Minimum footplate distance 35cm. ☆ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External		11523	no list
		○ Aluminium footplate EVOLUTION - rake and depth adjustable Minimum footplate distance 35cm. ☆ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External		07649	126,00
		○ One piece Carbon fibre footplate JOKER Minimum footplate distance 35cm. not available for seat width 48		03348	252,00
		○ One piece Carbon fibre footplate, rake adjustable Minimum footplate distance 35cm. ☆ Internal ☐ External		03346	298,00
JOKER ZERO		○ Plastic footrest, rake adjustable Minimum footplate distance 15cm. ☆ 2/3 Internal ☐ 2/3 External		06814	189,00
		○ One piece aluminium footplate - Flip-up!, rake and depth adjustable Minimum footplate distance 15cm. ☆ Internal ☐ 2/3 Internal ☐ External ☐ 2/3 External		06856	189,00
		○ One piece Carbon fibre footplate rake adjustable, swing away, flip-up. Minimum footplate distance 15cm. ☆ Internal ☐ External		06857	487,00

Date Order User Address

Front frame

☆ Aluminium standard frame



Colour aluminium frame

☆ Bright silver

☐ Polished black
☐ Opaque black☐ Bright blue
☐ Bright anthracite☐ Bright orange
☐ Bright fuchsia
☐ Bright red☐ Bright lime
☐ Bright blue china
☐ Pearl white

Seat

☆ Seat cloth with pocket (under seat)

☐ Seat with tension straps

08883 224,00

☐ Carbon rigid seat

08336 475,00

☐ Pocket (under seat) [25x11x4,5]

08281 74,00

Rear wheels axle

☆ Aluminum rear wheels axle

☐ Carbon fiber rear wheels axle

(Lightest configuration)

07647 148,00


Backrest adjustment and fixing system

☆ Backrest with double locking (upright and folded)
with rake adjustable from -4° to +16°☐ Backrest with upright locking system ONLY (no locking when folded)
with rake adjustable from -4° to +16°☐ Fixed back(non-folding back)
with rake adjustable from -4° to +16°

Backrest

☆ Adjustable velcro straps and padded transpiring upholstery

☆ Aluminium backrest

☐ Titanium backrest☐ Carbon fibre backrest

06811 189,00

04389 749,00

☐ **AIR** Integrated posture backrest (NOT REMOVABLE)

Carbon shell - Lateral retaining 7,5cm

Depth adjustment ±1cm. - Angle adjustment ±10°

not available for seat width 48

1.047,00

H 28		H 34		H 40	
Backrest Height (TS) min 34 (.....)		Backrest Height (TS) min 38 (.....)		Backrest Height (TS) min 44 (.....)	
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
☐	☐	☐	☐	☐	☐

Installation on:

☆ Backrest tubes

☐ Aluminium backrest☐ Titanium backrest☐ Carbon fibre backrest

06811 189,00

04389 749,00

Date	Order	User	Address																																																																						
Handle tube ☆ Tubes with standard handle ☐ Tubes with short handle ☐ SPORT tubes Without push handles ☐ Carbon fibre tubes with short handle (Lightest configuration) 04390 148,00 Not available with: Posture backrest PHYSIO - Posture backrest A I R ☐ Folding handles 02958 252,00 																																																																									
Casters ☆ 4" Soft ☐ 3" Sport (Lightest configuration) ☐ 5" Soft ☐ 6" Soft (Available only with 100° angle) ☐ 6" Pneumatic (Available only with 100° angle) ☐ 3" Active Soft (wide) (Aluminium alloy rim) 06824 126,00 ☐ 4" Active Sport (narrow) (Aluminium alloy rim) 01398 126,00 ☐ 4" Active Soft (wide) (Aluminium alloy rim) 03943 126,00 ☐ 5" Active Soft (wide) (Aluminium alloy rim) 07518 126,00 ☐ 5" Pneumatic 10844 63,00																																																																									
Colour set: Rear wheels and forks - To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24) ☆ Black 22" 24" 25" 26" 00833 126,00 ☐ Bright silver 22" 24" 00833 126,00 ☐ Bright anthracite 24" 00833 126,00 ☐ Bright orange 22" 24" 00833 126,00 ☐ Bright fuchsia 22" 24" 00833 126,00 ☐ Bright blue 22" 24" 00833 126,00 ☐ Bright red 22" 24" 00833 126,00 ☐ Bright lime 22" 24" 00833 126,00 ☐ Pearl white 22" 24" 00833 126,00																																																																									
Rear wheels ☆ 24" Ultra-light Progeo - Standard! ☐ 22" Ultra-light Progeo - Standard! ☐ 25" Ultra-light Progeo - Standard! (colour: black) ☐ 26" Progeo - Standard! (colour: black) ☐ 24" ULTRA-LIGHT PROGEO KEY SPORT Weight reduction approx. 450gr. 04383 876,00 ☐ 24" Ultra-light "Spinergy LX" 01685 876,00 ☐ 24" Ultra-light "Spinergy LXL" ☐ Black hub and white spokes ☐ Hub and spokes black 07846 876,00 ☐ 24" Ultra-light "Spinergy Flex rim" Integrated aluminium/rubber handrim 02270 1.374,00 ☐ 25" Ultra-light "Spinergy LX" 02101 876,00 ☐ 25" Ultra-light "Spinergy LXL" ☐ Black hub and white spokes 07847 876,00 ☐ 26" Ultra-light "Spinergy LX" 02102 876,00 ☐ 24" Off-road (aluminium handrim black) 01392 224,00																																																																									
Rear wheel tyre ☆ High pressure tyre (22"-24"-25"-26") (Lightest configuration) ☐ High pressure profile tyre (24") ☐ Marathon Plus (24"-25") with anti-puncture protection 01669 63,00 ☐ Anti-puncture - Colour Grey (22") 04181 98,00 ☐ Anti-puncture - Colour Black (24") 11083 126,00																																																																									
<table border="1"> <thead> <tr> <th>Handrim</th> <th>Wheel distance ☆ 3 cm</th> <th>□ 4 cm</th> <th></th> <th></th> </tr> </thead> <tbody> <tr> <td>☆ Aluminium</td> <td>(22"-24"-25"-26")</td> <td>(Lightest configuration)</td> <td></td> <td></td> </tr> <tr> <td>○ Aluminium black SOFT-TOUCH</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>09170</td> <td>126,00</td> </tr> <tr> <td>○ Chromed Steel</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>01684</td> <td>148,00</td> </tr> <tr> <td>○ Titanium</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>00016</td> <td>275,00</td> </tr> <tr> <td>○ Ultra Grepp</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>00017</td> <td>315,00</td> </tr> <tr> <td>○ Max Grepp</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>00834</td> <td>349,00</td> </tr> <tr> <td>○ Ergo Para</td> <td>(24")</td> <td></td> <td>00632</td> <td>349,00</td> </tr> <tr> <td>○ TOP FRICTION</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>11023</td> <td>399,00</td> </tr> <tr> <td>○ NOVA H</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>11025</td> <td>263,00</td> </tr> <tr> <td>○ NOVA D</td> <td>(24"-25")</td> <td>(24" handrim, on 25" wheel is installed with spacers)</td> <td>11026</td> <td>286,00</td> </tr> <tr> <td>○ SURGE LT</td> <td>(25")</td> <td>Oval countour with a rubber insert</td> <td>06427</td> <td>360,00</td> </tr> <tr> <td>○ Quad (24")</td> <td>(24")</td> <td></td> <td>00693</td> <td>200,00</td> </tr> <tr> <td>○ Silicon cover for standard aluminium handrim</td> <td></td> <td>Available for measures (22" - 24" - 26")</td> <td>04050</td> <td>160,00</td> </tr> </tbody> </table>				Handrim	Wheel distance ☆ 3 cm	□ 4 cm			☆ Aluminium	(22"-24"-25"-26")	(Lightest configuration)			○ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	126,00	○ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	148,00	○ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	275,00	○ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	315,00	○ Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	349,00	○ Ergo Para	(24")		00632	349,00	○ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	399,00	○ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	263,00	○ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	286,00	○ SURGE LT	(25")	Oval countour with a rubber insert	06427	360,00	○ Quad (24")	(24")		00693	200,00	○ Silicon cover for standard aluminium handrim		Available for measures (22" - 24" - 26")	04050	160,00
Handrim	Wheel distance ☆ 3 cm	□ 4 cm																																																																							
☆ Aluminium	(22"-24"-25"-26")	(Lightest configuration)																																																																							
○ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	126,00																																																																					
○ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	148,00																																																																					
○ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	275,00																																																																					
○ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	315,00																																																																					
○ Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	349,00																																																																					
○ Ergo Para	(24")		00632	349,00																																																																					
○ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	399,00																																																																					
○ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	263,00																																																																					
○ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	286,00																																																																					
○ SURGE LT	(25")	Oval countour with a rubber insert	06427	360,00																																																																					
○ Quad (24")	(24")		00693	200,00																																																																					
○ Silicon cover for standard aluminium handrim		Available for measures (22" - 24" - 26")	04050	160,00																																																																					

Date	Order	User	Address
------------	-------------	------------	---------------

Wheel locks

☆ Bent lever push to lock 

☐ Straight lever push to lock 

☐ Straight lever pull to lock 

☐ Sport brake (Lightest configuration)  06428 **126,00**

☐ Brake with extension lever push to lock  01949 **63,00**

☐ Brake with extension lever pull to lock  03372 **63,00**

Armrests

☐ Tubular armrests (Swing away, removable and height adjustable armrest)  00694 **252,00**

Accessories

☐ Pad for front frame protectors 07058 **52,00**

☐ Backrest-seat connection cloth 10241 **17,00**

☐ Removable storage bag (with quick coupling) 09811 **126,00**

SAFETY

☐ 4 points tie-down system rings  08654 **189,00**

Optional 1

Only one optional devices for side can be chosen

☐ Transit wheels Exclude the installation of Anti-tippers  03290 **148,00**

☐ Anti-tipper BASE Exclude the installation of Transit wheels

☐ Right 06819 **98,00**

☐ Left 06820 **98,00**

☐ Anti-tipper with rapid hookup (inside rotation to reduce the encumbrance) Exclude the installation of Transit wheels

☐ Right 06821 **126,00**

☐ Left 06822 **126,00**

Optional 2

☐ Push handles for companion (Removable and height adjustable)  03645 **177,00**

☐ Tetraclip (Quick release axle for quadriplegic) 03947 **177,00**

☐ Safety pelvic belt 03972 **63,00**

☐ Harness 04103 **126,00**

☐ Lumbar support (to insert in the backrest upholstery) 00013 **40,00**

☐ Swing away lateral support Not available with Posture backrest PHYSIO - Posture backrest A I R

☐ Right 02097 **252,00**

☐ Left 02098 **252,00**

☐ Spoke guards Black (22"-24") 01043 **148,00**

☐ Spoke guards Transparent (22"-24") 00835 **148,00**

☐ Spoke guards Progeo design (22"-24") 00019 **177,00**

☐ Spoke guards Dragon (22"-24") 02106 **177,00**

Additional Kit

☐ 24" Off-road rear wheels Off-Road tyre - aluminium handrim (black colour) 03378 **498,00**

☐ 24" FAT rear wheels FAT tyre 3" - aluminium handrim (black colour) 09357 **876,00**

Date Order User Address

PHYSIO - Posture backrest (Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment ($\pm 18^\circ$)

- ☐ **AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7.5 cm		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.047,00

- ☐ **POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape	Lateral retaining wing B=regular		
	R	R	R
7.5 cm			
	☐	☐	☐

1.013,00

Configuration DEEP with side support PAD (height adjustable)

- ☐ Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)


☐ Right
10272

29,00


☐ Left
10272

29,00

Large (h 20cm)


☐ Right
10273


40,00


☐ Left
10273

40,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.



- ☐ Large (18x12)

05647 183,00

- ☐ Small (14x9)

05643 183,00

- ☐ Lateral support extension

Quantity (number pcs) (.....)

10271 69,00

- ☐ Lateral support extension

Quantity (number pcs) (.....)

- ☐ Small (14x9)

05642 183,00

- ☐ Large (18x12)

05646 183,00

- ☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio
Not available for backrest Physio AIR carbon shell



- ☐ Large (34x14)
☐ Small (28x14)

05648 252,00

Permobil LimitedUnit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts. AL7 1ANUnited Kingdom - Phone: +44 1484 722888

Date	Order	User	Address
------	-------	------	---------

Note

Wheelchair : Progeo JOKER - Rev. 2023.1	Prices - VAT excluded		Send this page only if any optional have been chosen	
	to be added to the price list		☐	standard (supplied if no requests are indicated)
			☐	on request (without extra charge)
			☐	optional (cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

REFERENCE TABLES

1	Front height				
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster					
FORK		Casters			
		3"	4"	5"	6"
SPORT		max 52	max 53	max 54	
SMALL		max 54	max 55	max 56	max 57
MEDIUM		max 56	max 57	max 58	max 59
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.					

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 36	min 38	min 39	min 40	
max 44	max 46	max 47	max 48	

3	Footplate option	Footplate distance
	Aluminium footplate rake adjustable	min 35
	One piece Carbon fibre footplate	min 35
	One piece Carbon fibre footplate rake adjustable	min 35
	Plastic footrest, rake adjustable	min 35
	One piece aluminium footplate - Flip-up!, rake and depth adjustable	min 35
	One piece Carbon fibre footplate rake adjustable, swing away, flip-up.	min 35
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm.)		
Calculated with the best configuration and with the footplates positioned on inside.		
NOTE! Lower footplate distances are possible only after our valuation		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all couloers)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)