



Date	UK (English)	
User	Company stamp	
End user weight (kg)		Rigid frame wheelchair
Language user manual	Responsible	£ 3205,00
		Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm	
Seat width (LS) ☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48	The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1.5 cm wider.
Seat depth (PS) ☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50	The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.
Footplate distance (DP) (cm) ± 0.5 cm - adjustment steps 1.5cm	(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)
Setting ☒ Standard ☐ Prudential ☐ Active ☐ Extreme	
Backrest Height (TS) Std.Value 33,5 (min 26-max 47) adjustment steps 1.5cm	
Camber ☒ 0° ☐ 2° ☐ 4°	

☒ Standard rear adjustment	☒ Carbon fibre mudguard	
	☐ Carbon fibre Straight sideguard	
	☐ Removable version	
		 08921 130,00
☐ FULL CARBON rear adjustment plate with integrated sideguard Lightest configuration		07646 967,00
	☒ Version mudguard (Available only with rear wheels 24")	
	☐ Version straight sideguard	

Front height - Tick the box with the required height

The possible difference in height between HIGH frame and LOW frame is 3cm (with same fork, caster and assembly hole).

	 SPORT fork (H95)	 Small size fork (H115)
LOW FRAME 90° Front height (min 44,5 - max 50) Footplate distance (min 33) U-TUBE Aluminium footplate Footplate distance (min 18 ÷ max 35) DPmax = Front height - 10cm	Casters 3" ☐ ~45 ☐ ~46 4" ☐ ~46 ☐ ~47 5" ☐ ~48,5	Casters 3" ☐ ~45,5 ☐ ~47 4" ☐ ~46,5 ☐ ~48 ☐ ~49 5" ☐ ~49 ☐ ~50
HIGH FRAME 90° Front height (min 47,5 - max 53) Footplate distance (min 36) U-TUBE Aluminium footplate Footplate distance (min 21 ÷ max 38) DPmax = Front height - 10cm	Casters 3" ☐ ~48 ☐ ~49 4" ☐ ~49 ☐ ~50 5" ☐ ~51,5	Casters 3" ☐ ~48,5 ☐ ~50 4" ☐ ~49,5 ☐ ~51 ☐ ~52 5" ☐ ~52 ☐ ~53

The real front-seat height could be slightly different from the chosen one considering it is affected by seat depth and rear seat height.

Rear height

Rear wheels 22" (min 36 - max 44)
Rear wheels 24" (min 38 - max 46)
Rear wheels 25" (min 39 - max 47)
Rear wheels 26" (min 41 - max 48)

The difference between Front Height and Rear Height must be in the 3 - 12 cm range.

Date Order User Address

Seat

- ☆ Seat cloth with pocket (under seat)
- ☐ Seat with tension straps
- ☐ Carbon rigid seat
- ☐ Pocket (under seat) [25x11x4,5]

08883	231,00
08336	489,00
08281	76,00

CUSTOMIZE COLORS

- ① **Frame colour**
(without extra charge)
- ☆ **Opaque black**
- ☐ Pearl white
- ☐ Matt orange
- ☐ Tekno anthracite
- ☐ Matt Midnight blue
- ☐ Matt grey
- ☐ Matt red

Frame colour on request

- ☐ Bright silver
- ☐ Bright red
- ☐ Bright blue
- ☐ Bright fuchsia
- ☐ Bright lime

12579 105,00

- ② **Fork colour**
- ☆ **Opaque black**
- ☐ Bright silver
- ☐ Pearl white
- ☐ Bright red
- ☐ Tekno anthracite
- ☐ Matt grey

08342 65,00


- ☐ Bright blue
- ☐ Bright fuchsia
- ☐ Bright lime
- ☐ Matt Midnight blue
- ☐ Matt red

- ③ **Colour rear wheels**
To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)
- ☆ **Black**

- ☐ Bright silver
- ☐ Pearl white
- ☐ Bright red
- ☐ Tekno anthracite
- ☐ Matt grey

04692 65,00

- ☐ Bright blue
- ☐ Bright fuchsia
- ☐ Bright lime
- ☐ Matt Midnight blue
- ☐ Matt red



Rear wheels axle

- ☆ Aluminum rear wheels axle
- ☐ Carbon fiber rear wheels axle

Lightest configuration

07647 152,00

Backrest adjustment and fixing system

- ☆ Backrest with double locking (upright and folded)
- ☐ Backrest with upright locking system ONLY (no locking when folded)
- ☐ Fixed back(non-folding back)

with rake adjustable from -4° to +16°

with rake adjustable from -4° to +16°

with rake adjustable from -4° to +16°

Backrest

- ☆ Adjustable velcro straps and padded transpiring upholstery

- ☆ Aluminium backrest

- ☐ Titanium backrest

- ☐ Carbon fibre backrest

06811 195,00

04389 771,00

- ☐ **AIR** Integrated posture backrest (NOT REMOVABLE)

Carbon shell - Lateral retaining 7,5cm

Depth adjustment ±1cm. - Angle adjustment ±10°

not available for seat width 48

1.078,00

H 28		H 34		H 40	
Backrest Height (TS) min 34		Backrest Height (TS) min 38		Backrest Height (TS) min 44	
(.....)		(.....)		(.....)	
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
☐	☐	☐	☐	☐	☐

Installation on:

- ☆ Backrest tubes


- ☐ Aluminium backrest

- ☐ Titanium backrest

- ☐ Carbon fibre backrest

06811 195,00

04389 771,00



Date Order User Address

Handle tube

☆ Tubes with standard handle

- ☐ Tubes with short handle
☐ SPORT tubes Without push handles

☐ Carbon fibre tubes with short handle *Lightest configuration*
Not available with: Posture backrest PHYSIO - Posture backrest A I R

04390 **152,00**

☐ Folding handles **Minimum backrest height 30,5cm**
Not available with: Push handles for companion

02958 **260,00****Casters**

☆ 4" Soft

- ☐ 3" Sport
☐ 5" Soft

Lightest configuration

- ☐ 3" Active Soft (wide) (Aluminium alloy rim)
☐ 4" Active Sport (narrow) (Aluminium alloy rim)
☐ 4" Active Soft (wide) (Aluminium alloy rim)
☐ 5" Active Soft (wide) (Aluminium alloy rim)
☐ 5" Pneumatic

06824 **130,00**01398 **130,00**03943 **130,00**07518 **130,00**10844 **65,00****Rear wheels**

☆ 24" Ultra-light Progeo - Standard!

- ☐ 22" Ultra-light Progeo - Standard!
☐ 25" Ultra-light Progeo - Standard!
☐ 26" Progeo - Standard!

(colour: black)
(colour: black)

- ☐ 24" ULTRA-LIGHT PROGEO KEY SPORT *Lightest configuration*
☐ 24" Ultra-light "Spinergy LX"
☐ 24" Ultra-light "Spinergy LXL"
☐ 24" Ultra-light "Spinergy Flex rim"
☐ 25" Ultra-light "Spinergy LX"
☐ 25" Ultra-light "Spinergy LXL"
☐ 26" Ultra-light "Spinergy LX"
☐ 24" Off-road (aluminium handrim black)

*Lightest configuration*04383 **902,00**01685 **902,00**07846 **902,00**02270 **1.415,00**02101 **902,00**07847 **902,00**02102 **902,00**01392 **231,00**☐ Black hub and white spokes☐ Hub and spokes black

Integrated aluminium/rubber handrim

☐ Black hub and white spokes**Rear wheel tyre**☆ High pressure tyre (22"-24"-25"-26") *Lightest configuration*☐ High pressure profile tyre (24")

☐ Marathon Plus (24"-25") with anti-puncture protection

01669 **65,00**

☐ Anti-puncture - Colour Grey (22")

04181 **101,00**

☐ Anti-puncture - Colour Black (24")

11083 **130,00****Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"****Handrim****Wheel distance ☆ 3 cm**☐ 4 cm

☆ Aluminium

(22"-24"-25"-26")

Lightest configuration☐ Aluminium black SOFT-TOUCH

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

09170 **130,00**☐ Chromed Steel

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

01684 **152,00**☐ Titanium

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00016 **283,00**☐ Ultra Grepp

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00017 **324,00**☐ Max Grepp

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00834 **359,00**☐ TOP FRICTION

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11023 **411,00****Oval countour with a rubber insert**☐ NOVA H

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11025 **271,00****Ergonomic black profile**☐ NOVA D

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11026 **295,00****Ergonomic profile with rubber coating (ideal for quadriplegics)**☐ SURGE LT

(25")

Oval countour with a rubber insert

06427 **371,00**☐ Quad (24")

(24")

00693 **206,00**☐ Silicon cover for standard aluminium handrim

Available for measures (22"- 24" - 26")

04050 **165,00****Wheel locks**

☆ Bent lever push to lock

- ☐ Straight lever push to lock
☐ Straight lever pull to lock

☐ Sport brake☐ Brake with extension lever push to lock☐ Brake with extension lever pull to lock06428 **130,00**01949 **65,00**03372 **65,00**

Date Order User Address

Footplate

☆ Aluminium footplate. Rake not adjustable

☆ Aluminium footplate


☐ with plastic plate☐ One piece aluminium footplate - Flip-up! FLIP-BACK! - rake and depth adjustable

IMPORTANT NOTICE! From the footplate to the floor minimum 6cm !! (Minimum Front height = DP + 4cm)

With 5" front wheel the difference between Front Height and Rear Height must be greater than 5cm

☐ Aluminium footplate - rake and depth adjustable

☆ Internal

☐ 2/3 Internal☐ External☐ 2/3 External

09420 195,00




04456 260,00

☐ U-TUBE - Aluminium footplate, rake and depth adjustable**LOW FRAME: minimum footplate distance 18 - maximum footplate distance 35****HIGH FRAME: minimum footplate distance 21 - maximum footplate distance 38**

09853 231,00

☐ One piece Carbon fibre footplate. Rake adjustable

☆ Internal

☐ External

04110 383,00

Armrests☐ Tubular armrests

(Swing away, removable and height adjustable armrest)



00694 260,00

☐ Removable armrest with carbon-fiber straight sideguard

Available only with: Standard rear adjustment



09079 231,00

Accessories☐ Pad for front frame protectors

07058 54,00

☐ Backrest-seat connection cloth

10241 18,00

☐ Removable storage bag (with quick coupling)

09811 130,00

Optional 1Only one optional devices for side can be chosen☐ Transit wheels

Exclude the installation of Anti-tippers



03290 152,00

☐ Anti-tipper with rapid hookupExclude the installation of Transit wheels
(inside rotation to reduce the encumbrance)☐ Right☐ Left

06821 130,00

06822 130,00

Optional 2☐ Push handles for companion

(Removable and height adjustable)



03645 182,00

☐ Tetraclick

(Quick release axle for quadriplegic)

03947 182,00

☐ Safety pelvic belt

03972 65,00

☐ Harness

04103 130,00

☐ Lumbar support

(to insert in the backrest upholstery)

00013 41,00

☐ Swing away lateral support☐ Right☐ Left

02097 260,00

02098 260,00

Not available with:

Posture backrest PHYSIO - Posture backrest A I R

☐ Spoke guards Black

(22"-24")

01043 152,00

☐ Spoke guards Transparent

(22"-24")

00835 152,00

☐ Spoke guards Progeo design

(22"-24")

00019 182,00

☐ Spoke guards Dragon

(22"-24")

02106 182,00

Additional Kit☐ 24" Off-road rear wheels

Off-Road tyre - aluminium handrim (black colour)

03378 513,00

☐ 24" FAT rear wheels

FAT tyre 3" - aluminium handrim (black colour)

09357 902,00

Date.....Order.....User.....Address.....

PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- ☐ *AIR* Carbon shell
- not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape 7,5 cm 		Lateral retaining wing A=high - R=regular					
		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.078,00

- ☐ *POSTURE* Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 7.5 cm	Lateral retaining wing R=regular		
	R	R	R

1.043,00

Configuration DEEP with side support PAD (height adjustable)


- ☐ Aluminum lateral with padded, breathable upholstery.



Small (h 13cm)		Large (h 20cm)	
	☐ Right 10272		☐ Right 10273
	30,00		41,00
	☐ Left 10272		☐ Left 10273
	30,00		41,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.



☐ Large (18x12)	05647	188,00
☐ Small (14x9)	05643	188,00
☐ Lateral support extension	Quantity (number pcs).....	
☐ Lateral support extension	10271	71,00
☐ Small (14x9)	05642	188,00
☐ Large (18x12)	05646	188,00

- ☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio Posture



☐ Large (34x14)	05648	295,00
☐ Small (28x14)		

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields Welwyn Garden City, Herts. AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo JOKER R2 - Rev. 2024.1	Prices - VAT excluded to be added to the price list	Send this page only if any optional have been chosen	
		☐ standard	(supplied if no requests are indicated)
		☐ on request	(without extra charge)
		☐ optional	(cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb