

Deliver to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:



Deliver to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Order date:
 Purchase order:
 Number required:

Marked for:
 Sunrise order no:

☐ Order ☐ Calculation ☐ Built 4 me (longer leadtime)

☒ = Standard ☐ = Cost option
☐ = No cost option

Reimbursement code:

Max. user weight: 125 kg

Part No. 770 00 001
 KRY01

USER WEIGHT kg

Mandatory information to ensure correct configuration

Base Price £4,250

KRYPTON R: STANDARD CONFIGURATION

KRYPTON R

Base specification includes:

- Tubular Footrest
- Aluminium backrest
- Removable side guards
- Carbon axle tube
- Exo back upholstery with colours
- Carbotechnology fork in colours
- Aluminium rim castor with colours
- Lightweight wheels in anodised colours
- Folding backrest
- Stickers in colour options



CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - in mm; measured from outside frame tube to outside frame tube

320	340	360	380	400	420	440	460
☐	☐	☐	☐	☐	☐	☐	☐

SEAT DEPTH SELECTION - in mm

340	360	380	400	420	440	460	480
☐	☐	☐	☐	☐	☐	☐	☐

FRONT SEAT HEIGHT SELECTION - in dependance of front castor size, fork, frame height version

The high frame is 40 mm higher than the low frame.

		Front seat height in mm														
		Fork (length)	Frame:	430	440	450	460	470	480	490	500	510	520	530	540	550
3"	111 mm	low	☐	☐	☐											
	123 mm				☐	☐	☐									
	Frog Leg						☐	☐								
	One Arm		☐	☐	☐											
	111 mm	high					☐	☐	☐							
	123 mm								☐	☐	☐					
	Frog Leg										☐	☐				
	One Arm						☐	☐	☐							
4"	111 mm	low		☐	☐	☐	☐									
	123 mm				☐	☐	☐	☐								
	Frog Leg						☐	☐	☐							
	One Arm			☐	☐	☐	☐									
	111 mm	high						☐	☐	☐	☐					
	123 mm								☐	☐	☐	☐				
	Frog Leg										☐	☐	☐			
	One Arm							☐	☐	☐	☐					
5"	111 mm	low				☐	☐	☐								
	123 mm						☐	☐	☐							
	Frog Leg							☐	☐	☐						
	One Arm				☐	☐	☐									
	111 mm	high								☐	☐	☐				
	123 mm										☐	☐	☐			
	Frog Leg											☐	☐	☐		
	One Arm									☐	☐	☐				
6"	123 mm	high										☐	☐	☐		
	Frog Leg												☐	☐	☐	

REAR SEAT HEIGHT SELECTION - in dependance of rear wheel size

The rear seat height is adjustable ± 10 mm.

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

Caution: with rear seat heights 460 mm and above, the width between the sideguards increases by 35 mm

Rear seat height in mm

Rear Wheels	380	390	400	410	420	430	440	450	460	470	480	490	500
24"	☐ ⁽¹⁾	☐ ⁽²⁾	☐ ⁽²⁾	☐	☐	☐	☐	☐	☐	☐	☐	☐ ⁽³⁾	
25"		☐ ⁽¹⁾	☐ ⁽²⁾	☐ ⁽²⁾	☐	☐	☐	☐	☐	☐	☐	☐	☐ ⁽³⁾

(1): No adjustment of rear seat-height possible

(2): No downwards-adjustment of rear-seat-height possible. Only +1 cm adjustable.


(3): No upwards-adjustment of rear-seat-height possible. Only -1 cm adjustable.

										£
SEAT UPHOLSTERIES AND CUSHIONS										
Seat Sling										
KRR020012	☐ Seat sling "airflow" with 1 utility bag									
KRR020019	☐ Seat sling "airflow" superlight. Without bags, max user weight 85kg (1) (L)									
KRR020017	☐ Seat sling "airflow" with 2 utility bags									345
KRR020018	☐ Seat sling "airflow" with utility bag and catheter bag (from SW 340mm)									410
KRR020000	☐ Without seat sling									
Seat cushion black cover										
Foam: (medium soft)		Thickness:		KRR020101	☐ 30 mm	KRR020102	☐ 50 mm			85
Latex: (very soft)		Thickness:		KRR020103	☐ 50 mm					85
JAY Cushions - for complete product range see Jay order forms										
JAY020002	☐ JAY Basic									75
JAY020001	☐ JAY Easy Visco, Curved Base									180
JAY020007	☐ JAY Easy Fluid, Curved Base									415
JAY020004	☐ JAY Lite									290
Seat accessories										
KRR090018	☐ Positioning belt (only with Standard backrest)									70
BACKREST										
KRR030001	☐ Standard backrest (aluminium, black) Choose back height factory setting (mm). Height adjustments possible in the range 250-300; 300-400; 400-450 mm. 250 - 300 300 - 350 350 - 400 400 - 450 ☐ 250 ☐ 300 ☐ 350 ☐ 400 ☐ 275 ☐ 325 ☐ 375 ☐ 425 ☐ 300 ☐ 350 ☐ 400 ☐ 450									
KRR030006	☐ Carbon backrest - fixed height, crosstube very close to the upholstery, "Krypton" rubber grip (2) (L) Choose back height (mm) no height adjustment possible. ☐ 250 ☐ 325 ☐ 400 ☐ 275 ☐ 350 ☐ 425 ☐ 300 ☐ 375 ☐ 450									300
KRR030102	☐ Backrest, angle-adjustable, fold-down, locked when folded, colour: (-31° to +15°) Colour: ☒ Black ☐ Orange ☐ Silver ☐ Red ☐ Blue KRR030036 KRR030037 KRR030038 KRR030039									90
KRR030012	☐ Backrest, angle-adjustable (non foldable) (-31° to +15°) (only in black anodizing)									
KRR030015	☐ Backrest, angle-adjustable, fold-down (-31° to +15°) (only in black anodizing)									
Push handles: choose only one option										
KRR030200	☐ Omit push handles (L)									
KRR030203	☐ Fold-down push handles (from 325mm backrest-height)									165
KRR030201	☐ Push handles long									140
KRR030204	☐ Height adj. push handles (+) (only with Standard backrest; from 325mm backrest height)									220
KRR030205	☐ Height adj. push handles for Jay backs (+) (only with Standard backrest; from 325mm backrest height)									220
Backrest upholstery										
KRR030005	☐ EXO superlight Backrest sling, synthetic fibre, padded, adjustable, black (Evo Exo Version) (L) Colour, binding tape: ☒ Black ☐ Orange ☐ Silver ☐ Red ☐ Blue									
KRR030316	☐ EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black (Evo Exo Version) Colour, binding tape: ☒ Black ☐ Orange ☐ Silver ☐ Red ☐ Blue									
KRR030317	☐ EXO PRO Backrest, synthetic fibre, padded, adj., breathable, spacer fabric, 1 accessory bag, black (Evo Exo Version) Colour, binding tape: ☒ Black ☐ Orange ☐ Silver ☐ Red ☐ Blue									90
KRR030000	☐ Without backrest upholstery									
(1): Max. user weight 85 kg, only available up to seat-width 400 mm (2): Not with height adj push handles, not with Jay J3 backs, not with positioning-belt (L)= lightest option (+)= adds 500 g or more to the wheelchair-weight										

				£
JAY 3 BACKS AND ACCESSORIES - only with Standard backrest (+)				
Jay 3 Carbon Back - Shallow Contour				
Lower Thoracic:	☐ JAY060001	Short - 170 mm	☐ JAY060002	Tall - 240 mm 820
Mid Thoracic:	☐ JAY060003	Short - 300 mm	☐ JAY060004	Tall - 380 mm 820
Accessories for JAY 3 Backs				
WHT030415	☐	Whitmyer Headrest, COBRA XTRA, flip back	(only with Jay 3 backs)	345
WHT030416	☐	Whitmyer Headrest, AXYS, Standard	(only with Jay 3 backs)	325
JAY090058	☐	JAY Accessory bag	(only with Jay 3 Carbon backs)	60
SIDEGUARDS AND ARMRESTS				
KRR040107	☐	Sideguard, composite, removable, 21 cm high (31 cm long), black		
KRR040101	☐	Sideguard, carbon fibre, slim, black	(1) (L)	260
KRR040102	☐	Sideguard, aluminium, cold weather protection (colour: black, matt)		230
KRR040103	☐	Sideguard, aluminium, fender, cold weather protection (colour: black, matt)	(min rear seat height 390 mm)	315
KRR040104	☐	Sideguard, carbon fibre, fender, cold weather protection, black	(min rear seat height 390 mm)	395
KRR040052	☐	Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable	(2) (+)	180
KRR040117	☐	Tubular armrest, Alu, swing away (sideways), padded, removable	(+)	160
KRR040000	☐	w/o Sideguard		
FOOTREST AND FOOTPLATES				
KRR050029	☒	Footrest, tubular, lightweight, aluminium	(L)	
KRR050026	☐	Platform footplate lightweight angle & depth adjustable carbon		395
KRR050125	☐	Performance footplate angle adjustable	(3)	285
KRR050030	☐	Footrest, tubular, lightweight, aluminium, with cover		95
KRR050025	☐	Platform footplate lightweight angle & depth adjustable composite		195
KRR050024	☐	Platform footplate lightweight angle & depth adjustable aluminium (colour: black, matt)		195
KRR050128	☐	Footboard, platform, plastic, angle and depth adjustable, flip-up (to the back)	(4)	210
KRR050127	☐	Adjustable side protection, to create a side-stop on the footboard	(5)	30
LOWER LEG LENGTH (LLL) Attention: LLL has to be 30 mm less than seat height front; all measurements in mm				
☐ 360 ☐ 370 ☐ 380 ☐ 390 ☐ 400 ☐ 410 ☐ 420 ☐ 430 ☐ 440 ☐ 450 ☐ 460 ☐ 470 ☐ 480 ☐ 490 ☐ 500 ☐ 510				
FRAME VERSIONS				
Frame length				
KRR010206	☒	Frame standard		
KRR010207	☐	Frame long (20 mm additional frame length)	(not in combination with SD 480 mm)	
Frame angle				
KRR010012	☒	100° frame (UK 80°)		
KRR010013	☐	92° Performance frame (UK 88°)		
Frame inset				
KRR010101	☒	Inset: Small (10 mm)		
KRR010103	☐	Inset: Performance (30 mm)	(not available for seat width 320 mm)	
Castor position (see explanation on back pages)				
KRR010301	☒	Castor house position wide		
KRR010302	☐	Castor house position narrow		
(1) With 24" rear wheel min. rear seat height 410 mm; with 25" rear wheel min. rear seat height 430 mm (2) Only with COG 56 mm (3) Not with SW 320 mm and 30 mm inset (4) User weight limitation 75 kg, not with 4" castor (5) Only for platform footplates lightweight angle & depth adjustable (L)= lightest option (+)= adds 500 g or more to the wheelchair-weight				

CENTRE OF GRAVITY (COG) - all measurements are in mm; measured from front of back post; choose only one option

COG56	☐	56 mm
COG68	☐	68 mm
COG80	☐	80 mm
COG92	☐	92 mm
COG104	☐	104 mm

**FRAME ACCESSORIES**

KRR090031	☐	Suspension up to 70 kg	(1) (+)	(chair will be built using aluminium axle-stems and axle-tube)	585
KRR090032	☐	Suspension up to 120 kg	(1) (+)	(chair will be built using aluminium axle-stems and axle-tube)	585

Anti tip flip up (user operated) (2)

KRR090033	☐	Left	KRR090008	☐	Right	(+)	(Please choose only one)	175
-----------	-----------------------	------	-----------	-----------------------	-------	-----	--------------------------	-----

Tip assist (step tube) (2)

KRR090001	☐	Left	KRR090002	☐	Right			35
KRR090003	☐	Pair				(+)	(not available with anti-tip)	70
KRR090009	☐	Crutch holder					(only in combination with tip assist)	80

Anti-tip tubes swing-away (2)

KRR090004	☐	Left	KRR090005	☐	Right	(+)		125
KRR090006	☐	Pair				(+)	(not available with tip-assist)	195
KRR090007	☐	Quickie style anti-tip tubes pluggable	(2) (+)				(Pair)	95

KRF090016	●	Frame protection, black					(textile, with velcro)	
KRR090577	☐	Frame protection, Plastic, glued to the frame, black						85
KRR090010	☐	Transit wheels	(2) (+)					235

FORKS

KRR080013	☐	Castor fork, Carbotecture®, 111mm	(L)								
Colour:	☒	Black	☐	Silver	☐	Orange	☐	Red (3)	☐	Blue (3)	
KRR080029	☐	Castor fork, Aluminium, 123mm									
Colour:	☒	Black	☐	Silver	☐	Orange					
KRR080008	☐	Castor fork, aluminium, one-arm								(only with 4" & 5" Alu soft, Alu solid)	215
Colour:	☒	Black	☐	Silver	☐	Orange	☐	Red	☐	Blue	
KRR080007	☐	Castor fork, aluminium, suspension (Frog Legs®), black	(+)								365

CASTORS

KRR080306	☐	Castor, solid (soft), rim: aluminium								(only 4" and 5")	
Colour, rim:	●	Silver	☐	Black	☐	Orange	☐	Red	☐	Blue	
KRR080304	☐	Castor, solid (hard), rim: aluminium	Colour, rim: Silver							(only 4" and 5")	
KRR080301	☐	Castor, solid (hard), rim: composite, black								(3"-6", not with One arm fork)	
KRR080302	☐	Castor, solid (soft), rim: composite, black								(only 5" and 6", not with One arm fork)	

REAR WHEELS

KRR070003	☐	Rear wheel, Lightweight, 24 spokes (silver), straight spokes								(Std.: black hub/rim; colour options last page)	
KRR070004	☐	Rear wheel, Proton, 24 spokes (black), straight spokes	(L)							(user weight limit 100kg)	250
KRR070005	☐	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver									495
KRR070007	☐	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black									595
KRR070000	☐	Without rear wheel									

(1):	Not possible with carbon axle-stems and axle: chair will be built with aluminium axle-stems and axle-tube
(2):	Not with KRR07049: camber adapter 0°, superlight.
(3):	Red and blue only in combination with soft castor
	(L)= lightest option (+)= adds 500 g or more to the wheelchair-weight

			£
TYRES			(choose only one tyre type)
KRR070107	☐ Tyre, pneumatic, slick, Schwalbe One	(L)	115
KRR070102	☐ Tyre, puncture proof, tread	(+)	
KRR070104	☐ Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	(+)	95
KRR070101	☐ Tyre, pneumatic, fine profile, Schwalbe Right Run	(not in 25")	0
HANDRIMS			
KRR070201	☐ Handrims, aluminium, powder-coated, clear	(L)	
KRR070202	☐ Handrims, aluminium, anodised, silver	(L)	
KRR070203	☐ Handrims, aluminium, powder-coated, coloured	(L)	(select colour on back pages, charged once)
KRR070212	☐ Handrims, Ellipse 3R, black	(+)	(ergonomic handrims with rubber-strip)
KRR070204	☐ Handrims, stainless steel, silver	(+)	
KRR070205	☐ Handrims, titanium, untreated		240
KRR070206	☐ Handrims, titanium, powder-coated, clear		240
KRR070207	☐ Handrims, Supergrip ®	(+)	250
KRR070208	☐ Handrims, Maxgrepp ®	(+)	(only 24") 350
KRR070250	☐ Handrim cover, rubber	(not with 25")	110
Hand rim mounting			
KRR070317	☒ Handrim mount, narrow (17 mm)		
KRR070307	☐ Handrim mount, super narrow (7 mm)	(only with stainless steel)	
KRR070327	☐ Handrim mount, wide (27 mm)		
WHEEL CAMBER			(choose only one option)
KRR070049	☐ 0° superlight, no width-adjustment	(1) (L)	
KRR070600	☐ 0° width-adjustment +/-15mm		
KRR070603	☐ 3° width-adjustment +/-15mm		
WHEEL LOCKS (Choose only one option)			
KRR060003	☐ Compact wheel lock	(2) (L)	
KRR060001	☐ Standard wheel lock		
KRR060024	☐ Lightweight compact brake EVO		200
KRR060004	☐ Lightweight Compact wheel lock		185
KRR060002	☐ Knee lever wheel lock		95
KRR060010	☐ Extension for knee lever wheel lock, slip-on, swing-away	(only with KRR060002)	50
KRR060006	☐ One arm wheel lock (left hand actuation)	(+)	200
KRR060007	☐ One arm wheel lock (right hand actuation)	(+)	200
ACCESSORIES			
KRR090000	☒ Tool kit		
KRR090024	☒ High pressure pump		90
KRR090025	☐ Air pump		20
KRR090026	☐ Mobile phone pocket		55
KRR090030	☐ Back pack		75
(1): Camber adapter glued into the axle-tube. No adjustment of the track-width possible.			
(2): Adjustment of COG ± 2 cm on a chair in the field.			
(L)= lightest option (+)= adds 500 g or more to the wheelchair-weight			

				£
COLOURS FOR HANDRIMS				
	Handrim (1) £75	Hub/Rim £200		
01 Matte Black	☐	☐	12 Royal Blue	☐
02 Gloss Black	☐	☐	13 Brilliant Ocean	☐
03 Metallic Shadow	☐		14 Scarlet Noir	☐
04 Starlight White	☐		15 Red Pepper	☐
05 Silver Ore	☐		16 Cherry Lustre	☐
06 Storm Gray	☐		17 Turbo Green	☐
07 Desert Tan	☐		18 Pink Power	☐
08 Desert Sage	☐		19 Satin Rose	☐
09 Salute Green	☐		20 Mystic Purple	☐
10 Velvet Sapphire	☐		21 Solar Orange	☐
11 Midnight Blue	☐		22 Velocity Yellow	☐
FRAME-COLOURS (Not including axle or backrest)				
KRR100139	☐ Natural Carbon, no coating			
KRR100137	☐ Transparent, gloss			
KRR100138	☐ Transparent, matt			
KRR101101	☐ Matte Black			
KRR101102	☐ Gloss Black			
FRAME-GRAPHICS "Flying Pixel" (scratch resistant, non-removeable, includes the "Krypton" sticker)				
KRR100590	☐ Frame graphics Flying Pixel, silver			
KRR100591	☐ Frame graphics Flying Pixel, red			
KRR100592	☐ Frame graphics Flying Pixel, blue			
COLOUR STICKERS "QUICKIE" and "Krypton" (3)				
KRR100550	☐ Black			
KRR100551	☐ Silver			
KRR100552	☐ Orange			
KRR100553	☐ Red			
KRR100554	☐ Blue			
COLOURS FOR LIGHTWEIGHT WHEELS 24" (4)				
KRR100336	☐ Rear wheel, lightweight anodised, orange (4)			
KRR100337	☐ Rear wheel, lightweight anodised, blue (4)			
KRR100338	☐ Rear wheel, lightweight anodised, red (4)			
Comments for special build:				
				
Krypton R - Base Price				£4,250
Option price				
Total Price				
<i>(1) only with KRR070203, handrim aluminium, powder coated</i> <i>(3) No choice possible when frame graphics "Flying Pixel" are chosen.</i> <i>(4) Anodised colours for hub and rim, black spokes. Only with light weight wheel in 24"</i>				

WEIGHT

We want to be very transparent on weight. The following information will be a guide to the lightest possible Krypton R: how to achieve a 6.2 kg wheelchair, and how to order the lightest possible configuration for your Krypton R

Krypton R: 6.2 kg configuration Transport weight: 3.7 kg

Frame: high frame, 92°, 30 mm inset, seat width: 360 mm, seat depth: 400 mm. Seat height: front 470 mm, rear 420 mm. Seat-sling Airflow superlight.

Wheel-lock : lightweight compact wheel-lock

Footrest: tubular footrest



Backrest: Carbon backrest, EXO superlight backrest-upholstery, without push-handles

Sideguards: Carbon fibre, slim-style

Rear wheels: 24" Proton wheels with Aluminium handrims (coated or anodized) and Schwalbe One tyres

Forks and castors: 111 mm carbotechnology fork, 3" castor

Please mind that there is a 5% tolerance on all weight specifications due to weight deviations on mounting parts (tyres, castor-wheels etc.)

LIGHTWEIGHT-GUIDE FOR ORDERING THE KRYPTON R

Sunrise Medical offers a wide variety of options for the Krypton R to suit every user's needs. However, most options add significant weight to the wheelchair. To make this more transparent when filling in the order-form, and to assist you in ordering the lightest possible wheelchair, we indicate the approximate weight-increase of the options on the order form:

- (L)** indicates the lightest possible option within the range of options
- (+)** indicates that this option adds more than 500 g to the wheelchair's weight compared to the lightest option within the range

LIGHTWEIGHT OPTIONS FOR THE KRYPTON R

We developed several options for the Krypton R to save weight on the chair without sacrificing daily-use suitability, comfort or durability



Carbon backrest (KRR030006): full carbon backrest with very close cross-tube (20mm cross-tube to backrest-upholstery) (-100 g compared to Standard backrest)



Seat sling "airflow" superlight (KRR020019): fully adjustable seat-sling with composite buckles and lighter straps (max. user weight 85 kg) (-100 g compared to seat sling "airflow")



EXO superlight (KRR030005): lightweight version of the EXO backrest upholstery. EVA-foam padding, adjustable straps, wider cut of the upholstery to cover backrest tubes. (-160 g compared to Std. EXO backrest-upholstery)



Axle-plug superlight (KRR070049): (only 0° camber): glued into the axle-tube, no adjustability of the track-width possible (-160 g compared to width-adjustable axle-plug KRR070600 / 603)

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) *For the user:* Allow space between sideguard and hip on each side of the chair
(approx. 10 - 20 mm altogether)
- b) *On the wheelchair:* From outside frame tube to outside frame tube.

Seat depth (SD):

- a) *For the user:* When sitting upright: leave space for 2 - 3 fingers between
front edge of the seat upholstery and back of knee
(approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) *On the wheelchair:* From back tube to front of seat upholstery.

Front seat height:

- a) *For the user:* When sitting with 90° knee angle:
From back of knee to bottom of shoe sole.
- b) *On the wheelchair:* From front of seat tube to the ground.
(min. 30 mm longer than user measurement)
Remember to include the cushion height.

Rear seat height:

- b) *On the wheelchair:* From back of seat tube to the ground.
Note: Max. seat height difference rear to front is 130 mm

Back height:


- a) *For the user:* From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) *On the wheelchair:* From back of seat tube to upper edge of backrest upholstery.
Note: The user's ability to sit upright has a major influence on the back height
and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability.
0° is the angle 90° to the seat. Negative angles are <90°
and positive angles are >90°.

Low Seat Frame vs. High Seat Frame:

The high frame is 40 mm higher than the low frame. So you can choose shorter and lighter forks to minimise the weight of your wheelchair.



Technical Data Overview:

Max. overall height:	950 mm	Max. overall width:	Camber	24" rear wheel	25" rear wheel
Max. overall length:	930 mm		For 0° camber:	SW + 170 mm	SW + 170 mm
Max. turning circle:	700 mm		For 3° camber:	SW + 230 mm	SW + 240 mm
Max. safe slope:	7°				
Wheelchair weight:	from 6.2 kg				
Dismantled weight:	from 3.7 kg				
Weight of heaviest part:	2.1 kg				

All changes in the configuration might change the seat heights.
Manufacturing tolerances are +/- 10 mm for all measurements.

© Sunrise Medical 2017

SEAT UPHOLSTERIES AND CUSHIONS



KRR020012 Seat sling, airflow



KRR020019 Seat sling, airflow, superlight

BACKREST



*KRR030006
Carbon Backrest*



*KRR030102
Backrest, angle-adjustable, lock when folded*



*KRR030201
Push handles long*



*KRR030203
Push handles, fold-down*



*KRR030204
Push handles, height-adjustable*



*KRR030316
EXO backrest sling*




*KRR030317
EXO PRO backrest sling*



*KRR030005
EXO superlight backrest sling*

FOOTREST AND FOOTPLATES



*KRR050029
Tubular footrest lightweight*



*KRR050025
Platform footrest LW composite plate*



*KRR050026
Platform footrest LW carbon plate*



*KRR050125
Performance footrest alu angle adj*



*KRR050127
Foot positioning plate*

SIDEGUARD & ARMREST



KRR040101 Slim style sideguard



KRR040104 Carbon fibre, cold weather protection, fender



KRR040107 Sideguard, composite, removable (black)




KRR040117 Tubular armrest, swing away,



KRR040052 Single-post sideguard, armpad (33 cm)


FRAME VERSIONS: CASTOR HOUSE POSITIONS

Position wide:



Further to the front and wider position of the castors for optimised rolling resistance

Position narrow:



20 mm further rearwards and 20 mm per side more narrow position of the castors for a more narrow front-end and improved manoeuvrability in tight spaces



KRR010301 Castor house position wide



KRR010302 Castor house position narrow



Castor house position wide vs. narrow

FRAME ACCESSORIES



KRR090032
Suspension



KRR090001
Tip assist



KRR090004
Anti-tip tubes SA



KRR090007
Quickie style anti tip tubes pluggable



KRR090033
Anti tip flip up (user operated)



KRR090010
Transit Wheels



KRR090010
Frame protection, Plastic, glued to the frame, black

FORKS



KRR080013 & KRR080306
Castor fork, Carbotechnology & Castor solid (soft) aluminum



KRR080008
Castor fork, aluminum, one-arm, red

REAR WHEELS & HANDRIMS



KRR070001



KRR070002



KRR070003
Lightweight wheel



KRR070004
Proton wheel



KRR070211
Handrims Surge®



KRR070208
Handrims Maxgrepp

TYRES



KRR070101
Schwalbe Right Run



KRR070107
Schwalbe One



KRR070104
Schwalbe Marathon Plus



KRR070102
Tyre, puncture proof

WHEEL LOCKS



KRR060001
Standard brake



KRR060002
Knee lever brake



KRR060003
Compact brake



KRR060004
Lightweight Compact brake

FRAME COLOURS



1. *Transparent, gloss*
2. *Transparent, matt*
3. *Black, gloss*
4. *Black, matt*

FRAME GRAPHIC "FLYING PIXELS"



From left to right:

- *Frame graphics Flying Pixel, red*
- *Frame graphics Flying Pixel, silver*
- *Frame graphics Flying Pixel, blue*

24 " LIGHT WEIGHT WHEEL COLOURS



*KRR100336
Lightweight wheel, anodized orange*



*KRR100337
Lightweight wheel, anodized blue*



*KRR100338
Lightweight wheel, anodized red*