

M5 Corpus®

More adventure

The M5 Corpus is made for more adventure, which gives it a capability all of its own for users who require more speed, stronger motors, higher elevation, better reach and exceptional comfort, stability and traction in varying driving situations. With powerful LED lights and increased driving range and battery life, the colourful new M5 Corpus is designed to give you much more.

Active Reach

Tilts the seat forward by up to 20° and raises it to extend functional reach, aid in transfers and your daily activities.

Active Height

Enables you to drive at up to 5 km/h while fully elevated and combines seamlessly with *Active Reach* to give you the ultimate advantage in accessing your environment.

Corpus: the ultimate seat system

Improved comfort, pelvic positioning, lateral stability and immersion are all key objectives of the new Corpus cushions.

Corpus cushions feature softer dual-density foam on the backrest, a re-designed pelvic well and Stretch-Air™, a breathable, water-resistant cushion fabric that is machine-washable at 60°C.

More control

Agile independent suspension, works on all 6 wheels, offering superior ride comfort, traction and stability to give you more control when you are out exploring.

More drivability

The new Permobil battery and Permobil VoltPro charger work together to provide you with up to 25% greater driving distance, optimized charging times and longer battery life.







The softer dual-density backrest cushion and re-designed pelvic well increase immersion into the seat, giving more comfort *and* stability

MORE SUPPORT

BodiLink positioning accessories specially designed for the Corpus seating system (option)

MORE COLOUR

A fresh new aesthetic allows you to mix & match shroud/hub insert colour options to match your style.

MORE VISIBILITY

LED lights front and back are twice as bright for improved visibility.

MORE FREEDOM

to move in a greater variety of drive positions with enhanced stability.

MORE CARRYING CAPACITY
Optional luggage/essentials carrier

ARCTIC





M5 Corpus® Specifications

Max. speed	
Max. Driving Range (G24: 85 Ah batteries)	
Max. user weight	
Obstacle climbing (with run up)	
Min. turning diameter	
Pivot width	
Electronics R-Net	
Crash test ISO 7176-19 – 4-point tie-down	
Crash test ISO 7176-19 - Dahl docking station	
Tested according to EN12182/EN12184 CE	
Max. Length incl. seat	
Base width	
Width incl. seat and armrest	
Height	
Stowage length/ height	880 mm/870 mm
Weight incl. batteries 85 Ah	
Seat to floor height (excl. cushion)	
Seat elevator height (excl. cushion)	450-800 mm
Seat depth	
Seat width	420–570 mm (50 mm increments)
Backrest width	
Backrest height (excl. seat cushion)	470; 545–670 mm (25 mm increments
Armrest height	185–320 mm
Distance between armrests	
Backrest angle adjustment (manual)	
Seat plane adjustment – Anterior/posterior tilt	
*Range may vary depending on user weight and if the wheelchair is used frequently on slopes, rough ground or to climb kerbs often.	
5	



















