

Karma Mobility Ltd

Email: info@karmamobility.co.uk

Website: www.karmamobility.co.uk

Tel 0845 630 3436

Fax 0845 630 3736

Dealer Prescription Form

VIP\_16081\_en\_rev10

User Information	Dealer Information
Order No.	Order No.
Order Date:	Order Date:
Client Name:	Client Name:
Phone:	Phone:



### Standard Features

- 115 kg Maximum Weight Capacity
- Self Propulsion with 22" Q-Release Rear Wheels
- Transit Propulsion with 14" Q-Release Rear Wheel
- 7" Caster Wheels
- 90° to 157° Backrest Angle
- Patented V Shaped Seating System
- Seat Depths: 17", 18"
- Detachable Headrest with Neck Pillow
- Elevating Legrests
- Detachable Seat and Back Cushion
- Slowing Brakes

◆ Standard feature    ○ Free of Charge

Seat																																											
<b>Seat Width and Seat Depth</b>																																											
<table border="1"> <thead> <tr> <th>Depth \ Width</th> <th>380mm(16")</th> <th>405mm(17")</th> <th>430mm(18")</th> <th>455mm(19")</th> <th>480mm(20")</th> <th>depth</th> </tr> </thead> <tbody> <tr> <td>405mm(15")</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td>16" STD</td> </tr> <tr> <td>405mm(16")</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td>16" STD</td> </tr> <tr> <td>405mm(17")</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td></td> <td></td> <td>16" STD</td> </tr> <tr> <td>405mm(18")</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td>18" STD</td> </tr> <tr> <td>405mm(20")</td> <td></td> <td></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td>18" STD</td> </tr> </tbody> </table>	Depth \ Width	380mm(16")	405mm(17")	430mm(18")	455mm(19")	480mm(20")	depth	405mm(15")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			16" STD	405mm(16")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			16" STD	405mm(17")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			16" STD	405mm(18")			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	18" STD	405mm(20")			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	18" STD	
Depth \ Width	380mm(16")	405mm(17")	430mm(18")	455mm(19")	480mm(20")	depth																																					
405mm(15")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			16" STD																																					
405mm(16")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			16" STD																																					
405mm(17")	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			16" STD																																					
405mm(18")			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	18" STD																																					
405mm(20")			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	18" STD																																					
<b>Seat Width</b>																																											
<input type="checkbox"/> Padded Upholstery	<input type="radio"/>																																										
<b>Seat Depth</b>																																											
<input type="checkbox"/> STD (Standard 15" to 17" = 16" & 18", 20" = 18" deep)	<input type="radio"/>																																										
<input type="checkbox"/> +1"	£10.00																																										
<input type="checkbox"/> +2"	£10.00																																										
<input type="checkbox"/> Footrest Hanger Bracker L/H	£40.00																																										
<input type="checkbox"/> Footrest Hanger Bracker R/H	£40.00																																										
<b>Backrest</b>																																											
<b>Backrest Height</b>																																											
<input type="checkbox"/> Standard Backrest Upholstery	<input type="radio"/>																																										
<b>Backrest Upholstery</b>																																											
<input type="checkbox"/> Strap Tension Adjustable Backrest (Including Detachable Padded Cover)	<input type="radio"/>																																										
<b>Legrests/Footplates</b>																																											
<b>Legrest RH</b>																																											
<input type="checkbox"/> Elevating Legrest (0° to 80°) RH Swing in/Swing out	<input type="radio"/>																																										
<input type="checkbox"/> 80° RH Swing out Legrest	◆																																										

<input type="checkbox"/> Stump Support (If instead of ELR/STD-->No upgrade charge, just swap to stump suport)	£40.00
---	--------

<b>Legrests/Footplates (cont)</b>	
<b>Legrest LH</b>	
<input type="checkbox"/> Elevating Legrest (0° to 80°) RH Swing in/Swing out	○
<input type="checkbox"/> 80° RH Swing out Legrest	◆
<input type="checkbox"/> Stump Support	£40.00
<b>Footplate RH</b>	
<input type="checkbox"/> Flip-up Tube-in-Centre Footplate	○
<input type="checkbox"/> Flip-up LH Angle Adjustable Footplate (Not available for 14" seat width)	£20.00
<input type="checkbox"/> If Standard Footplates Not Required Deduct	-£11.00
<input type="checkbox"/> Flip-up Tube-in-Centre Footplate for ELR	◆
<input type="checkbox"/> Angle Adjustable Footplate for ELR	£20.00
<b>Footplate LH</b>	
<input type="checkbox"/> Flip-up Tube-in-Centre Footplate	○
<input type="checkbox"/> Flip-up LH Angle Adjustable Footplate (Not available for 14" seat width)	£20.00
<input type="checkbox"/> If Standard Footplates Not Required Deduct	-£11.00
<input type="checkbox"/> Flip-up Tube-in-Centre Footplate for ELR	◆
<input type="checkbox"/> Angle Adjustable Footplate for ELR	£20.00
<b>Armrests</b>	
<b>Armrest RH</b>	
<input type="checkbox"/> Flip-back PU Armpad RH 220~240mm Height	○
<b>Armrest LH</b>	
<input type="checkbox"/> Flip-back PU Armpad LH 220~240mm Height	○
<b>Front Castors/Forks</b>	
<b>Castors</b>	
<input type="checkbox"/> 175mm x 30mm Solid Tire (7")	○
<input type="checkbox"/> 175mm x 45mm solid Tyre (7" x 1.75")	£24.00
<input type="checkbox"/> 175mm x 45mm Pneumatic Tyre (7" x 1.75")	£35.00
<input type="checkbox"/> If Standard Castor not Required Deduct	-£12.00
<b>Forks</b>	
<input type="checkbox"/> Standard Fork	○
<input type="checkbox"/> 2" Wide Fork	£12.00
<input type="checkbox"/> If Standard Castor not Required Deduct	-£10.00
<b>Rear Wheels</b>	
<b>Rear Wheel</b>	
<input type="checkbox"/> 14" Solid Tire Plastic Rim	○
<input type="checkbox"/> 14" Pneumatic Tire Plastic Rim	◆
<input type="checkbox"/> 22" Solid Tire Spoked Rim with Ergo Handrim	○
<input type="checkbox"/> 22" Pneumatic Tire Spoked Rim with Ergo Handrim	◆
<b>Brakes</b>	
<input type="checkbox"/> Pull to Lock Brakes (Only available for 20", 22" and 24" rear wheel)	○
<input type="checkbox"/> Attendant Brakes	○
<b>Safety Options</b>	
<input type="checkbox"/> Anti-tips - pair	○
<b>Push handles</b>	
<input type="checkbox"/> Without headrest brackets - for 3rd party backrests	◆

<b>Accessories</b>	
<input type="checkbox"/> Lap Belt (Safety Belt)	£6.40
<input type="checkbox"/> Heel loops - pair	£6.60
<input type="checkbox"/> Lap Tray	£25.00
<input type="checkbox"/> Cane Holder	£37.50
<input type="checkbox"/> Footboard Standard	£25.00
<input type="checkbox"/> Footboard Deep	£25.00
<input type="checkbox"/> Foldable Headrest	£65.00
<b>Price</b>	
	<b>Total price</b>