

Permobil

Progeo Joker



Aluminum



Rigid



Adjustable

permobil

Fully personalized frame and easy adjustments

The Joker range combines elegant rigid frame design with practical detailing. Form and function are harnessed to meet your mobility needs without compromising your sense of style.

For active manual wheelchair users who may need to make adjustments as they gain driving experience or progress on their rehabilitation journey.



Progeo Joker

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37.5 - 40 - 42.5 - 45 - 47.5 - 50 cm
Seat height	Front from 40 to 59 cm (adjustable) Rear from 38 to 46 cm (adjustable)
Footsupport distance	From 35 to 54 cm (adjustable)
Backrest height	From 24.5 to 47 cm (adjustable)
Front frame angle	85° - 90° - 100°
Total length	88 cm approx.
Total width	Seat width + 17 cm (camber 0°)
Backrest angle	From 74° to 94° (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 5.9 kg approx.
Wheelchair weight	Starting from 8.5 kg approx.
Maximum load	125 kg

BLACK

BRIGHT LIME

BRIGHT FUCHSIA

BRIGHT BLUE

BRIGHT RED

Personalise your ride*

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Actual product colors may vary.