

permobil Order Form - Progeo noir 2.0

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	UK (English)	
User	Company stamp	Rigid frame wheelchair
End user weight (kg)		£ 8071,00
Language user manual	Responsible	Retail price - VAT not included

Carbon fibre frame made-to-measure.

NOIR is a wheelchair without possibilities of adjustments after its manufacturing.

See reference tables

Take care to write all the requested measures (in centimetres)

Frame Angle

85° 90° 95°

Seat depth (PS)

35 37,5
 40 42,5
 45 47,5

See note (1)

Backrest rake adjustable from -4° to +16°

Seat width (LS)

33 36
 39 42
 45

Front height (min 44-max 54)

Rear height (min 37-max 48)
calculated with rear wheels 24"

Setting C

Extremely prudential 7,5
 Prudential 9
 Standard 10,5
 Active 12
 Extreme 13,5

camber

2° 0° 4°

Footplate option

Carbon fibre frame with integrated footplate, not adjustable
(Lightest configuration)

Aluminium footplate adjustable distance +/- 1,5cm (rake and depth adjustable)
Minimum footplate distance 35cm.

Internal 2/3 Internal 2/3 External External

11531 **no list**

Seat

Seat cloth with pocket (under seat)

Seat with tension straps 08883 **224,00**

Carbon rigid seat 08336 **475,00**

Pocket (under seat) [25x11x4,5] 08281 **74,00**

Pag. 1 / 6

UK_202304

Rev. 2023.2_Mar


customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Date Order User Address

Since that the paint is sensitive to the surface material and texture, the colour selected for the wheelchair frame can also differ considerably from what is displayed on a RAL colour chart or also compared to what is displayed on a monitor.

CARBON

☆ Polished transparent painting




Opaque transparent painting 08925 624,00

ESSENCE

Polished transparent painting

<input type="radio"/> ANTHRACITE Essence	10351	876,00
<input type="radio"/> BLUE Essence	08928	876,00
<input type="radio"/> GREEN Essence	11163	876,00
<input type="radio"/> LIME Essence	08926	876,00
<input type="radio"/> ORANGE Essence	10349	876,00
<input type="radio"/> RED Essence	08927	876,00
<input type="radio"/> SAND Essence	11165	876,00
<input type="radio"/> VIOLET Essence	11167	876,00
<input type="radio"/> WHITE Essence	10350	876,00



Opaque transparent painting

<input type="radio"/> ANTHRACITE Opaque Essence	08929	1.121,00
<input type="radio"/> BLUE Opaque Essence	10354	1.121,00
<input type="radio"/> GREEN Opaque Essence	11162	1.121,00
<input type="radio"/> LIME Opaque Essence	10352	1.121,00
<input type="radio"/> ORANGE Opaque Essence	10355	1.121,00
<input type="radio"/> RED Opaque Essence	10353	1.121,00
<input type="radio"/> SAND Opaque Essence	11164	1.121,00
<input type="radio"/> VIOLET Opaque Essence	11166	1.121,00
<input type="radio"/> WHITE Opaque Essence	10356	1.121,00

my ESSENCE colour RAL _____

Glossy Opaque

Note

11186 1.293,00

RACE

Polished transparent painting

<input type="radio"/> ANTHRACITE Race	11168	876,00
<input type="radio"/> BLUE Race	11170	876,00
<input type="radio"/> GREEN Race	11172	876,00
<input type="radio"/> LIME Race	11174	876,00
<input type="radio"/> ORANGE Race	11176	876,00
<input type="radio"/> RED Race	11178	876,00
<input type="radio"/> SAND Race	11180	876,00
<input type="radio"/> VIOLET Race	11182	876,00
<input type="radio"/> WHITE Race	11184	876,00



Opaque transparent painting

<input type="radio"/> ANTHRACITE Opaque Race	11169	1.121,00
<input type="radio"/> BLUE Opaque Race	11171	1.121,00
<input type="radio"/> GREEN Opaque Race	11173	1.121,00
<input type="radio"/> LIME Opaque Race	11175	1.121,00
<input type="radio"/> ORANGE Opaque Race	11177	1.121,00
<input type="radio"/> RED Opaque Race	11179	1.121,00
<input type="radio"/> SAND Opaque Race	11181	1.121,00
<input type="radio"/> VIOLET Opaque Race	11183	1.121,00
<input type="radio"/> WHITE Opaque Race	11185	1.121,00

my RACE colour RAL _____

Glossy Opaque

Note

11187 1.293,00

FLASH

Polished transparent painting

<input type="radio"/> BLUE Flash	11190	876,00
<input type="radio"/> WHITE Flash	11192	876,00
<input type="radio"/> RED Flash	11188	876,00



Polished transparent painting

<input type="radio"/> BLUE Opaque Flash	11191	1.121,00
<input type="radio"/> WHITE Opaque Flash	11193	1.121,00
<input type="radio"/> RED Opaque Flash	11189	1.121,00

Fork colour

- ☆ Opaque black
- Bright anthracite
- Bright orange
- Bright silver
- Bright blue
- Pearl white
- Bright fuchsia
- Bright lime
- Bright red

Date	Order	User	Address
------------	-------------	------------	---------------

Backrest adjustment and fixing system

Backrest with double locking (upright and folded) with rake adjustable from **-4° to +16°**
 Backrest with upright locking system ONLY (no locking when folded) with rake adjustable from **-4° to +16°**
 Fixed back(non-folding back) with rake adjustable from **-4° to +16°**

Integrated backrest (Banded, not folding back) 10841 **1.248,00**
 Exclude installation of the titanium backrest

☆ 0° (90° to the floor)
 (.....) to the ground to the ground

Backrest

Adjustable velcro straps and padded transpiring upholstery ☆ Carbon fibre backrest
 Titanium backrest
 Not available with: Integrated backrest (Banded, not folding back)

AIR Integrated posture backrest (NOT REMOVABLE) 1.047,00
 Carbon shell - Lateral retaining 7,5cm
 Depth adjustment ±1cm. - Angle adjustment ±10°
 not available for seat width 48

H 28 Backrest Height (TS) min 34 (.....)	H 34 Backrest Height (TS) min 38 (.....)	H 40 Backrest Height (TS) min 44 (.....)			
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Installation on:
 Backrest tubes
 (Carbon tubes with integrated backrest)

Carbon fibre backrest
 Titanium backrest
 Not available with: Integrated backrest (Banded, not folding back)

Handle tube

Tubes with short handle
 SPORT tubes Without push handles (Lightest configuration)
 Folding handles

02958 **252,00**







Casters

3" Sport (Lightest configuration)
 4" Soft
 3" Active Soft (wide) (Aluminium alloy rim)
 4" Active Sport (narrow) (Aluminium alloy rim)
 4" Active Soft (wide) (Aluminium alloy rim)
 5" Active Soft (wide) (Aluminium alloy rim)
 5" Pneumatic

Rear wheels - See chart N°3

24" ULTRA-LIGHT PROGEO KEY SPORT (Lightest configuration)
 24" Ultra-light "Spinergy LX"
 24" Ultra-light "Spinergy LXL"
 24" Ultra-light "Spinergy Flex rim" Black hub and white spokes Hub and spokes black
 25" Ultra-light "Spinergy LX" **Integrated aluminium/rubber handrim**
 25" Ultra-light "Spinergy LXL" Black hub and white spokes
 26" Ultra-light "Spinergy LX"

04353 **527,00**


Date	Order	User	Address
Rear wheel tyre			
<input checked="" type="checkbox"/> High pressure tyre	(24"-25"-26")	(Lightest configuration)	
<input type="checkbox"/> High pressure profile tyre	(24")		
<input type="checkbox"/> Marathon Plus	(24"-25")	with anti-puncture protection	01669 63,00
<input type="checkbox"/> Anti-puncture - Colour Black	(24")		11083 126,00
Handrim			
	Wheel distance <input checked="" type="checkbox"/> 3 cm	<input type="checkbox"/> 4 cm	
<input checked="" type="checkbox"/> Aluminium	(24"-25"-26")	(Lightest configuration)	
<input type="checkbox"/> Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170 126,00
<input type="checkbox"/> Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684 148,00
<input type="checkbox"/> Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016 275,00
<input type="checkbox"/> Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017 315,00
<input type="checkbox"/> Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834 349,00
<input type="checkbox"/> Ergo Para	(24")		00632 349,00
<input type="checkbox"/> TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023 399,00
<input type="checkbox"/> Oval countour with a rubber insert			
<input type="checkbox"/> NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025 263,00
<input type="checkbox"/> Ergonomic black profile			
<input type="checkbox"/> NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026 286,00
<input type="checkbox"/> Ergonomic profile with rubber coating (ideal for quadriplegics)			
<input type="checkbox"/> SURGE LT	(25")	Oval countour with a rubber insert	06427 360,00
<input type="checkbox"/> Quad	(24")		00693 200,00
<input type="checkbox"/> Silicon cover for standard aluminium handrim		Available for measures (22" - 24" - 26")	04050 160,00
Wheel locks			
<input checked="" type="checkbox"/> Sport brake			
<input type="checkbox"/> Bent lever push to lock			
<input type="checkbox"/> Straight lever push to lock			
<input type="checkbox"/> Straight lever pull to lock			
<input type="checkbox"/> Brake with extension lever push to lock			
<input type="checkbox"/> Brake with extension lever pull to lock			
Mudguard			
<input checked="" type="checkbox"/> Carbon fibre mudguard			
<input type="checkbox"/> Carbon fibre Straight sideguard	<input checked="" type="checkbox"/> Aligned wheel	<input type="checkbox"/> Under wheel cm	<input type="checkbox"/> Above wheel cm
<input type="checkbox"/> Removable carbon fibre mudguard			04109 252,00
<input type="checkbox"/> Removable Carbon fibre straight sideguard	<input checked="" type="checkbox"/> Aligned wheel	<input type="checkbox"/> Under wheel cm	<input type="checkbox"/> Above wheel cm
			04109 252,00
Armrests			
<input type="checkbox"/> Tubular armrests	(Swing away, removable and height adjustable armrest)		 00694 252,00
Accessories			
<input type="checkbox"/> Pad for front frame protectors			07058 52,00
<input type="checkbox"/> Backrest-seat connection cloth			10241 17,00
<input type="checkbox"/> Removable storage bag (with quick coupling)			09811 126,00
Optional 1			
<input type="checkbox"/> Transit wheels			 03290 148,00
<input type="checkbox"/> Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance)		
	<input type="checkbox"/> Right		06821 126,00
	<input type="checkbox"/> Left		06822 126,00
Optional 2			
<input type="checkbox"/> Safety pelvic belt			03972 63,00
<input type="checkbox"/> Harness			04103 126,00
<input type="checkbox"/> Lumbar support	(to insert in the backrest upholstery)		00013 40,00
<input type="checkbox"/> Tetraclip	(Quick release axle for quadriplegic)		03947 177,00
Additional Kit			
<input type="checkbox"/> 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)		03378 498,00
<input type="checkbox"/> 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)		09357 876,00

Date Order User Address


PHYSIO - Posture backrest (Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT	28	34	40	1.047,00		
Shell's shape	Lateral retaining wing A=high - R=regular					
	A	R	A	R	A	R
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52	1.013,00
Shell's shape	Lateral retaining wing B=regular			
	R	R	R	
	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

Configuration DEEP with side support PAD (height adjustable)







- Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)	Large (h 20cm)
	
<input type="checkbox"/> Right 10272 29,00	<input type="checkbox"/> Right 10273 40,00
<input type="checkbox"/> Left 10272 29,00	<input type="checkbox"/> Left 10273 40,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.

	<input type="radio"/> Large (18x12)	05647	183,00
	<input type="radio"/> Small (14x9)	05643	183,00
	<input type="radio"/> Lateral support extension	Quantity (number pcs).....	10271 69,00
	<input type="radio"/> Lateral support extension	Quantity (number pcs).....	
	<input type="radio"/> Small (14x9)	05642	183,00
	<input type="radio"/> Large (18x12)	05646	183,00

- Headrest adjustable - Double ball mounting system
Available only with backrest Physio
Not available for backrest Physio AIR carbon shell



- Large (34x14) **05648** **252,00**
- Small (28x14)

Date	Order	User	Address
------	-------	------	---------

REFERENCE TABLES

Progeo (Español) NOIR is a wheelchair without possibilities of adjustments after its manufacturing.

(1) The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

1		A Footplate width		
	DP	35÷39	40÷45	46÷50
	LS			
	33	21,5	21	20,5
	36	24,5	24	23,5
	39	27,5	27	26,5
	42	30,5	30	29,5
	45	33,5	33	32,5
	48	36,5	36	35,5

2	Setting	C :
	Extremely prudential	7,5
	Prudential	9
	Standard	10,5
	Active	12
	Extreme	13,5

3	Maximum width
Camber	
0°	LS + 18
2°	LS + 21
4°	LS + 24

Maximum width calculated including 3 cm. distance between tyre and handrim

Send this page only if any optional have been chosen	
☆	standard (supplied if no requests are indicated)
<input type="checkbox"/>	on request (without extra charge)
<input type="radio"/>	optional (cost to be added at the basic price)