




Date	UK (English)	
User	Company stamp	
End user weight (kg)		Folding frame wheelchair
Language user manual	Responsible	£ 3950,00
		Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm						
Seat width (LS)	☐ 33	☐ 36	☐ 39	☐ 42	☐ 45	☐ 48
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1 cm wider. For seat enlargement see section enlargement of mudguard and armrest with integrated mudguard						
Seat depth (PS)	☐ 35	☐ 37,5	☐ 40	☐ 42,5	☐ 45	☐ 47,5
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.						
Footplate distance (DP)	 (cm) ±0.5cm - adjustment steps 1.5cm (Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)					
Setting	☆ Standard ☐ Prudential ☐ Active					
Backrest Height (TS)	Std.Value 33,5(min 26-max 47) adjustment steps 1.5cm					
camber	☆ 0° ☐ 2°					
Front height STD 49						
Rear height STD 42						

Front frame		
☆ Aluminium standard frame		
		
☆	☐ (max 5" caster)	☐ (max 4" caster)

Frame colour		
Frame colour (without extra charge)		
☆ Opaque black	☐ POLAR white	☐ Matt orange
☐ Tekno anthracite	☐ Matt Midnight blue	☐ Matt grey
		☐ Matt red
Frame colour on request		
☐ Bright silver	☐ Bright fuchsia	☐ Bright blue china
☐ Bright blue	☐ Bright lime	☐ Bright red
		12579 109,00

Backrest
☆ Adjustable velcro straps and padded transpiring upholstery

Handle tube	
☆ Tubes with short handle	
☐ Tubes with standard handle	
☐ SPORT tubes Without push handles	
Lightest configuration	
☐ Folding handles	
Not available with: Push handles for companion	
	02958 270,00

Casters		
☆ 4" Soft		
☐ 3" Sport	Lightest configuration	
☐ 5" Soft		
☐ 6" Soft	(Available only with 100° angle)	
☐ 6" Pneumatic	(Available only with 100° angle)	
☐ 3" Active Soft (wide)	(Aluminium alloy rim)	06824 95,00
☐ 4" Active Sport (narrow)	(Aluminium alloy rim)	01398 95,00
☐ 4" Active Soft (wide)	(Aluminium alloy rim)	03943 95,00
☐ 5" Active Soft (wide)	(Aluminium alloy rim)	07518 95,00
☐ 5" Pneumatic		10844 68,00

Support front fork	
☆ Standard support front fork	
☐ Front fork support V Design 2.0	
Exclude the installation of : Quick release front forks	

Date Order User Address

Special Options

☐ Rear wheels plate extension	01950	189,00
☐ lockingband/strip to close/fold the wheelchair	06833	25,00
☐ Drawstring to unlock the back with a singleaction	06834	135,00

Colour set: includes rear wheel hub and rim and front forks

Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black	22" 24" 25" 26"		
☐ Bright silver	22" 24"	00833	135,00
☐ Bright lime	22" 24"	00833	135,00
☐ Bright blue	22" 24"	00833	135,00
☐ Bright red	22" 24"	00833	135,00
☐ Bright fuchsia	22" 24"	00833	135,00
☒ POLAR white	22" 24"	00833	135,00
☐ Tekno anthracite	22" 24"	00833	135,00
☐ Matt Midnight blue	22" 24"	00833	135,00
☐ Matt grey	22" 24"	00833	135,00
☐ Matt red	22" 24"	00833	135,00

Rear wheels

☆ 24" Ultra-light Progeo - Standard!			
☐ 22" Ultra-light Progeo - Standard!			
☐ 25" Ultra-light Progeo - Standard!	(colour: black)		
☐ 26" Progeo - Standard!	(colour: black)		
☐ 24" ULTRA-LIGHT PROGEO KEY SPORT	Lightest configuration	04383	695,00
☐ 24" Ultra-light "Spinergy LX"		01685	695,00
☐ 24" Ultra-light "Spinergy LXL"	☐ Black hub and white spokes	07846	695,00
☐ 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1.150,00
☐ 25" Ultra-light "Spinergy LX"		02101	695,00
☐ 25" Ultra-light "Spinergy LXL"	☐ Black hub and white spokes	07847	695,00
☐ 26" Ultra-light "Spinergy LX"		02102	695,00
☐ 24" Off-road (aluminium handrim black)		01392	270,00

Rear wheel tyre

☆ High pressure tyre	(22"-24"-25"-26")	Lightest configuration	
☐ High pressure profile tyre	(24")		
☐ Marathon Plus	(24"-25")	with anti-puncture protection	01669 95,00
☐ Anti-puncture - Colour Black	(24")		11083 95,00
Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"			
☐ Anti-puncture - Colour Grey	(22")		04181 105,00

Handrim**Wheel distance****☆ 3 cm****□ 4 cm**

☆ Aluminium	(22"-24"-25"-26")	Lightest configuration	
☐ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170 160,00
☐ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684 160,00
☐ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016 290,00
☐ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017 290,00
☐ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023 290,00
Oval countour with a rubber insert			
☐ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025 290,00
Ergonomic black profile			
☐ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026 290,00
Ergonomic profile with rubber coating (ideal for quadriplegics)			
☐ Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050 172,00

Wheel locks

☆ Bent lever push to lock			
☐ Straight lever push to lock			
☐ Straight lever pull to lock			
☐ Sport brake	Lightest configuration	06428	135,00
☐ Brake with extension lever push to lock		01949	68,00
☐ Brake with extension lever pull to lock		03372	68,00



Date Order User Address

Automatic Closure footplate


☆ Footplate COMPACT, automatic Closure

*Lightest configuration***11529**


Minimum footplate distance 35cm.

☐ Internal☐ External**RAKE AND DEPTH ADJUSTABLE****Footplate - See reference tables N°3**


Warning!! All flip-up footrests, need a stretch calf band and standard fork-support

Weight increase approx. 500gr.☐ Plastic footrest. rake adjustable

Minimum footplate distance 35cm.

☆ 2/3 Internal ☐ 2/3 External**11525**☐ One piece footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal☆ 2/3 Internal ☐ External**11527**☐ One piece aluminium footplate - Flip-up!. rake and depth adjustable

Minimum footplate distance 35cm.

☐ Internal☆ 2/3 Internal ☐ 2/3 External**11528**☐ One piece Carbon fibre footplate, rake adjustable, swing away, flip-up.

Minimum footplate distance 35cm.


☆ Internal

☐ External**03291****294,00****Mudguard****Mudguard - FIXED**

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard**Mudguard enlargement**

☆ + 0 cm

☐ + 1 cm☐ + 2 cm**12648**☐ Extra charge for : REMOVABLE VERSION

NO → Mudguard enlargement

**08921****135,00****Armrests**☐ Tubular armrests

(Swing away, removable and height adjustable armrest)

**00694****270,00****Removable armrest with integrated mudguard**☐ Desk armrests PROGEO

(Removable and height adjustable)

**09956****240,00****Enlargement of Desk armrests PROGEO**

☆ + 0 cm

☐ + 2 cm

Available only with Ultra-light Progeo standard wheels and camber 0°

12648☐ Removable armrest with carbon-fiber straight sideguard**09079****240,00****Accessories**☐ Pad for front frame protectors**07058****56,00**☐ Backrest-seat connection cloth**10241****25,00**☐ Removable storage bag (with quick coupling)**09811****135,00**

Date Order User Address

Optional 1

☐ Transit wheels			00836	158,00
☐ Anti-tipper with rapid hookup	(inside rotation to reduce the encumbrance)		00023	145,00
	☐ Right		00024	145,00
	☐ Left			
☐ Tipping aid	☐ Right		00025	79,00
	☐ Left		01046	79,00
☐ Crutches support	☐ Right		00838	79,00
	☐ Left		00837	79,00

Optional 2

☐ Tetraclip	(Quick release axle for quadriplegic)		03947	189,00
☐ Safety pelvic belt			03972	68,00
☐ Harness			04103	145,00
☐ Quick release front forks	Available with: Standard support front fork		00030	158,00
☐ Push handles for companion	(Removable and height adjustable)		01048	203,00
☐ Lumbar support	(to insert in the backrest upholstery)		00013	43,00
☐ Swing away lateral support	Not available with Posture backrest PHYSIO		02097	270,00
Prevents backrest folding	☐ Right		02098	270,00
	☐ Left			
☐ Spoke guards Black	(22"-24")		01043	158,00
☐ Spoke guards Transparent	(22"-24")		00835	158,00
☐ Spoke guards Progeo design	(22"-24")		00019	189,00

Additional Kit

☐ 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378	534,00
☐ 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357	895,00

Date.....

Order.....


User.....

Address.....

PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)


Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

☐  Carbon shell
not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7,5 cm 		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.121,00

☐  Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 	Lateral retaining wing R=regular		
	R	R	R
			

1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



☐ Right
10272

49,00

☐ Left
10272

49,00

Large (h 20cm)



☐ Right
10273


49,00

☐ Left
10273

49,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Aluminum lateral with padded, breathable upholstery.



☐ Large (18x12)

05647 196,00

☐ Small (14x9)

05643 196,00

☐ Lateral support extension

Quantity (number pcs)..... 10271 74,00

☐ Lateral support extension

Quantity (number pcs).....

☐ Small (14x9)

05642 196,00

☐ Large (18x12)

05646 196,00

☐ Headrest adjustable - Double ball mounting system
Available only with backrest Physio Posture



☐ Large (34x14) 05648 385,00

☐ Small (28x14)

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts. AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo EGO - Rev. 2025.1	Prices - VAT excluded to be added to the price list	Send this page only if any optional have been chosen	
		☒ standard	(supplied if no requests are indicated)
		☐ on request	(without extra charge)
		☐ optional	(cost to be added at the basic price)

REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK	Casters			
	3"	4"	5"	6"
SMALL 	max 46	max 47	max 48	
MEDIUM 	max 50	max 51	max 52	max 54
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.				

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 35	min 38	min 39	min 40	
max 43	max 45	max 46	max 48	

3	Footplate option	Footplate distance
	Footplate COMPACT, automatic Closure	min 36
	Plastic footrest	min 36
	One piece footplate - Flip-up!	min 36
	One piece aluminium footplate - Flip-up!	min 36
	One piece Carbon fibre footplate	min 36
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		
NOTE! Lower footplate distances are possible only after our valuation		
Calculated with the best configuration and with the footplates positioned on inside.		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all colours)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)