

Date

UK (English)

User

Company stamp

End user wheight (kg)

Folding tilting frame wheelchair

Language user manual

Responsible

£ 4250,00

Retail price - VAT not included

Information on end user

End user pathology :

Presence of dystonia or important spasms : ☐ YES ☐ NO

End user wheight (kg):.....Kg

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm

MODEL ADULT

☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48☐ kit to increase seat width (LS + 3cm)

09876

98,00

☆ Not mounted

☐ Mounted

Seat width (LS)

MODEL JUNIOR

☐ 27 ☐ 30 ☐ 33 ☐ 36 ☐ 39☐ kit to increase seat width (LS + 3cm)

11476


98,00

☆ Not mounted

☐ Mounted


Seat depth (PS) - adjustable!

Short frame

☐ 35 ☐ 37,5 ☐ 40PS 37,5
PS 35


(PS 40)

Medium frame

☐ 40 ☐ 42,5 ☐ 45PS 42,5
PS 40


(PS 45)

Long frame

☐ 45 ☐ 47,5 ☐ 50PS 47,5
PS 45

(PS 50)

Backrest tube position



Footplate distance (DP) - (Real distance between seat and footplate)

IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)

Minimum footplate distance ~22cm

.....(cm)

 ± 0.5 cm - adjustment steps 1.5cm

Seat height (parallel to the ground)


CASTERS 6"

Rear wheels			
12"	16"	22"	24"
☐ 42		☐ 42	
☆ 43,3		☐ 43,3	☐ 43,3
☐ 44,5		☐ 44,5	☐ 44,5
☐ 45,8		☐ 45,8	☐ 45,8
☐ 47	☐ 47	☐ 47	☐ 47
	☐ 48,3	☐ 48,3	☐ 48,3


CASTERS 7"

Rear wheels			
12"	16"	22"	24"
☐ 44,5		☐ 44,5	☐ 44,5
☐ 45,8		☐ 45,8	☐ 45,8
☐ 47	☐ 47	☐ 47	☐ 47
	☐ 48,3	☐ 48,3	☐ 48,3



Front frame

☆ Straight frame
(not available for seat width 27)


Front frame, swing away, and removable




☆ 110°

☐ 120°☐ 95°

(Exclude the installation of casters 7")

☐ Frame with abduction 6cm
not available for seat width 42 - 45 - 48
Maximum footplate distance 42

Front frame, swing away, and removable

☆ 95°
(Not available casters 7")☐ Elevating front frame, swing away and removable
not available for seat width 27 - 30 - 33
Minimum footplate distance 25cm.Include calf supports and
plastic footrest with heel
band

09730

270,00

Date..... Order..... User..... Address.....

Frame colour (without extra charge)

☆ Opaque black

☐ AVIO blue☐ POLAR white☐ Matt orange☐ Tekno anthracite☐ Matt grey☐ Matt red**Frame colour on request**

○ Bright red

○ Bright fuchsia


○ Bright lime

○ Bright blue china


09729

270,00**One piece Push handle for companion**One piece Push handle with tilting control -
REMOVABLE

adjustable angle

With fixed back(non-folding back)
adjustable distance 30÷70cm
(to the seat)With folding back backrest
adjustable distance 20÷60cm
(to the seat)**Backrest**☆ Fixed back(non-folding back)
adjustable angle from -4° to +12°○ Folding back backrest
adjustable angle from -4° to +12°

09728

294,00○ Reclining backrest (Non-folding back)
dynamic reclining with gas spring

10090

337,00The real seat depth, from the front part of the
seat cloth to the backrest tube will be 1,5 cm
more than the declared one.**Backrest upholstery**

☆ WITH NO backrest upholstery, NOR tension adjustment straps

Backrest Height (TS) (std 42cm)

○ WITH backrest upholstery and tension adjustment straps. Complete with backrest-seat connection cloth.

04693

135,00

With fixed back(non-folding back)->minimum height 39cm

Backrest Height (TS) (std 42cm)**Casters**

☆ 6" Soft

☐ 6" Pneumatic☐ 7" Soft**Not available with: Front frame 95°****Rear wheels**

☆ 12" transit rear wheels - Anti-puncture

☐ 12" transit rear wheels - Pneumatic☐ 16" transit rear wheels - Pneumatic

○ 24" Ultra-light Progeo - Standard!

04696

135,00

○ 22" Ultra-light Progeo - Standard!

04696

135,00

○ 24" Ultra-light "Spinergy LX"

04691

695,00

○ 24" Ultra-light "Spinergy LXL"

☐ Black hub and white spokes☐ Hub and spokes black

07848

695,00

○ 24" Ultra-light "Spinergy Flex rim"

Integrated aluminium/rubber handrim

02270

1.150,00

○ 24" Off-road (aluminium handrim black)

04690

337,00○ 12" with brake handles for companion (drum brake) ☐ Anti-puncture☐ Pneumatic

04105

349,00

○ 16" with brake handles for companion (drum brake) - Pneumatic

04107

349,00

○ 24" with brake handles for companion (drum brake) (aluminium handrim - Wheel distance 3cm)

03384

386,00**Colour rear wheels**

Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black 22" 24"

00833

135,00

○ Bright silver 22" 24"

00833

135,00

○ Bright lime 22" 24"

00833

135,00

○ Bright blue 22" 24"

00833

135,00

○ Bright red 22" 24"

00833

135,00

○ Bright fuchsia 22" 24"

00833

135,00

○ POLAR white 22" 24"

00833

135,00

○ Tekno anthracite 22" 24"

00833

135,00

○ Matt Midnight blue 22" 24"

00833

135,00

○ Matt grey 22" 24"

00833

135,00

○ Matt red 22" 24"

00833

135,00**Handrim****Wheel distance ☆ 3 cm**☐ 4 cm

☆ Aluminium (22-24")

○ Aluminium black SOFT-TOUCH (24")

09170

160,00

○ Chromed Steel (24")

01684

160,00

○ Titanium (24")

00016

290,00

○ Ultra Grepp (24")

00017

290,00

○ TOP FRICTION (24")

11023

290,00**Oval countour with a rubber insert**

○ NOVA H (24")

11025

290,00**Ergonomic black profile**

○ NOVA D (24")

11026

290,00**Ergonomic profile with rubber coating (ideal for quadriplegics)**

○ Silicon cover for standard aluminium handrim

04050

172,00

Date Order User Address

Rear wheel tyre

- ☆ High pressure tyre (22-24")
- ☐ High pressure profile tyre (24")
- ☐ Marathon Plus (24") with anti-puncture protection 01669 **95,00**
- ☐ Anti-puncture - Colour Grey (22") 04181 **105,00**
- ☐ Anti-puncture - Colour Black (24") 11083 **95,00**
- Not available with Rear wheels : 24" Ultra-light "Spinergy LXL"**

Wheel locks

- ☆ Bent lever push to lock
- ☐ Straight lever push to lock
- ☐ Straight lever pull to lock
- ☐ Brake with extension lever push to lock 01949 **68,00**
- ☐ Brake with extension lever pull to lock 03372 **68,00**

Footplate (Not available with Elevating front frame, swing away and removable)

- ☆ Plastic footrest, rake adjustable
Minimum footplate distance ~22cm
Not available with: Front frame with levating footrest, swing away and removable
- ☐ One piece aluminium footplate - Flip-up! rake and depth adjustable
Not available for seat width 27-30 with straight front frame
Minimum footplate distance ~21cm
Not available with: Front frame with levating footrest, swing away and removable
- ☐ Internal ☆ 2/3 Internal ☐ 2/3 External ☐ External

Removable armrest with integrated mudguard

- ☆ Desk armrests PROGEO (Removable and height adjustable)
- ☐ Removable armrest with ABS straight sideguard (Removable and height adjustable)
- Real distance between seat and armrest pads (cm) ☐ 21 ☐ 23 ☐ 25 ☆ 27 ☐ 29

SAFETY

- ☐ 4 points tie-down system rings
(installation must be carried out by the manufacturer)
- A head restraint suitable for transportation must be available and fitted accordingly during transportation.

08654 **203,00****Accessories**

- ☐ Pad for front frame protectors 07058 **56,00**
- ☐ Removable storage bag (with quick coupling) 09811 **135,00**

Optional 1

- ☆ Anti-tipper with rapid hookup (inside rotation to reduce the encumbrance)
INCLUDED ON THE BASIC PRICE
- ☐ Super Head rest (Folding with the wheelchair) 11477 **488,00**
- ☐ Transit wheels 00836 **158,00**
- ☐ Tipping aid ☐ Right 00025 **79,00**
☐ Left 01046 **79,00**
- ☐ Crutches support ☐ Right 00838 **79,00**
☐ Left 00837 **79,00**

Optional 2

- ☐ Harness 04103 **145,00**
- ☐ Safety pelvic belt 03972 **68,00**
- ☐ Tray table (removable) 00508 **294,00**
- ☐ Lumbar support (to insert in the backrest upholstery) 00013 **43,00**
- ☐ Spoke guards Black (22-24") 01043 **158,00**
- ☐ Spoke guards Transparent (22-24") 00835 **158,00**
- ☐ Spoke guards Progeo design (22-24") 00019 **189,00**

Date.....Order.....User.....Address.....

PHYSIO - Std. Backrest upholstery not provided

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

Posture backrest KID JUNIOR VERSION

- ☐ *KID* Aluminum shell

WARNING: with this backrest the seat depth is reduced by 5 cm approx.

BACKREST HEIGHT	29	34	40
Shell's shape 	Lateral retaining wing R=regular		
	R ☐ 	R ☐ 	R ☐ 

1.085,00

Posture backrest ADULT VERSION

- ☐ *AIR* Carbon shell
not available for seat width 27 - 30

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7,5 cm 		A	R	A	R	A	R
							
		☐	☐	☐	☐	☐	☐

1.121,00

- ☐ *POSTURE* Aluminum shell
not available for seat width 27 - 30

BACKREST HEIGHT	34	40	52
Shell's shape 	Lateral retaining wing R=regular		
	R ☐ 	R ☐ 	R ☐ 

1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



- ☐ Right
10272

49,00



- ☐ Left
10272

49,00

Large (h 20cm)



- ☐ Right
10273

49,00



- ☐ Left
10273

49,00

Date.....

Order.....


User.....

Address.....

Accessories

KID AIR POSTURE

- ☐ SWING AWAY - Lateral support
- Not available for backrest Physio AIR carbon shell
- Fixing and adjustment hardware complete with protection sleeve
- Alluminum lateral with padded, breathable upholstery.



- ☐ Headrest adjustable - Double ball mounting system
- Available only with backrest Physio Posture and Kid

- ☐ Large (34x14)
- ☐ Small (28x14)



05647 196,00

05643 196,00

Quantity (number pcs).....
10271 74,00

Quantity (number pcs).....

05642 196,00

05646 196,00

05648 385,00

00064

Send this page only if any optional have been chosen	
☒ standard	(supplied if no requests are indicated)
☐ on request	(without extra charge)
☐ optional	(cost to be added at the basic price)