


permobil Order Form - Progeo DUKE

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	UK (English)	
User	Company stamp	
End user weight (kg)		Rigid frame wheelchair
Language user manual	Responsible	£ 4363,00
		Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm		
Seat width (LS)	☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48	
The seat width is measured from the outer sides of the frame, while the distance between the two side guards is 1.5 cm wider.		
Seat depth (PS)	☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5 ☐ 50	
The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.		
Footplate distance (DP)	 (cm) ±0.5cm - adjustment steps 1.5cm	
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		
Setting	★ Standard ☐ Prudential ☐ Active	
Backrest Height (TS)	Std.Value 33,5 (min 26-max 47) adjustment steps 1.5cm	
Camber	★ 0° ☐ 2° ☐ 4°	

Front frame	
★ Carbon fibre frame, Angle 90° to the seat	


Front height - Tick the box with the required height		
The possible difference in height between HIGH frame and LOW frame is 3cm (with same fork, caster and assembly hole).		
	 SPORT fork (H95)	 Small size fork (H115)
LOW FRAME 90° Front height (min 44,5 - max 50) Footplate distance (min 34,5) U-TUBE Aluminium footplate Footplate distance (min 18 ÷ max 35) DPmax = Front height - 10cm	Casters 3" ☐ ~45 ☐ ~46 4" ☐ ~46 ☐ ~47 5" ☐ ~48,5	Casters 3" ☐ ~45,5 ☐ ~47 4" ☐ ~46,5 ☐ ~48 ☐ ~49 5" ☐ ~49 ☐ ~50
HIGH FRAME 90° Front height (min 47,5 - max 53) Footplate distance (min 37,5) U-TUBE Aluminium footplate Footplate distance (min 21 ÷ max 38) DPmax = Front height - 10cm	Casters 3" ☐ ~48 ☐ ~49 4" ☐ ~49 ☐ ~50 5" ☐ ~51,5	Casters 3" ☐ ~48,5 ☐ ~50 4" ☐ ~49,5 ☐ ~51 ☐ ~52 5" ☐ ~52 ☐ ~53
The real front-seat height could be slightly different from the chosen one considering it is affected by seat depth and rear seat height.		
Rear height The difference between Front Height and Rear Height must be in the 3 - 12 cm range.		
Rear wheels 22" (min 36 - max 44) Rear wheels 24" (min 38 - max 46) Rear wheels 25" (min 39 - max 47) Rear wheels 26" (min 41 - max 48)		

Seat		
★ Seat cloth with pocket (under seat)		
☐ Seat with tension straps	08883	231,00
☐ Carbon rigid seat	08336	489,00
☐ Pocket (under seat) [25x11x4,5]	08281	76,00

Date Order User Address

CUSTOMIZE COLORS

- ① **Fork colour**
☆ Opaque black
☐ Bright silver
☐ Pearl white
☐ Bright red
☐ Tekno anthracite
☐ Matt grey
☐ Bright blue
☐ Bright fuchsia
☐ Bright lime
☐ Matt Midnight blue
☐ Matt red
- 08342 65,00
- ② **Colour rear wheels**
To be matched only with: Rear wheels Ultra-light Progeo - Standard! (22 - 24)
☆ Black
☐ Bright silver
☐ Pearl white
☐ Bright red
☐ Tekno anthracite
☐ Matt grey
☐ Bright blue
☐ Bright fuchsia
☐ Bright lime
☐ Matt Midnight blue
☐ Matt red
- 04692 65,00



Rear wheels axle

- ☆ Aluminum rear wheels axle
☐ Carbon fiber rear wheels axle *Lightest configuration*
- 07647 152,00

Backrest adjustment and fixing system

- ☆ Backrest with double locking (upright and folded) **with rake adjustable from -4° to +16°**
☐ Backrest with upright locking system ONLY (no locking when folded) **with rake adjustable from -4° to +16°**
☐ Fixed back(non-folding back) **with rake adjustable from -4° to +16°**

Backrest

- ☆ Adjustable velcro straps and padded transpiring upholstery

☆ Aluminium backrest

☐ Titanium backrest

☐ Carbon fibre backrest




06811 195,00

04389 771,00

- ☐ **AIR** Integrated posture backrest (NOT REMOVABLE)

Carbon shell - Lateral retaining 7,5cm
Depth adjustment ±1cm. - Angle adjustment ±10°
not available for seat width 48

H 28		H 34		H 40	
Backrest Height (TS) min 34 (.....)		Backrest Height (TS) min 38 (.....)		Backrest Height (TS) min 44 (.....)	
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
					
☐	☐	☐	☐	☐	☐


Installation on:

☆ Backrest tubes

☐ Aluminium backrest

☐ Titanium backrest

☐ Carbon fibre backrest



06811 195,00

04389 771,00

Handle tube

- ☆ Tubes with standard handle

☐ Tubes with short handle

☐ SPORT tubes Without push handles

- ☐ Carbon fibre tubes with short handle **Weight reduction approx. 110gr.**
Not available with: Posture backrest PHYSIO - Posture backrest A I R

04390 152,00

- ☐ Folding handles **Minimum backrest height 30,5cm**
Not available with: Push handles for companion




02958 260,00

Date Order User Address

Footplate


☆ Aluminium footplate. Rake not adjustable

☆ WITHOUT plastic plate

☐ with plastic plate☐ One piece aluminium footplate - Flip-up! FLIP-BACK! - rake and depth adjustable

IMPORTANT NOTICE! From the footplate to the floor minimum 6cm !! (Minimum Front height = DP + 4cm)

With 5" front wheel the difference between Front Height and Rear Height must be greater than 5cm




09420 195,00

☐ Aluminium footplate, rake and depth adjustable

☆ Internal

☐ 2/3 Internal☐ External☐ 2/3 External


04456 260,00

☐ U-TUBE - Aluminium footplate, rake and depth adjustable**LOW FRAME: minimum footplate distance 18 - maximum footplate distance 35****HIGH FRAME: minimum footplate distance 21 - maximum footplate distance 38**

09853 231,00

☐ One piece Carbon fibre footplate. Rake adjustable

☆ Internal

☐ External

04110 383,00

Mudguard

☆ Carbon fibre mudguard

☐ Carbon fibre Straight sideguard☐ Removable version

08921 130,00

Casters

☆ 4" Soft

☐ 3" Sport*Lightest configuration*☐ 5" Soft☐ 3" Active Soft (wide)

(Aluminium alloy rim)

06824 130,00

☐ 4" Active Sport (narrow)

(Aluminium alloy rim)

01398 130,00

☐ 4" Active Soft (wide)

(Aluminium alloy rim)

03943 130,00

☐ 5" Active Soft (wide)

(Aluminium alloy rim)

07518 130,00

☐ 5" Pneumatic

10844 65,00

Rear wheels

☆ 24" Ultra-light Progeo - Standard!

☐ 22" Ultra-light Progeo - Standard!☐ 25" Ultra-light Progeo - Standard!

(colour: black)

☐ 26" Progeo - Standard!

(colour: black)

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT*Lightest configuration*

04383 902,00

☐ 24" Ultra-light "Spinergy LX"

01685 902,00

☐ 24" Ultra-light "Spinergy LXL"☐ Black hub and white spokes☐ Hub and spokes black

07846 902,00

☐ 24" Ultra-light "Spinergy Flex rim"

Integrated aluminium/rubber handrim

02270 1.415,00

☐ 25" Ultra-light "Spinergy LX"

02101 902,00

☐ 25" Ultra-light "Spinergy LXL"☐ Black hub and white spokes

07847 902,00

☐ 26" Ultra-light "Spinergy LX"

02102 902,00

☐ 24" Off-road (aluminium handrim black)

01392 231,00

Rear wheel tyre

☆ High pressure tyre

(22"-24"-25"-26")

Lightest configuration☐ High pressure profile tyre

(24")

☐ Marathon Plus

(24"-25")

with anti-puncture protection

01669 65,00

☐ Anti-puncture - Colour Grey

(22")

04181 101,00

☐ Anti-puncture - Colour Black

(24")

11083 130,00

Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"

Date Order User Address


Handrim	Wheel distance ★ 3 cm	□ 4 cm		
☆ Aluminium	(22"-24"-25"-26")	Lightest configuration		
○ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	130,00
○ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	152,00
○ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	283,00
○ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	324,00
○ Max Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00834	359,00
○ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	411,00
○ Oval countour with a rubber insert				
○ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	271,00
○ Ergonomic black profile				
○ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	295,00
○ Ergonomic profile with rubber coating (ideal for quadriplegics)				
○ SURGE LT	(25")	Oval countour with a rubber insert	06427	371,00
○ Quad (24")	(24")		00693	206,00
○ Silicon cover for standard aluminium handrim		Available for measures (22"- 24" - 26")	04050	165,00

Wheel locks			
☆ Bent lever push to lock			
□ Straight lever push to lock			
□ Straight lever pull to lock			
○ Sport brake		06428	130,00
○ Brake with extension lever push to lock		01949	65,00
○ Brake with extension lever pull to lock		03372	65,00

Armrests			
Tubular armrests	(Swing away, removable and height adjustable armrest)		00694 260,00
○ Removable armrest with carbon-fiber straight sideguard			09079 231,00

Accessories			
○ Pad for front frame protectors		07058	54,00
○ Backrest-seat connection cloth		10241	18,00
○ Removable storage bag (with quick coupling)		09811	130,00

Optional 1			
<u>Only one optional devices for side can be chosen</u>			
○ Transit wheels	Exclude the installation of Anti-tippers		03290 152,00
○ Anti-tipper with rapid hookup	Exclude the installation of Transit wheels (inside rotation to reduce the encumbrance)		06821 130,00
	○ Right		06822 130,00
	○ Left		

Optional 2			
○ Push handles for companion	(Removable and height adjustable)		03645 182,00
○ Tetraclick	(Quick release axle for quadriplegic)		03947 182,00
○ Safety pelvic belt			03972 65,00
○ Harness			04103 130,00
○ Lumbar support	(to insert in the backrest upholstery)		00013 41,00
○ Swing away lateral support	○ Right		02097 260,00
Not available with:	○ Left		02098 260,00
Posture backrest PHYSIO - Posture backrest A I R			
○ Spoke guards Black	(22"-24")		01043 152,00
○ Spoke guards Transparent	(22"-24")		00835 152,00
○ Spoke guards Progeo design	(22"-24")		00019 182,00
○ Spoke guards Dragon	(22"-24")		02106 182,00


Additional Kit			
○ 24" Off-road rear wheels	Off-Road tyre - aluminium handrim (black colour)	03378	513,00
○ 24" FAT rear wheels	FAT tyre 3" - aluminium handrim (black colour)	09357	902,00

Date.....Order.....User.....Address.....

PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- ☐ **AIR** Carbon shell
- not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape 7,5 cm 		Lateral retaining wing A=high - R=regular					
		A	R	A	R	A	R
		 ☐	 ☐	 ☐	 ☐	 ☐	 ☐

1.078,00

- ☐ **POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 7.5 cm 	Lateral retaining wing R=regular		
	R	R	R
	 ☐	 ☐	 ☐

1.043,00

Configuration DEEP with side support PAD (height adjustable)

- ☐ Alluminum lateral with padded, breathable upholstery.



- ☐ Right
10272
- ☐ Left
10272

30,00
30,00




- ☐ Right
10273
- ☐ Left
10273

41,00
41,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



- ☐ Large (18x12)
- ☐ Small (14x9)
- ☐ Lateral support extension
- ☐ Lateral support extension
- ☐ Small (14x9)
- ☐ Large (18x12)

05647	188,00
05643	188,00
Quantity (number pcs).....	10271 71,00
Quantity (number pcs).....	
05642	188,00
05646	188,00

- ☐ Headrest adjustable - Double ball mounting system

Available only with backrest Physio Posture



- ☐ Large (34x14)
- 05648 295,00
- ☐ Small (28x14)

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts. AL7 1AN United Kingdom - Phone: +44 1484 722888

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo DUKE - Rev. 2024.1	Prices - VAT excluded	Send this page only if any optional have been chosen	
	to be added to the price list	☐ standard	(supplied if no requests are indicated)
		☐ on request	(without extra charge)
		☐ optional	(cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb