


Date	UK (English)	Folding frame wheelchair
User	Company stamp	
End user wheight (kg)		£ 1900,00
Language user manual	Responsible	Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm	
Seat width (LS)	☐ 33 ☐ 36 ☐ 39 ☐ 42 ☐ 45 ☐ 48
The seat width is measured from the outer sides of the frame	
For seat enlargement see section enlargement of mudguard and armrest with integrated mudguard	
Seat depth (PS)	☐ 35 ☐ 37,5 ☐ 40 ☐ 42,5 ☐ 45 ☐ 47,5
Footplate distance (DP)	 (cm) ± 0.5 cm - adjustment steps 1.5cm
(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)	
Setting	☒ Standard ☐ Prudential ☐ Active ☐ Extreme
Backrest Height (TS)	Std.Value 36(min 30-max 54) adjustment steps 1.5cm
Crossbrace (black)	☒ Single ☐ Double (Not seat width 33-36)
Front height (STD 49) - Rear height (STD 44)	

Rear frame	
Aluminium frame	☒ Aluminium standard frame - Same colour of front frame  The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. camber ☒ 0° ☐ 1,5° ☐ 3°
	☐ Aluminium frame DYNAMIC 2.0 - (COLOUR: BLACK) 09955 337,00  The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one. With rake adjustable backrest ($\pm 9^\circ$) Exclude the installation of : Brake handles for companion (Drum Brake) camber ☒ 0° ☐ 1,5° ☐ 2° ☐ 4°
Titanium frame TMAX	 Front frame + Titanium frame T MAX 06858 736,00 The real seat depth, from the front part of the seat cloth to the backrest tube will be 1 cm more than the declared one. <u>Brushed Titanium frame - not varnished</u> <u>Weight reduction approx. 450gr</u> camber ☒ 0° ☐ 1,5° ☐ 3°

Front frame	
☐ Short frame - Ideal with Seat depth (PS) 35 - 37,5 Not available for depth 45 - 47,5 ☐ Medium frame - Ideal with Seat depth (PS) 40 - 42,5 Not available for depth 47,5 ☐ Long frame - Ideal with Seat depth (PS) 45 - 47,5	Frame Angle ☒ 110° ☐ 100°  

Frame colour	
Frame colour (without extra charge) ☒ Opaque black ☐ Tekno anthracite ☐ Matt grey ☒ POLAR white ☐ Matt Midnight blue ☐ Matt red	
Frame colour on request ☐ Bright silver ☐ Bright fuchsia ☐ Bright blue china ☐ Bright blue ☐ Bright lime ☐ Bright red	12579 109,00

Backrest
☒ Adjustable velcro straps and padded transpiring upholstery

Handle tube
☒ Tubes with standard handle ☐ Folding handles Not available with: Push handles for companion, Brake handles for companion (Drum Brake)
 02958 270,00

Date	Order	User	Address
------	-------	------	---------

Colour set: Includes rear wheel hub and rim, rear wheel plates and front forks

Only standard Progeo 22"/24" wheels have coloured hubs and rims

☆ Black	22" 24" 25" 26"		
☐ Bright silver	22" 24"	00833	135,00
☐ Bright lime	22" 24"	00833	135,00
☐ Bright blue	22" 24"	00833	135,00
☐ Bright red	22" 24"	00833	135,00
☐ Bright fuchsia	22" 24"	00833	135,00
☒ POLAR white	22" 24"	00833	135,00
☐ Tekno anthracite	22" 24"	00833	135,00
☐ Matt Midnight blue	22" 24"	00833	135,00
☐ Matt grey	22" 24"	00833	135,00
☐ Matt red	22" 24"	00833	135,00

Casters

☆ 5" Soft			
☐ 4" Soft			
☐ 3" Sport	Lightest configuration		
☐ 6" Soft	(Available only with 110° angle)		
☐ 6" Pneumatic	(Available only with 110° angle)		
☐ 7" Soft	(Available only with 110° angle)		
☐ 3" Active Soft (wide)	(Aluminium alloy rim)	06824	95,00
☐ 4" Active Sport (narrow)	(Aluminium alloy rim)	01398	95,00
☐ 4" Active Soft (wide)	(Aluminium alloy rim)	03943	95,00
☐ 5" Active Soft (wide)	(Aluminium alloy rim)	07518	95,00
☐ 5" Pneumatic		10844	68,00

Rear wheels

☆ 24" Ultra-light Progeo - Standard!			
☐ 22" Ultra-light Progeo - Standard!			
☐ 25" Ultra-light Progeo - Standard!	(colour: black)		
☐ 26" Progeo - Standard!	(colour: black)		
☐ 24" ULTRA-LIGHT PROGEO KEY SPORT	Lightest configuration	04383	695,00
☐ 24" Ultra-light "Spinergy LX"		01685	695,00
☐ 24" Ultra-light "Spinergy LXL"	☐ Black hub and white spokes	07846	695,00
☐ 24" Ultra-light "Spinergy Flex rim"	Integrated aluminium/rubber handrim	02270	1.150,00
☐ 25" Ultra-light "Spinergy LX"		02101	695,00
☐ 25" Ultra-light "Spinergy LXL"	☐ Black hub and white spokes	07847	695,00
☐ 26" Ultra-light "Spinergy LX"		02102	695,00
☐ 24" Off-road (aluminium handrim black)		01392	270,00
☐ 24" with brake handles for companion (drum brake) (aluminium handrim - Wheel distance 3cm)		03384	386,00
Available only with camber 0° Not available with Frame Dynamic 2.0°			
☐ One arm drive - (24" rear wheels Aluminium handrim High pressure Tyre) - camber 0°			
	☐ Right ☐ Left	02103	430,00

Rear wheel tyre

☆ High pressure tyre	(22"-24"-25"-26")	Lightest configuration	
☐ High pressure profile tyre	(24")		
☐ Marathon Plus	(24"-25")	with anti-puncture protection	01669 95,00
☐ Anti-puncture - Colour Black	(24")		11083 95,00
Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"			
☐ Anti-puncture - Colour Grey	(22")		04181 105,00

Date Order User Address

Handrim	Wheel distance ☆ 3 cm	□ 4 cm		
☆ Aluminium	(22"-24"-25"-26")	Lightest configuration		
○ Aluminium black SOFT-TOUCH	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	09170	160,00
○ Chromed Steel	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	01684	160,00
○ Titanium	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00016	290,00
○ Ultra Grepp	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	00017	290,00
○ TOP FRICTION	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11023	290,00
Oval countour with a rubber insert				
○ NOVA H	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11025	290,00
Ergonomic black profile				
○ NOVA D	(24"-25")	(24" handrim, on 25" wheel is installed with spacers)	11026	290,00
Ergonomic profile with rubber coating (ideal for quadriplegics)				
○ Silicon cover for standard aluminium handrim	Available for measures (22"- 24" - 26")		04050	172,00

Wheel locks				
☆ Bent lever push to lock				
□ Straight lever push to lock				
□ Straight lever pull to lock				
○ Sport brake				
			06428	135,00
○ Brake with extension lever push to lock				
○ Brake with extension lever pull to lock				
			01949	68,00
			03372	68,00

Footplate - See reference tables N°3				
☆ Plastic footrest, rake adjustable (Minimum footplate distance 33cm.)	☆ 2/3 Internal	□ 2/3 External		11525
□ Automatic Closure footplate , rake and depth adjustable (Minimum footplate distance 30cm.)	□ Internal	☆ 2/3 Internal	□ External	
○ Minimum footplate distance 25cm. Available only with Multi-function and rotative system - 03945 203,00				
□ One piece footplate - Flip-up! rake and depth adjustable (Minimum footplate distance 30cm.)	□ Internal	☆ 2/3 Internal	□ External	
○ Minimum footplate distance 25cm. Available only with Multi-function and rotative system - 03945 203,00				
□ One piece aluminium footplate - Flip-up! rake and depth adjustable (Minimum footplate distance 30cm.)	□ Internal	☆ 2/3 Internal	□ 2/3 External	□ External
○ One piece Carbon fibre footplate rake adjustable, swing away, flip-up. (Minimum footplate distance 30cm.)	☆ Internal	□ External		03291 294,00
○ Footplate COMPACT, automatic Closure, rake and depth adjustable Minimum footplate distance 35cm.	☆ Internal	□ External		09957 270,00

Mudguard				
Mudguard - FIXED				
☆ ABS mudguard				
□ ABS Straight sideguard				
○ Carbon fibre mudguard				03946 337,00
○ Carbon fibre Straight sideguard				08918 337,00
Mudguard enlargement				
☆ + 0 cm	Not available with Rear wheels with brake handles for companion (Drum Brake)			
□ + 1 cm □ + 2 cm				12648
Mudguard - REMOVABLE				
The distance between the two side guards is 1 cm wider than the seat width				
○ Removable ABS mudguard				08921 135,00
○ Removable ABS fibre straight sideguard				08921 135,00
○ Removable carbon fibre mudguard				08919 466,00
○ Removable Carbon fibre straight sideguard				08920 466,00

Date	Order	User	Address
Armrests <i>Tubular armrests with standard support combined with anti-tipper wheels exclude the possibility of installing the transit wheels.</i>			
☐ Tubular armrests		(Swing away, removable and height adjustable armrest) ☐ Standard Support (Not with frame Dynamic) ☐ Steel support with safety locking system	 04828 189,00 00694 270,00
Removable armrest with integrated mudguard			
☐ Desk armrests PROGEO		(Removable and height adjustable) Not available with Frame Dynamic 2.0 + camber 4°	 09956 240,00
Enlargement of Desk armrests PROGEO			
☒ + 0 cm ☐ + 2 cm		Available only with Ultra-light Progeo standard wheels and camber 0°	12648
☐ Flip-up armrests PROGEO		(Removable and height adjustable) Not available with Frame Dynamic 2.0	 12277 240,00
☐ Removable armrest with ABS straight sideguard <i>The distance between the two side guards is 1 cm wider than the seat width</i>			 08881 240,00
☐ Removable armrest with carbon-fiber straight sideguard <i>The distance between the two side guards is 1 cm wider than the seat width</i>			08882 486,00
Accessories			
☐ Pad for front frame protectors			07058 56,00
☐ Backrest-seat connection cloth			10241 25,00
☐ Removable storage bag (with quick coupling)			09811 135,00
Optional 1			
☐ Transit wheels			00836 158,00
☐ Anti-tipper with rapid hookup		(inside rotation to reduce the encumbrance) ☐ Right ☐ Left	 00023 145,00 00024 145,00
☐ Tipping aid		☐ Right ☐ Left	 00025 79,00 01046 79,00
☐ Crutches support		☐ Right ☐ Left	 00838 79,00 00837 79,00
Optional 2			
☐ Tetraclick		(Quick release axle for quadriplegic)	03947 189,00
☐ Rear wheels plate extension ☐ Handybike adaptor		(Not with frame Dynamic 2.0)	01950 189,00 00506 96,00
☐ Push handles for companion		(Removable and height adjustable)	 01048 203,00
☐ Harness ☐ Safety pelvic belt			04103 145,00 03972 68,00
☐ Tray table (removable)			00508 294,00
☐ Lumbar support		(to insert in the backrest upholstery)	00013 43,00
☐ Super Head rest		(Folding with the wheelchair)	 04049 486,00
☐ Swing away lateral support Not available with Posture backrest PHYSIO		☐ Right ☐ Left	 02097 270,00 02098 270,00
☐ Spoke guards Black		(22"-24")	01043 158,00
☐ Spoke guards Transparent		(22"-24")	00835 158,00
☐ Spoke guards Progeo design		(22"-24")	00019 189,00
Additional Kit			
☐ 24" Off-road rear wheels		Off-Road tyre - aluminium handrim (black colour)	03378 534,00
☐ 24" FAT rear wheels		FAT tyre 3" - aluminium handrim (black colour)	09357 895,00

Date.....

Order.....

User.....

Address.....

PHYSIO - Posture backrest - Std. Backrest upholstery not provided




Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

☐ **AIR** Carbon shell
not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7.5 cm 	A	R	A	R	A	R	
							
	☐	☐	☐	☐	☐	☐	

1.121,00

☐ **POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape	Lateral retaining wing R=regular		
			
	☐	☐	☐

1.085,00

Configuration DEEP with side support PAD (height adjustable)

Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)



☐ Right
10272

49,00

☐ Left
10272

49,00

Large (h 20cm)



☐ Right
10273


49,00

☐ Left
10273

49,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



☐ Large (18x12)

05647 196,00

☐ Small (14x9)

05643 196,00

☐ Lateral support extension

Quantity (number pcs).....
10271 74,00

☐ Lateral support extension

Quantity (number pcs).....

☐ Small (14x9)

05642 196,00

☐ Large (18x12)

05646 196,00

☐ Headrest adjustable - Double ball mounting system
Available only with backrest Physio Posture

☐ Large (34x14)
☐ Small (28x14)



05648 385,00

Permobil LimitedUnit 7, Polaris Centre, 41 Brownfields,Welwyn Garden City, Herts. AL7 1ANUnited Kingdom - Phone: +44 1484 722888

Date	Order	User	Address
------------	-------------	------------	---------------

Note

Wheelchair : Progeo EXELLE - Rev. 2025.1	Prices - VAT excluded to be added to the price list	Send this page only if any optional have been chosen	
		☒ standard	(supplied if no requests are indicated)
		☐ on request	(without extra charge)
		☐ optional	(cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

REFERENCE TABLES

1	Front height	
With closed frame (95°) max 5" caster		
Casters	Front height	
	min.	max
3"	41	52
4"	42	53
5"	44	54
6"	48	55
7"	50	56
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.		

2	Rear height				
STANDARD	Rear frame	Rear wheels			
		22"	24"	25"	26"
	Straight plate 	min. 37 max 45	min. 39 max 47	min. 40 max 48	min. 41 max 49
	Shaped plate 	min. 36 max 49	min. 38 max 51	min. 39 max 52	min. 40 max 53
	Rear wheels plate extension with Tubular armrests 	min. 37 min. 39 max 46	min. 39 min. 41 max 48	min. 40 min. 43 max 50	min. 41 min. 44 max 51
DYNAMIC 2.0		min. 36 max 49	min. 38 max 51	min. 39 max 52	min. 40 max 53

3	Footplate option	Footplate distance	
		Frame Angle 110°	Frame Angle 100°
	Plastic footrest	min. 33	min. 31
	Automatic Closure footplate	min. 30	min. 28
	Automatic Closure footplate , with multifunction rotative system	min. 25	min. 23
	One piece footplate - Flip-up!	min. 30	min. 28
	One piece footplate - Flip-up! , with multifunction rotative system	min. 25	min. 23
	One piece aluminium footplate - Flip-up!	min. 30	min. 28
	One piece Carbon fibre footplate	min. 30	min. 28
Calculated with the best configuration and with the footplates positioned on inside.			
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)			

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all coulors)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)

Wheelchair : Progeo EXELLE - Rev. 2025.1