

permobil Order Form - Progeo JOKER

Permobil Limited Unit 7, Polaris Centre, 41 Brownfields, Welwyn Garden City, Herts, AL7 1AN United Kingdom - Phone: +44 1484 722888

Date

User

End user wheight (kg)

Language user manual

UK (English)

Company stamp

Responsible

Rigid frame wheelchair
CRASH TEST APPROVED

£ 2592,00

Retail price - VAT not included

Wheelchair measures - The tolerance on all the dimensions indicated and/or required is ± 10 mm

Seat width (LS)

☐ 33

☐ 36

☐ 39

☐ 42

☐ 45

☐ 48

Seat depth (PS)

☐ 35

☐ 37,5

☐ 40

☐ 42,5

☐ 45

☐ 47,5

☐ 50

The real seat depth, from the front part of the seat cloth to the backrest tube will be 1,5 cm more than the declared one.

Footplate distance (DP)

..... (cm) ±0.5cm - adjustment steps 1.5cm

(Real distance between seat and footplate) - IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)

Setting

☒ Standard

☐ Prudential

☐ Active

☐ Extreme

Backrest Height (TS)

Std.Value 33,5(min 26-max 47) adjustment steps 1.5cm

camber

☒ 0°

☐ 2°


☐ 4°

Front height STD 49

.....

Rear height STD 42

.....



☆ Aluminum rear adjustment plate

☆ Removable mudguard ABS

☐ Removable ABS fibre straight sideguard

☐ Removable carbon fibre mudguard


☐ Removable Carbon fibre straight sideguard

03946

324,00

08918

324,00



☐ FULL CARBON rear adjustment plate with integrated sideguard


☆ Version mudguard (Available only with rear wheels 24")

☐ Version straight sideguard

07646

967,00

JOKER



☆ Aluminium footplate, rake adjustable

Minimum footplate distance 35cm.

☆ Internal

☐ 2/3 Internal

☐ External

☐ 2/3 External

☐ Aluminium footplate EVOLUTION - rake and depth adjustable

Minimum footplate distance 35cm.

☆ Internal

☐ 2/3 Internal

☐ External

☐ 2/3 External

☐ One piece Carbon fibre footplate JOKER

Minimum footplate distance 35cm.

not available for seat width 48


Not available with: Frame Angle 85°

☐ One piece Carbon fibre footplate, rake adjustable

Minimum footplate distance 35cm.

☆ Internal

☐ External



11523



07649

130,00



03348


260,00



03346

307,00

JOKER ZERO



☐ Plastic footrest, rake adjustable

Minimum footplate distance 15cm.

☆ 2/3 Internal

☐ 2/3 External

☐ One piece aluminium footplate - Flip-up!, rake and depth adjustable

Minimum footplate distance 15cm.

☆ Internal

☐ 2/3 Internal

☐ External

☐ 2/3 External

☐ One piece Carbon fibre footplate rake adjustable, swing away, flip-up.

Minimum footplate distance 15cm.


☆ Internal

☐ External




06814

195,00



06856

195,00



06857

502,00

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Pag. 1 / 7



UK_202404

Rev. 2024.1_Apr

Date Order User Address

Front frame

☆ Aluminium standard frame



Frame colour

Frame colour (without extra charge)

☆ Opaque black
☐ Tekno anthracite☐ Pearl white
☐ Matt Midnight blue☐ Matt orange
☐ Matt grey☐ Matt red

Frame colour on request

☐ Bright silver
☐ Bright blue☐ Bright fuchsia
☐ Bright lime☐ Bright blue china
☐ Bright red

12579

105,00

Seat

☆ Seat cloth with pocket (under seat)

☐ Seat with tension straps

08883

231,00

☐ Carbon rigid seat

08336

489,00

☐ Pocket (under seat) [25x11x4,5]

08281

76,00

Rear wheels axle

☆ Aluminum rear wheels axle

☐ Carbon fiber rear wheels axle

Lightest configuration

07647

152,00

Backrest adjustment and fixing system

☆ Backrest with double locking (upright and folded)
with rake adjustable from -4° to +16°☐ Backrest with upright locking system ONLY (no locking when folded)
with rake adjustable from -4° to +16°☐ Fixed back(non-folding back)
with rake adjustable from -4° to +16°

Backrest

☆ Adjustable velcro straps and padded transpiring upholstery

☆ Aluminium backrest

☐ Titanium backrest

06811

195,00

☐ Carbon fibre backrest

04389

771,00

☐ **AIR** Integrated posture backrest (NOT REMOVABLE)

Carbon shell - Lateral retaining 7,5cm

Depth adjustment ±1cm. - Angle adjustment ±10°

not available for seat width 48

1.078,00

H 28	H 34	H 40			
Backrest Height (TS) min 34 (.....)	Backrest Height (TS) min 38 (.....)	Backrest Height (TS) min 44 (.....)			
Lateral retaining wing A=high - R=regular					
A	R	A	R	A	R
					
☐	☐	☐	☐	☐	☐

Installation on:

☆ Backrest tubes

☐ Aluminium backrest☐ Titanium backrest☐ Carbon fibre backrest

06811

195,00

04389

771,00

Date Order User Address

Handle tube

★ Tubes with standard handle

- ☐ Tubes with short handle
☐ SPORT tubes Without push handles

☐ Carbon fibre tubes with short handle *Lightest configuration*
Not available with: Posture backrest PHYSIO - Posture backrest A I R

04390 152,00

☐ Folding handles **Minimum backrest height 30,5cm**
Not available with: Push handles for companion



02958 260,00

Casters

★ 4" Soft

☐ 3" Sport☐ 5" Soft☐ 6" Soft☐ 6" Pneumatic*Lightest configuration*

(Available only with 100° angle)

(Available only with 100° angle)

☐ 3" Active Soft (wide)

(Aluminium alloy rim)

06824 130,00

☐ 4" Active Sport (narrow)

(Aluminium alloy rim)

01398 130,00

☐ 4" Active Soft (wide)

(Aluminium alloy rim)

03943 130,00

☐ 5" Active Soft (wide)

(Aluminium alloy rim)

07518 130,00

☐ 5" Pneumatic

10844 65,00

Colour set: includes rear wheel hub and rim and front forks

Only standard Progeo 22"/24" wheels have coloured hubs and rims

★ Black

22" 24" 25" 26"

☐ Bright silver

22" 24"

00833 130,00

☐ Bright lime

22" 24"

00833 130,00

☐ Bright blue

22" 24"

00833 130,00

☐ Bright red

22" 24"

00833 130,00

☐ Bright fuchsia

22" 24"

00833 130,00

☐ Pearl white

22" 24"

00833 130,00

☐ Tekno anthracite

22" 24"

00833 130,00

☐ Matt Midnight blue

22" 24"

00833 130,00

☐ Matt grey

22" 24"

00833 130,00

☐ Matt red

22" 24"

00833 130,00

Rear wheels

★ 24" Ultra-light Progeo - Standard!

☐ 22" Ultra-light Progeo - Standard!☐ 25" Ultra-light Progeo - Standard!☐ 26" Progeo - Standard!

(colour: black)

(colour: black)

☐ 24" ULTRA-LIGHT PROGEO KEY SPORT*Lightest configuration*

04383 902,00

☐ 24" Ultra-light "Spinergy LX"

01685 902,00

☐ 24" Ultra-light "Spinergy LXL"☐ Black hub and white spokes☐ Hub and spokes black

07846 902,00

☐ 24" Ultra-light "Spinergy Flex rim"

Integrated aluminium/rubber handrim

02270 1.415,00

☐ 25" Ultra-light "Spinergy LX"

02101 902,00

☐ 25" Ultra-light "Spinergy LXL"☐ Black hub and white spokes

07847 902,00

☐ 26" Ultra-light "Spinergy LX"

02102 902,00

☐ 24" Off-road (aluminium handrim black)

01392 231,00

Rear wheel tyre

★ High pressure tyre

(22"-24"-25"-26")

Lightest configuration☐ High pressure profile tyre

(24")

☐ Marathon Plus

(24"-25")

with anti-puncture protection

01669 65,00

☐ Anti-puncture - Colour Grey

(22")

04181 101,00

☐ Anti-puncture - Colour Black

(24")

11083 130,00

Not available with Rear wheels : 24" ULTRA-LIGHT PROGEO KEY SPORT, 24" Ultra-light "Spinergy LXL"

Handrim

Wheel distance ★ 3 cm

☐ 4 cm

★ Aluminium

(22"-24"-25"-26")

Lightest configuration☐ Aluminium black SOFT-TOUCH

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

09170 130,00

☐ Chromed Steel

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

01684 152,00

☐ Titanium

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00016 283,00

☐ Ultra Grepp

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00017 324,00

☐ Max Grepp

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

00834 359,00

☐ TOP FRICTION

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11023 411,00

Oval countour with a rubber insert

☐ NOVA H

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11025 271,00

Ergonomic black profile

☐ NOVA D

(24"-25")

(24" handrim, on 25" wheel is installed with spacers)

11026 295,00

Ergonomic profile with rubber coating (ideal for quadriplegics)

☐ SURGE LT

(25")

Oval countour with a rubber insert

06427 371,00

☐ Quad (24")

(24")

Available for measures (22"- 24" - 26")

00693 206,00

☐ Silicon cover for standard aluminium handrim

04050 165,00

Date Order User Address

Wheel locks

☆ Bent lever push to lock

☐ Straight lever push to lock

☐ Straight lever pull to lock

☐ Sport brake

Lightest configuration

☐ Brake with extension lever push to lock

☐ Brake with extension lever pull to lock



06428

130,00

01949

65,00

03372

65,00

Armrests

☐ Tubular armrests

(Swing away, removable and height adjustable armrest)



00694

260,00

Accessories

☐ Pad for front frame protectors

07058

54,00

☐ Backrest-seat connection cloth

10241

18,00

☐ Removable storage bag (with quick coupling)

09811

130,00

SAFETY

☐ 4 points tie-down system rings

The use of a headrest during transport is mandatory
(installation must be carried out by the manufacturer)



08654

195,00

Optional 1**Only one optional devices for side can be chosen**

☐ Transit wheels

Exclude the installation of Anti-tippers



03290

152,00

☐ Anti-tipper BASE

Exclude the installation of Transit wheels

☐ Right

06819

101,00

☐ Left

06820

101,00

☐ Anti-tipper with rapid hookup

(inside rotation to reduce the encumbrance)

Exclude the installation of Transit wheels

☐ Right

06821

130,00

☐ Left

06822

130,00

**Optional 2**

☐ Push handles for companion

(Removable and height adjustable)



03645

182,00

☐ Tetraclip

(Quick release axle for quadriplegic)

03947

182,00

☐ Safety pelvic belt

03972

65,00

☐ Harness

04103

130,00

☐ Lumbar support

(to insert in the backrest upholstery)

00013

41,00

☐ Swing away lateral support

Not available with Posture backrest PHYSIO - Posture backrest A I R

☐ Right

02097

260,00

☐ Left

02098

260,00



☐ Spoke guards Black

(22"-24")

01043

152,00

☐ Spoke guards Transparent

(22"-24")

00835

152,00

☐ Spoke guards Progeo design

(22"-24")

00019

182,00

☐ Spoke guards Dragon

(22"-24")

02106

182,00

Additional Kit

☐ 24" Off-road rear wheels

Off-Road tyre - aluminium handrim (black colour)

03378

513,00

☐ 24" FAT rear wheels

FAT tyre 3" - aluminium handrim (black colour)

09357


902,00

Date.....Order.....User.....Address.....

PHYSIO - Posture backrest
(Std. Backrest upholstery not provided)

Removable with quick release
Suitable for tubes with diameter 18÷26mm
Depth adjustment (4cm)
Angle adjustment (±18°)

- ☐ **AIR** Carbon shell
- not available for seat width 48

BACKREST HEIGHT		28		34		40	
Shell's shape		Lateral retaining wing A=high - R=regular					
7,5 cm 	A	R	A	R	A	R	
							
	☐	☐	☐	☐	☐	☐	

1.078,00

- ☐ **POSTURE** Aluminum shell

BACKREST HEIGHT	34	40	52
Shell's shape 7.5 cm 	Lateral retaining wing R=regular		
	R	R	R
			
	☐	☐	☐

1.043,00

Configuration DEEP with side support PAD (height adjustable)


- ☐ Alluminum lateral with padded, breathable upholstery.



Small (h 13cm)		Large (h 20cm)	
	☐ Right 10272		☐ Right 10273
	30,00		41,00
	☐ Left 10272		☐ Left 10273
	30,00		41,00

SWING AWAY - Lateral support

Fixing and adjustment hardware complete with protection sleeve
Alluminum lateral with padded, breathable upholstery.



☐ Large (18x12)	05647	188,00
☐ Small (14x9)	05643	188,00
☐ Lateral support extension	Quantity (number pcs) (.....)	
☐ Lateral support extension	10271	71,00
☐ Small (14x9)	05642	188,00
☐ Large (18x12)	05646	188,00

- ☐ Headrest adjustable - Double ball mounting system
- Available only with backrest Physio Posture



☐ Large (34x14)	05648	295,00
☐ Small (28x14)		

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

Permobil LimitedUnit 7, Polaris Centre, 41 Brownfields Welwyn Garden City, Herts. AL7 1ANUnited Kingdom - Phone: +44 1484 722888

Date Order User Address

Note

Wheelchair : Progeo JOKER - Rev. 2024.1		Prices - VAT excluded to be added to the price list		Send this page only if any optional have been chosen	
				☐	standard (supplied if no requests are indicated)
				☐	on request (without extra charge)
				☐	optional (cost to be added at the basic price)

customer_support_uk@permobil.com - https://www.permobil.com/en-gb

REFERENCE TABLES

1	Front height			
With closed frame (95°) max 5" caster With closed frame (90°) max 5" caster With closed frame (85°) max 4" caster				
FORK		Casters		
	3"	4"	5"	6"
SPORT 	max 52	max 53	max 54	
SMALL 	max 54	max 55	max 56	max 57
MEDIUM 	max 56	max 57	max 58	max 59
The size of the fork will be defined only by Progeo in accordance with wheelchair final configuration.				

2	Rear height			
Rear wheels				
22"	24"	25"	26"	
min 36	min 38	min 39	min 40	
max 44	max 46	max 47	max 48	

3	Footplate option	Footplate distance
	Aluminium footplate rake adjustable	min 35
	One piece Carbon fibre footplate	min 35
	One piece Carbon fibre footplate rake adjustable	min 35
	Plastic footrest, rake adjustable	min 35
	One piece aluminium footplate - Flip-up!, rake and depth adjustable	min 35
	One piece Carbon fibre footplate rake adjustable, swing away, flip-up.	min 35
IMPORTANT NOTICE! From the footplate to the floor minimum 4cm !! (Minimum Front height = DP + 4cm)		
Calculated with the best configuration and with the footplates positioned on inside.		
NOTE! Lower footplate distances are possible only after our valuation		

4	Rear wheels	
	Model	Features
	22" Ultra-light Progeo - Standard!	(not available in anthracite color)
	24" Ultra-light Progeo - Standard!	(available for all couloers)
	24" ULTRA-LIGHT PROGEO KEY SPORT	(silver hub, black rim and 16 spokes)
	24" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	24" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes) - (black hub, rim and 18 spokes)
	24" Ultra-light "Spinergy Flex rim"	(silver hub, black rim and spokes) - (integrated aluminium/rubber handrim - encumbrance reduction up to 3cm. - all tyre)
	24" Off-road	(colour: black) - (aluminium handrim - tyre Off-road)
	25" Ultra-light "Spinergy LX"	(silver hub, black rim and spokes)
	25" Ultra-light "Spinergy LXL"	(black hub and rim, 18 white spokes)
	26" Progeo - Standard!	(colour: black)
	26" Ultra-light "Spinergy LX"	(silver hub, black rim and 18 spokes)