

Q100 R & Q200 R

Retail



Billing Address Or Customer-ID	Shipping Address	Identical To Billing Address
C.-ID:	C.-ID:	
Name:	Name:	
Street:	Street:	
Town/PC:	Town/PC:	
Order Information	User Weight	
Purchase Order No.:	User Weight Capacity: Q100 R: 125 kg / Q200 R: 136 kg	
Art. No.: Q100 R : 22600100C & Q200 R - 22600000C-300		
Order Date:	Day	Month Year

Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	Sling Seat / Solid Seat Pan
Typical WC Weight Without Batteries	From 60 kg	Seat Size Range	Width 380 (with In-Rig Armrest/Legrest) - 540 mm
			Depth 410 - 510 mm

ORDER FORM

Drive Bases			£
QJR000019	Q100 R (RWD Base Chair) - 12" Drive Wheels, Solid & 8" Castors, Solid		£3,695
QJR000022	Q200 R (RWD Base Chair)		£3,795

Frame Width			
QJR020128	420 mm Frame Width – Seat Width Range With Armrest Adjustment 420-480 mm (380 mm Available With In-Rig Armrest)		0
QJR020129	480 mm Frame Width – Seat Width Range With Armrest Adjustment 480-540 mm (440 mm Available With In-Rig Armrest)		0

Seat Width		
SW38	380 mm (1)	
SW42	420 mm	
SW44	440 mm (2)	
SW46	460 mm	
SW48	480 mm	
SW52	520 mm	
SW54	540 mm	

Seat Depth		
SD41	410 mm	
SD435	435 mm	
SD46	460 mm	
SD485	485 mm	
SD51	510 mm	

(1) Supplied With In-Rig Armrest

(2) Supplied With In-Rig Armrest When 480 mm Frame Width Selected

Seat Type		£
QJR020087	Sling Seat	0
QJR020127	Solid Seat Pan (1)	140

Seat Cushion Type

QJR020098	Seat Cushion, Standard, 50 mm	0
-----------	-------------------------------	---

JAY Seat Cushion

JAY Seat Cushion		£
JAY020002	JAY Basic	80
JAY020006	JAY Easy Visco With Flat Base	180
JAY020015	JAY Balance, Fluid, Microclimatic, With Positioning Kit, Until SW 500 mm	755
JAY020009	JAY Ion	180
JAY020024	JAY Union	350
JAY090106	JAY Union, Reduced Profile	350

SEDEO Seat Cushion

QJR020073	SEDEO, Seat Cushion, Profiled, Height 70 mm	160
QJR020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm	260
QJR020078	SEDEO, Seat Cushion, Stabilizer, Height 90 mm	280

SEDEO Seat Cushion Cover

SEDEO Seat Cushion Cover		£
QJR020518	Seat Cushion Cover, Incontinence Protection	0
QJR020080	Seat Cushion Cover, Comfort	0

Power Seating Options

Seat Angle Adjustment (Adjustable To 0°, 3°, 6°) - Preset Selection

QJR020125	Seat Angle, 0° Preset	0
QJR020096	Seat Angle, 3° Preset	0
QJR020126	Seat Angle, 6° Preset (2)	0

Q200 R

QJR130013	Posterior Tilt, Powered, 0° - 22° (With 7 cm Centre Of Gravity Movement)	710
-----------	--	-----

Backrest Frame Type

QJR030147	Backrest, Tool Angle Adjustable 90°, 94°, 98°, 102°	0
-----------	---	---

Q200 R

QJR030160	Backrest, Powered Recline, 90° - 120°, Flip Down	560
-----------	--	-----

COG Adjustment

COG075	Centre Of Gravity, Adjustable, +0 mm to -75 mm (3)	340
--------	--	-----

Seat Heights

SH45	45cm, Seat to Floor Height (To pan 3° Fixed Tilt) - Range - 43 - 47cm (4)(5)	130
SH47	47cm, Seat to Floor Height (To pan 3° Fixed Tilt) - Range - 45 - 49cm (4)	0
SH49	49m, Seat Height Adjustment (To pan 3° Fixed Tilt) - Range - 49cm - 53cm (4) (5)	0

Backrest Height

BH50	500 mm	0
BH57	570 mm	100

Sling Back Size

QJR030150	Standard, 5 Tension Adj Backrest	0
-----------	----------------------------------	---

JAY Shape Back Size

JAY080006	JAY Shape Backrest - Oversized (Supplied 1" larger than seat width to allow greater immersion)	270
-----------	--	-----

(1) Recommended In Combination With JAY Seat Cushions
(2) Seat Angle 6° Not Available With Seat Height 490 mm
(3) Refer To Training Materials With COG Option Before Setup And Use
(4) 3° Seat Angle
(5) Not available with tilt

JAY Shape Laterals Thoracic Supports £

JAY080010	JAY SHAPE Lateral Thoracic supports (Pair) - Support 15 x 15 cm (Pair)	270
JAY080011	JAY SHAPE Lateral Thoracic supports (Pair) - Support 10 x 12 cm (Pair)	270

SEDEO Lite Back Cushion £

QJR030130	SEDEO LITE, Back Cushion, Profiled, Height 500 mm, Comfort Cover	210
QJR030131	SEDEO LITE, Back Cushion, Profiled, Adjustable Straps, Height 570 mm, Comfort Cover	295
QJR030132	SEDEO LITE, Back Cushion, Deep Contour, Passive, Adjustable Straps, Height 570 mm, Comfort Cover	350
QJR030133	SEDEO LITE, Back Cushion, Deep Contour, Active, Adjustable Straps, Height 570 mm, Comfort Cover	350

SEDEO Lateral Thoracic Supports £

QJR030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	160
QJR030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	160

Headrest Pad £

QJR030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm	95
QJR030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm	110
QJR030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm	170

Headrest Mount £

QJR030433	SEDEO, Headrest Mount	140
WHT030409	Whitmyer, Headrest Mount, AXYS	145

Headrest Receiver £

QJR030120	Headrest Receiver, Backrest Frame Mount	0
-----------	---	---

Armrest Type  £

QJR040047	Armrest, Single Post, Height-Adjustable, Removeable	0
QJR040048	Arm Pad, Integral Foam, Depth-Adjustable, 380 mm	0
QJR040049	Side Guard, Depth-Adjustable	0

UNI Swing Away Legrest £

QJR050115	Legrest, UNI, 80°, Swing Away - Select Footplate Below	0
QJR050138	Legrest, UNI, Manual Elevating, Swing Away, Pair	275

Legrest Positioning

QJR050150	In-Rig Legrest 20 mm Per Side	55
QJR050151	Out-Rig Legrest 20 mm Per Side	55

Footplates


QJR050135	Footplate, Fixed	0
QJR050137	Footplate, Angle Adjustable, Pair	80
QJR050142	Amputee Supports, Padded, Right	150
QJR050143	Amputee Supports, Padded, Left	150

SEDEO LITE Swing Away Legrest £

QJR050100	SEDEO LITE, Legrest, 70°, Swing Away	0
QJR050105	SEDEO LITE, Legrest, Manual Elevating, Swing Away	375
QJR050017	SEDEO LITE, Legrest, Power Elevating, Swing Away	695

Legrest Positioning

QJR050150	In-Rig Legrest 20 mm Per Side	55
QJR050151	Out-Rig Legrest 20 mm Per Side	55
QJR050013	Amputee Supports, Padded, Right	195
QJR050014	Amputee Supports, Padded, Left	195

SEDEO PRO Swing Away Legrest		£
QJR050015	SEDEO PRO, Legrest, 70°, Swing Away	130
Support Options		
QJR050082	Calf Support, Integral Foam, Depth-Adjustable, Pair	100
QJR050083	Calf Support, Padded, Depth-Adjustable, Pair	135
QJR050084	Heel Loop	0
QJR050085	Calf Strap	0
QJR050200	Knee Pad, Integral Foam	75
Centre Mount Legrest (Q100 R Only)		£
QJR050148	Centre Mount Legrest, Angle Adjustable, With Flip Up Footplates (1)	390
Remote Joystick – VR2 		£
QJR110036	Remote Joystick, VR2	0
QJR110037	Remote Joystick, VR2, 2 Actuators - Required if retro fitting tilt	0
QJR110040	Remote Joystick, VR2, 2 Actuators And Lights	0
Power Modules – VR2		£
QJR110085	Power Module, VR2, 60 A	0
QJR110001	Power Module, VR2, 70 A - Required if retro fitting tilt (2)	0
QJR110002	Power Module, VR2, 90 A (3)	135
Lighting Module – VR2		£
QJR110030	Lighting Module, VR2	0
Remote Joystick – R-Net		£
QJR110015	Remote Joystick, R-Net, LED, Actuators	485
QJR110016	Remote Joystick, R-Net, LED, Actuators And Lights	485
QJR110017	Remote Joystick, R-Net, Colour LCD	600
QJR110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	750
QJR110009	Traction Control System For R-Net (GYRO Technology)	400
Power Modules – R-Net		£
QJR110003	Power Module, R-Net, 90 A	0
Joystick Mounting		£
QJR110082	Joystick Mounting, Parallel, Swing Away	0
Joystick Mounting Location		£
QJR110052	Joystick Mounting Location, Right	0
QJR110053	Joystick Mounting Location, Left	0
Attendant Controls - Select Joystick Above		£
QJR110021	Attendant Control, Dual	425
QJR110158	Attendant Only Tiller Control	390
QJR110020	Joystick Mounted As Attendant Control	270
Attendant Mounting Location		
QJR110025	Attendant Mounting Location, Right	0
QJR110026	Attendant Mounting Location, Left	0
LED Lights & Indicators		£
QJR110201	LED Light And Indicator StVZO Approved	400

(1) Not available with kerb climber

(2) Required For Actuator Control

(3) With Q200 R 10 kph Motors Only

Remote Joystick – LinX		£
QJR110810	REM110, Drive Only	0
QJR110811	REM211, Drive & Actuators	105
QJR110818	Omit Joystick	0
Power Modules & ISM's – LinX		£
QJR110802	Power Module, LinX, 75 A, Drive Only	0
QJR110803	Power Module, LinX, 75 A , Actuators and Lights	105
Joystick Mounting - LiNX		£
QJR110836	Joystick Mounting, Parallel, Swing Away	0
Joystick Mount Location - LiNX		£
QJR110838	Joystick Mounting Location, Right	0
QJR110839	Joystick Mounting Location, Left	0
Attendant Controls - LiNX		£
QJR110840	Dual Attendant Control, LiNX	425
QJR110841	Omit Dual Control	0
QJR110842	Joystick Mounted As Attendant Control	270
QJR110843	Omit Attendant Control	0
Attendant Mounting Location		
QJR110844	Attendant Mounting Location, Right	0
QJR110845	Attendant Mounting Location, Left	0
Additional Modules- LiNX		£
QJR110830	Traction Control System, for LiNX Gyro Technology (1)	400
QJR110831	Omit Gyro	0
QJR110820	2 Actuator Module LiNX, ACT-200 (1)	105
QJR110822	Omit Actuator Box	0
Kerb Climber		£
QJR090038	Kerb Climber (2)	200

(1) Q200 R Only
(2) Not with centre mount

Drive Wheels - Q200-R Only		£
QJR070084	Drive Wheel, 12", Solid, Knobbly, Black	0
QJR070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
QJR070033	Drive Wheel, 14", Solid, Knobbly, Black	0
Castor Wheels - Q200-R Only		£
QJR080210	Castor Wheel, 9", Solid, Black	0
QJR080211	Castor Wheel, 9", Pneumatic, Black	0
Motor Packages		£
QJR180014	2-Pole Motor, 4 kph / 2.5 mph, Via Programming	0
QJR180011	2-Pole Motor, 6 kph / 4 mph, Standard	0
QJR180006	2-Pole Motor, 10 kph / 6 mph (1)(2)	310
Colour Accents		£
QJR140004	Black (RAL 9005)	0
QJR140002	Red (RAL 040 40 67)	0
QJR140005	Apple Green (RAL 110 70 77)	0
QJR140003	Blue (RAL 260 40 45)	0
Batteries & Chargers		
Batteries		
QJR160012	Battery, 40 Ah, Gel, Maintenance-Free	0
QJR160003	Battery, 50 Ah, Gel, Maintenance-Free - Required if retro fitting tilt	0
QJR160007	Battery, 55 Ah, Maintenance-Free (3)	0
Chargers		
QJR160100	Charger, Automatic, 220 V / 8 A, UK	0
QJR160109	Omit Charger	0
Positioning Supports & Accessories		£
QJR090128	Positioning Belt, Buckle, Short, 50 mm x 1700 mm	0
QJR090110	Rear View Mirror, Right	60
QJR090118	Rear View Mirror, Left	60
QJR090009	Crutch Holder	70
QJR030605	Abduction Pommel (4)	255
QJR090124	Tray Table, Swing Away	335
DAHL Cover Kit		£
QJR090164	DAHL Cover Kit (5)	130
Switches		£
QJR110601	Single Switch, Ø 64 mm, Black	120
QJR110602	Single Switch, Ø 64 mm, Red	120
QJR110603	Single Switch, Ø 64 mm, Green	120
Switch Mounting		
QJR110560	Switch Mounting, Link-It To The Hand Control (Available Post Launch)	325
QJR110566	Switch Mounting, Link-It To the Armrest (Opposite Site As The Hand Control)	325

(1) Q200 R Only. Not Available With Tiller Control

(2) Requires 90Amp VR2 power module

(3) Not With Power Tilt

(4) Requires Seat Pan

(5) Backrest Cover Only - Purchase Full DAHL Kit Separately From DAHL. Not Available With Power Recline

Comments For Personalisations



Q100 R & Q200 R

£

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year



Sunrise Medical

Thorns Road, Brierley Hill, West Midlands
DY5 2LD

Tel.: +44 (0) 845 605 66 88

Fax: +44 (0) 845 605 66 89

E-Mail: enquiries@sunmed.co.uk

www.SunriseMedical.co.uk