

Q500[®] M/F/R

with SEDEO[®] PRO



ORDER FORM

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address	
C.-ID:			C.-ID:			
Name:			Name:			
Street:			Street:			
Town/PC:			Town / PC:			
E-Mail:			E-Mail:			
Phone:			Phone:			
Order Information			User Weight			
Order	Calculation	Built-4-Me [®]	User Weight Capacity: 160 kg			
Purchase Order No.:			User Weight: kg			
Art. No.: Q500 M: 2250050M-000, Q500 F: 2250050F-000, Q500 R: 2250050R-000						
Order Date:	Day	Month	Year			

Seating Specifications

All Bases Prepared For DAHL Docking Station		Seat Type	SEDEO PRO
Typical WC Weight	Without Batteries 140 kg	Seat Size Range	Width 380 - 560 mm, Depth 400 - 560 mm
	With Batteries 165 kg	Seat To Floor Height Range	415 mm - 490 mm

MWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 620 mm	Tires / Castors	7" Front / 7" Rear (Black)
Drive Control	VR2, Option R-Net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	70 mm	Frame Colour	Black

FWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 610 mm	Tires / Castors	9" Rear (Black), Options 10"
Drive Control	VR2, Option R-net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	65 mm	Frame Colour	Black

RWD Base Specifications

User Weight Capacity	160 kg	Drive Wheels	13" / 14" Drive Wheels (Black)
Base Size	Width Min. 610 mm	Tires / Castors	9" Front (Black), Options 10"
Drive Control	VR2, Option R-net	Charger	8 A Dual Mode Charger
Driving Range	Up To 40 km	Battery Size	60 And 80 Ah Batteries
Drive Motors	4-Pole 6 kph, Options 8, 10 And 12.5 kph	Colour Accents	7 Colour Accents
Ground Clearance	65 mm	Frame Colour	Black

Q500 Drive Bases			£
QOP000023	Q500 M (MWD Base Chair)		6350
QOP000024	Q500 R (RWD Base Chair)		6750
QOP000029	Q500 F (FWD Base Chair)		6750

Seat Cushion Type			£
QOP020110	Omit Seat Cushion		0

SEDEO Seat Size				£
Seat Cushion Width		Seat Depth		
SW38	15" (385 mm)	SD40	15.5" (400 mm)	
SW42	16" (420 mm)	SD42	16" (420 mm)	
SW46	18" (460 mm)	SD44	17" (440 mm)	
SW50	20" (500 mm)	SD46	18" (460 mm)	
SW54	21" (540 mm)	SD48	19" (480 mm)	
SW56	22" (560mm)	SD50	19.5" (500 mm)	
		SD52	20" (520 mm)	
		SD54	21" (540 mm)	
		SD56	22" (560 mm)	

Seat Rail			£
QOP020085	Seat Rail, Short, For Seat Depth 400 - 520 mm		0
QOP020069	Seat Rail, Long, For Seat Depth 460 - 560 mm (1)		0

SEDEO Seat Cushion			£
QOP020073	SEDEO, Seat Cushion, Profiled, Height 70 mm		139
QOP020076	SEDEO, Seat Cushion, Support Memory, Height 90 mm		231
QOP020078	SEDEO, Seat Cushion, Stabilizer, Height 90 mm		263

SEDEO Seat Cushion Covers			£
QOP020518	Seat Cushion Cover, Incontinence Protection		0
QOP020080	Seat Cushion Cover, Comfort		0

SEDEO Seat Boards			£
QOP020082	SEDEO, Seat Board, Flat, Depth-Adjustable, Cushions Guide At Rear		122
QOP020083	SEDEO, Seat Pan, 50 mm, Depth-Adjustable, Cushions Guide At Rear		139

JAY Seat Size				£
Seat Cushion Width		Seat Depth		
SW40	15.5" (400 mm)	SD40	15.5" (400 mm)	
SW42	16" (420 mm)	SD42	16" (420 mm)	
SW44	17" (440 mm)	SD44	17" (440 mm)	
SW46	18" (460 mm)	SD46	18" (460 mm)	
SW48	19" (480 mm)	SD48	19" (480 mm)	
SW50	19.5" (500 mm)	SD50	19.5" (500 mm)	
SW52	20" (520 mm)	SD52	20" (520 mm)	
SW54	21" (540 mm)	SD54	21" (540 mm) (2)	
SW56	22" (560 mm)	SD56	22" (560 mm) (2)	

(1) Not with seat width 380 & 400 mm
(2) Only with seat width 560 mm

JAY Seat Cushions		£
JAY020005	JAY Easy Fluid With Flat Base	372
JAY020006	JAY Easy Visco With Flat Base	160
JAY020010	JAY Balance, Air, Single Valve, Microclimatic Cover, Until Seat Width 500 mm	607
JAY020012	JAY Balance, Air, Single Valve, Microclimatic Cover, From Seat Width 560 mm	607
JAY020014	JAY Balance, Air, Single Valve, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	686
JAY020016	JAY Balance, Air, Single Valve, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	686
JAY020011	JAY Balance, Fluid, Microclimatic Cover, Until Seat Width 500 mm	565
JAY020013	JAY Balance, Fluid, Microclimatic Cover, From Seat Width 560 mm	565
JAY020015	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, Until Seat Width 500 mm	644
JAY020017	JAY Balance, Fluid, Microclimatic Cover, With Positioning Kit, From Seat Width 560 mm	644
JAY020018	JAY Xtreme Active	611
JAY020023	JAY Leg Support Kit For The Xtreme Active (1)	47
JAY020019	JAY J2 Cushion, Until Seat Width 460 mm Or Seat Depth 480 mm	594
JAY020021	JAY J2 Cushion, From Seat Width 480 mm Or Seat Depth 500 mm	594
JAY020020	JAY J2 Cushion, Deep Contour, Until Seat Width 480 mm	634
JAY020022	JAY J2 Cushion, Deep Contour, From Seat Width 500 mm	634

Seat Heights (2)		£
SH415	415 mm (3) (4)	0
SH425	425 mm (5) (9)	0
SH44	440 mm (4)	0
SH45	450 mm (5) (9)	0
SH465	465 mm	0
SH49	490 mm	0

Powered Seating Options		£
QOP130011	Fixed Tilt Module (Adjustable 0°, 3°, 6°, 9°, Factory Setting 6°)	0
QOP130004	Posterior Tilt, Powered, 0° - 30° (5)	1003
QOP130012	Posterior Tilt, Powered, -5° - 25° (5) (6)	1150
QOP130007	Posterior Tilt, Powered, 0° - 50° (5) (6)	1279
QOP130005	Posterior Tilt, Powered, -5° - 45° (5) (6)	1404
QOP130009	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 30°), Powered (5)	1801
QOP130006	Seat Elevator (300 mm) / Tilt Module (Tilt 0° - 50°), Powered (5) (8)	2153

Backrest Frame Type		£
QOP030083	Backrest, Tool Angle Adjustable, 90° - 120°	0
QOP030084	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Right (7)	289
QOP030085	Backrest, Manual Recline, User Adjustable, 90° - 130°, Release Lever, Left (7)	289
QOP030086	Backrest, Powered Recline, 90° - 130° (5)	520
QOP030087	Backrest, Powered Recline, 90° - 150° With 100 mm Mechanical Anti-Shear (5)	936

Backrest Frame Option		£
QOP030098	Backrest, Flip Down, Quick Release	173

Push Handles		£
QOP030204	Push Handles, Height-Adjustable	187

(1) Only with seat width 560 mm

(2) Measured from front of the seat pan to the floor

(3) Not with 80Ah batteries

(4) Not With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered

(5) Not With QOP110036 Remote Joystick, VR2. In Case Of R-Net: QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L

(6) Not With Seat Height 415 Or 440 mm

(7) Not On Same Joystick Side

(8) Only with MWD base

(9) Only With QOP130009 Seat Elevator/Tilt Module (Tilt 0° - 30°), Powered or QOP130006 Seat Elevator/Tilt Module (Tilt 0° - 50°), Powered

SEDEO Pro Back Size			
Backrest Width		Backrest Height	
BW38	15" (385 mm) - Effective width 38cm	BH44	18" (440 mm) - Effective Height 44 - 50cm
BW42	16" (420 mm) - Effective width 38 - 42cm	BH50	19.5" (500 mm) - Effective Height 50 - 56cm
BW46	18" (460 mm) - Effective width 42 - 46cm	BH56	22" (560 mm) Effective Height 56 - 62cm
BW50	19.5" (500 mm) - Effective width 46 - 50cm		
BW54	21" (540 mm) - Effective width 50 - 54cm		
BW56	22" (560 mm) - Effective width 52 - 56cm		

SEDEO Pro Backrest		£
QOP030109	SEDEO PRO, Back Cushion, Profiled	130
QOP030113	SEDEO PRO, Back Cushion, Profiled, Adjustable Straps	244
QOP030114	SEDEO PRO, Back Cushion, Deep Contour, Passive	312
QOP030115	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Passive, Adjustable Straps	405
QOP030116	SEDEO PRO, Back Cushion, Flexible, Deep Contour, Active, Adjustable Straps	394
QOP030117	SEDEO PRO, Back Cushion, Flexible, Extra Deep Contour, Active, Adjustable Straps	437
	SEDEO PRO, Back Cushion, Adaptive (Choose B4Me Requirements Separately)	

SEDEO Pro Backrest Covers		£
QOP030059	Back Cushion Cover, Incontinence Protection	0
QOP030108	Back Cushion Cover, Comfort	0

SEDEO Lateral Thoracic Supports		£
QOP030607	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Right	150
QOP030608	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Long Pad, 150 x 130 mm, Left	150
QOP030609	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Right	150
QOP030610	SEDEO, Lateral, Thoracic Support, Swing Away, Height-Adjustable, Short Pad, 90 x 120 mm, Left	150

Headrest Receiver		£
QOP030119	Headrest Receiver, Back Cushion Mount	0
QOP030120	Headrest Receiver, Backrest Frame Mount	0

Headrest Pad		£
QOP030430	SEDEO, Headrest Pad, Small, (LxH) 240 x 150 mm (1)	84
QOP030431	SEDEO, Headrest Pad, Small, With Cover, (LxH) 240 x 150 mm (1)	99
QOP030432	SEDEO, Headrest Pad, Flexible, (LxH) 400 x 110 mm (1)	150
WHT030401	Whitmyer, Headrest Pad, Plush, (LxH) 203 x 127 mm (2)	130
WHT030402	Whitmyer, Headrest Pad, Plush, (LxH) 254 x 127 mm (2)	144
WHT030403	Whitmyer, Headrest Pad, Plush, (LxH) 356 x 127 mm (2)	209
WHT030404	Whitmyer, Headrest Pad, Adjust-A-Plush, Small, (LxH) 355 x 76 mm (2)	234
WHT030406	Whitmyer, Headrest Pad, Adjust-A-Plush, Large, (LxH) 460 mm x 127 mm (2)	235

Headrest Stem (3)		£
QOP030433	SEDEO, Headrest Mount	126
WHT030408	Whitmyer, Headrest Mount, LINX ²	158
WHT030409	Whitmyer, Headrest Mount, AXYS	130
WHT030410	Whitmyer, Headrest Mount, AXYS, Offset	144
WHT030411	Whitmyer, Headrest Mount, COBRA XTRA Flip Back	195

(1) Only use with SEDEO headrest mount

(2) Only use with Whitmyer headrest mount

(3) Choose a Whitmyer headrest mount and headrest receiver in case you would like to add a head control in the special controls section

JAY J3 Back Size		£
JAY030111	JAY J3 Back Size, 410 mm	0
JAY030112	JAY J3 Back Size, 460 mm	0
JAY030113	JAY J3 Back Size, 510 mm	0
JAY J3 Back Size, Shallow Contour		£
JAY030012	JAY J3 Back, Shallow Contour, Upper Thoracic, Short, 420 mm	691
JAY030019	JAY J3 Back, Shallow Contour, Upper Thoracic, Medium, 460 mm	691
JAY030025	JAY J3 Back, Shallow Contour, Upper Thoracic, Tall, 510 mm	
JAY030015	JAY J3 Back, Shallow Contour, Shoulder Height, Short, 530 mm	691
JAY030022	JAY J3 Back, Shallow Contour, Shoulder Height, Medium, 570 mm	691
JAY030028	JAY J3 Back, Shallow Contour, Shoulder Height, Tall, 610 mm	691
JAY J3 Back Size, Mid Contour		£
JAY030013	JAY J3 Back, Mid Contour, Upper Thoracic, Short, 420 mm	691
JAY030020	JAY J3 Back, Mid Contour, Upper Thoracic, Medium, 460 mm	691
JAY030026	JAY J3 Back, Mid Contour, Upper Thoracic, Tall, 510 mm	691
JAY030016	JAY J3 Back, Mid Contour, Shoulder Height, Short, 530 mm	691
JAY030023	JAY J3 Back, Mid Contour, Shoulder Height, Medium, 570 mm	691
JAY030029	JAY J3 Back, Mid Contour, Shoulder Height, Tall, 610 mm	691
JAY J3 Back Size, Mid Deep Contour		£
JAY030045	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Short, 420 mm	757
JAY030042	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Medium, 460 mm	757
JAY030046	JAY J3 Back, Mid Deep Contour, Upper Thoracic, Tall, 500 mm	757
JAY030043	JAY J3 Back, Mid Deep Contour, Shoulder Height, Short, 530 mm	757
JAY030047	JAY J3 Back, Mid Deep Contour, Shoulder Height, Medium, 570 mm	757
JAY030048	JAY J3 Back, Mid Deep Contour, Shoulder Height, Tall, 610 mm	757
JAY030050	JAY J3 Back, Mid Deep Contour, Extra Deep, Upper Thoracic, Medium, 460 mm	757
JAY J3 Back Size, Deep Contour		£
JAY030014	JAY J3 Back, Deep Contour, Upper Thoracic, Short, 420 mm	757
JAY030021	JAY J3 Back, Deep Contour, Upper Thoracic, Medium, 460 mm	757
JAY030027	JAY J3 Back, Deep Contour, Upper Thoracic, Tall, 500 mm	757
JAY030017	JAY J3 Back, Deep Contour, Shoulder Height, Short, 530 mm	757
JAY030024	JAY J3 Back, Deep Contour, Shoulder Height, Medium, 570 mm	757
JAY030030	JAY J3 Back, Deep Contour, Shoulder Height, Tall, 610 mm	757
JAY030051	JAY J3 Back, Deep Contour, Extra Deep, Upper Thoracic, Medium, 460 mm	757
JAY J3 Back, Spine Align Fit Kit		£
Spine Align Mini Fit Kit I		
JAY090059	JAY J3 Spine Align Mini Fit Kit I, 360 mm	38
JAY090029	JAY J3 Spine Align Mini Fit Kit I, 410 mm	38
JAY090030	JAY J3 Spine Align Mini Fit Kit I, 460 mm	38
JAY090031	JAY J3 Spine Align Mini Fit Kit I, 510 mm	38
Spine Align Mini Fit Kit II		
JAY090060	JAY J3 Spine Align Mini Fit Kit II, 360 mm	57
JAY090032	JAY J3 Spine Align Mini Fit Kit II, 410 mm	57
JAY090033	JAY J3 Spine Align Mini Fit Kit II, 460 mm	57
JAY090034	JAY J3 Spine Align Mini Fit Kit II, 510 mm	57
Positioning Boxes		
JAY090035	JAY J3 Positioning Box For Back	183
JAY090036	JAY J3 Hardware Box For Back	572

JAY J3 Back, Chest Strap (1)		£
JAY090040	JAY J3 Chest Strap, Centre Open Harness, 4 Point Strap, Medium	190
JAY090041	JAY J3 Chest Strap, Centre Open Harness, 4 Point Strap, Large	190
JAY090042	JAY J3 Chest Strap, Centre Open Harness, 4 Point Strap, Extra Large	190

JAY J3 Lateral Thoracic Pads, Right		£
JAY090020	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 100 mm, Right	56
JAY090021	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 127 mm, Right	56
JAY090022	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 152 mm, Right	56
JAY090023	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 100 mm, Right	56
JAY090024	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 127 mm, Right	56
JAY090025	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 152 mm, Right	56

JAY J3 Lateral Thoracic Pads, Left		£
JAY090011	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 100 mm, Left	56
JAY090012	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 127 mm, Left	56
JAY090013	JAY J3 Lateral Thoracic Pad, (LxH) 100 x 152 mm, Left	56
JAY090014	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 100 mm, Left	56
JAY090015	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 127 mm, Left	56
JAY090016	JAY J3 Lateral Thoracic Pad, (LxH) 152 x 152 mm, Left	56

JAY J3 Lateral Thoracic Brackets		£
JAY090026	JAY J3 Lateral Thoracic, Fixed Bracket, Right	56
JAY090017	JAY J3 Lateral Thoracic, Fixed Bracket, Left	56
JAY090028	JAY J3 Lateral Thoracic, Bracket, Swing Away, Right	105
JAY090019	JAY J3 Lateral Thoracic, Bracket, Swing Away, Left	105

(1) Only available with JAY Bracket

Armrest Types & Side Gaurds		£
QOP040019	Armrest, Single Post, Flip Back (1)	0
QOP040021	Armrest, Single Post, Removable (1)	0
QOP040123	Armrest, Reclining, Flip Back	173

Side Guard for Reclining Flip Back Armrest

QOP040036	Side Guard for reclining flip back arm (2)	0
-----------	--	---

Side Gaurds for Single Post Flip Back / Removable

QOP040033	Side Guard (3)	0
QOP040034	Side Guard, Integral Foam (4)	0
QOP040035	Side Guard, Wedge Shaped, Padded (5)	262

Arm Pads		£
QOP040023	Arm Pad, Integral Foam, 300 mm, Right	0
QOP040024	Arm Pad, Integral Foam, 300 mm, Left	0
QOP040025	Arm Pad, Integral Foam, 400 mm, Right	0
QOP040026	Arm Pad, Integral Foam, 400 mm, Left	0
QOP040027	Arm Pad, Integral Foam, 450 mm, Right (6)	0
QOP040028	Arm Pad, Integral Foam, 450 mm, Left (6)	0
QOP040204	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Right	69
QOP040203	Arm Pad, Comfort, Padded, PU Coated, 400 mm, Left	69
QOP040031	Arm Pad, Hemiplegia, Integral Foam, Right (7)	254
QOP040032	Arm Pad, Hemiplegia, Integral Foam, Left (8)	254

Thigh Supports & Pommel		£
QOP020342	Knee Support, Lateral, 120 x 120 mm, Swing Away, Right (9)	173
QOP020341	Knee Support, Lateral, 120 x 120 mm, Swing Away, Left (9)	173
QOP020113	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Right (10)	173
QOP020114	Thigh Support, Lateral, Small, 120 x 120 mm, Swing Away, Left (10)	173
QOP020121	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Right (10)	173
QOP020122	Thigh Support, Lateral, Large, 180 x 120 mm, Swing Away, Left (10)	173
QOP030605	Abduction Pommel (9)	237

Centre Mount Legrest (Single Piece Footplate)		£
QOP050086	Centre Mount Legrest, Tool Adjustable, With Flip Up And Angle-Adjustable Footboard (11)	0

Pad Options

QOP050087	Centre Calf Support, Integral Foam	0
QOP050092	Calf Pad, Integral Foam, Pair	0

Fixed Centre Mount Legrest (Split Footplate)		£
QOP050140	Fixed Centre Mount Legrest	116
QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
QOP050075	Calf Pad, Padded, Pair	0

(1) Not With Seat To Floor Height 415 mm

(2) Only With QOP040123 Armrest, Reclining, Flip Back

(3) Not With QOP040123 Armrest, Reclining, Flip Back

(4) Seat Width -20 mm Not With QOP040123 Armrest, Reclining, Flip Back

(5) Seat Width -80 mm, Only With QOP040021 Armrest, Single Post, Removable

(6) With QOP040019 Flip Back Armrest, Seat Width +20 mm

(7) Not With QOP110052 Joystick Mount, Right, QOP090040 Tray Table, Swing Away,

QOP090046 Tray Table, Half-Side, Swing Away, QOP090110 Rear View Mirror, Right

(8) Not With QOP110053 Joystick Mount, Left, QOP090040 Tray Table, Swing Away, QOP090046

Tray Table, Half-Side, Swing Away, QOP090118 Rear View Mirror, Left

(9) Not With QOP050015-17 SEDEO PRO, Legrest, Swing Away, QOP050080 SEDEO PRO, Hanger Receivers Only, QOP090122 Transfer Bar

(10) Not With Seat To Floor Height 415 mm, Not With QOP040035 Side Guard, Wedge Shaped, Padded

(11) Not With (Power) Centre Mount Legrests

(12) In Combination With Lift Or Tilt, Choose QOP050087 Centre Calf Support, Integral Foam Or

QOP050092 Calf Pad, Integral Foam, Pair

(13) Not With RWD

Powered Centre Mount		£
QOP050076	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension	1372
QOP050071	Centre Mount Legrest, Powered, 90° - 180°, 175 mm Extension, Footplate To Floor	1621
Footplates		
QOP050072	Footplate, Dual, Height- And Angle-Adjustable, 260 x 150 mm	0
QOP050073	Footplate, Dual, Height- And Angle-Adjustable, Large, 310 x 170 mm	0
Pad Options		
QOP050087	Central Calf Support, Integral Foam	0
QOP050092	Calf Pad, Integral Foam, Pair	0
QOP050075	Calf Pad, Padded, Pair	0

SEDEO Pro Swing Away Legrest		£
QOP050015	SEDEO PRO, Legrest, 70°, Swing Away	0
QOP050016	SEDEO PRO, Legrest, Manual Elevating, Swing Away	335
QOP050017	SEDEO PRO, Legrest, Power Elevating, Swing Away	624
Footplates		
QOP050088	Footplate, For Seat Width 385 - 400 mm	0
QOP050089	Footplate, For Seat Width 420 - 460 mm	0
QOP050090	Footplate, For Seat Width 480 - 520 mm	0
QOP050091	Footplate, For Seat Width 540 - 560 mm	0
QOP050081	Footplates, Coupled (1) (2)	0
Calf Pad / Strap		
QOP050084	Heel Loop	0
QOP050085	Calf Strap (2)	0
QOP050082	Calf Support, Integral Foam, Depth-Adjustable, Pair	104
QOP050083	Calf Support, Padded, Depth-Adjustable, Pair	128
QOP050200	Knee Pad, Integral Foam	68
Amputee Support (3)		
QOP050013	Amputee Support, Padded, Right	134
QOP050014	Amputee Support, Padded, Left	134

VR2 Controls		£
QOP110036	Remote Joystick, VR2	0
QOP110037	Remote Joystick, VR2, 2 Actuators	0
QOP110040	Remote Joystick, VR2, 2 Actuators And Lights	0
VR2 Power Module & Lighting Module		
QOP110002	Power Module, VR2, 90 A	0
QOP110030	Lighting Module, VR2 (4)	0

R-Net Controls		£
QOP110015	Remote Joystick, R-Net, LED, Actuators	434
QOP110016	Remote Joystick, R-Net, LED, Actuators And Lights	434
QOP110017	Remote Joystick, R-Net, Colour LCD	561
QOP110018	Remote Joystick, R-Net, Colour LCD, Heavy Duty	1120
QOP110013	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared Control Included	688
QOP110023	Remote Joystick, R-Net, High Resolution Colour LCD, Advanced, Infrared And Bluetooth Control Included	866
R-Net Power Module & ISM		
QOP110003	Power Module, R-Net, 90 A	0
QOP110004	Power Module, R-Net, 120 A	122
QOP110032	Actuator Control & Lights Module, R-Net, ISM-6L	0
R-Net Gyro		
QOP110009	Traction Control System For R-Net (GYRO Technology)	376

(1) Choose Also A Footplate

(2) Only For QOP050015 SEDEO PRO, Legrest, 70°, Swing Away

(3) Extra On A Set Of Swing Away Legrests, Attached Lose

(4) If LED Light And Indicator Is Chosen, Then This Is Required

Joystick Mount		£
QOP110050	Joystick Mounting, Fixed	0
QOP110056	Joystick Mounting, Parallel, Swing Away, Capsulated	74
Joystick Mounting Location		£
QOP110052	Joystick Mounting Location, Right	0
QOP110053	Joystick Mounting Location, Left	0
Attendant Controls		£
QOP110021	Dual Control	399
QOP110020	Joystick Mounted As Attendant Control (1)	243
Attendant Mounting Location		
QOP110025	Attendant Mounting Location, Right	0
QOP110026	Attendant Mounting Location, Left	0
Table Top Controls		£
QOP110143	Table Top Control, LCD, Centre Tray Mounted (3) (4) (5)	1996
Tray Table Mounting Location		
QOP110080	Tray Table Mounting Location, Left	0
QOP110081	Tray Table Mounting Location, Right	0
Joystick Knobs		£
QOP090505	Joystick Knob, Golf Ball, Ø 42 mm	37
QOP090506	Joystick Knob, T-Bar, 82 mm Wide	37
QOP090507	Joystick Knob, Soft Ball, Ø 70 mm	37
QOP090508	Joystick Knob, Ball, Ø 35 mm	37
QOP090045	Joystick Knob, Allround, Ball, Ø 48 mm	37
QOP090509	Joystick Knob, Chin Cup, Ø 37 mm	37
QOP090512	Joystick Knob, Stick, 80 x 15 mm	37
QOP090511	Joystick Knob, Mushroom, Ø 45 mm	37
LED Lights & Indicators		£
QOP110201	LED Light And Indicator StVZO Approved (2)	362
Drive Wheels		£
MWD, FWD & RWD Base		
QOP070065	Drive Wheel, 13", Pneumatic, Knobbly, Black	0
QOP070069	Drive Wheel, 13", Solid, Knobbly, Black	0
QOP070032	Drive Wheel, 14", Pneumatic, V-Type, Black	
QOP070031	Drive Wheel, 14", Solid, V-Type, Black	0
QOP070034	Drive Wheel, 14", Pneumatic, Knobbly, Black	0
QOP070033	Drive Wheel, 14", Solid, Knobbly, Black	

(1) Choose Also Remote Joystick
(2) QOP110030 Lighting Module VR2 Or QOP110032 Actuator Control & Lights Module, R-Net, ISM-6L Required
(3) Not with SW40, SW41, SW42, SW48, SW485

(4)Not with QOP110309 Hand Control, Allround Joystick, QOP110310 Hand Control, Allround Joystick Light, QOP110311 Hand Control, Allround Joystick Heavy Duty
(5) Only R-Net Controls, comes with joystick

Castor Wheels		£
MWD Base		
QOP080213	Castor Wheel, 7", Solid, Black	0
FWD/RWD		
QOP080210	Castor Wheel, 9", Solid, Black	0
QOP080211	Castor Wheel, 9", Pneumatic, Black	0
FWD/RWD		
QOP080210	Castor Wheel, 9", Solid, Black	
QOP080211	Castor Wheel, 9", Pneumatic, Black	
QOP080058	Castor Wheel, 10", Solid, Black	0
QOP080212	Castor Wheel, 10", Pneumatic, Black	0

Mudguards (FWD/RWD only)		£
QOP010007	Mudguard, Drive Wheel, 14"	112
QOP010008	Mudguard, Castor Wheel, 9"	88

Kerb Climber (RWD Only)		£
QOP090038	Kerb Climber	215

Suspension Setting		£
QOP150001	Standard Suspension Setting - Selected as default	0
QOP150002	Firm Suspension Setting (2) (3)	0

Motor Packages		£
QOP180012	2-Pole Motor, 6 kph/ 4 mph, Performance	0
QOP180009	4-Pole Motor, 6 kph / 4 mph	289
QOP180007	4-Pole Motor, 10 kph / 6.5 mph	363
QOP180008	4-Pole Motor, 13 kph / 8 mph	468

C-Me - Drive at 4.8kph in full Lift (4)		
QPL110034	C-Me Function, Enabled (5)	
QPL110035	C-Me Function, Disabled (6)	

Colour Accents		£
QOP140005	Apple Green (RAL 110 70 77)	0
QOP140004	Black (RAL 9005)	0
QOP140003	Blue (RAL 260 40 45)	0
QOP140006	Orange (RAL 060 60 80)	0
QOP140002	Red (RAL 040 40 67)	0
QOP140007	Silver (RAL 9006)	0
QOP140001	White (RAL 9010)	0

Batteries		£
QOP160003	Battery, 50 Ah, Gel, Maintenance-Free	0
QOP160001	Battery, 60 Ah, Gel, Maintenance-Free	0
QOP160005	Battery, 80 Ah, Gel, Maintenance-Free (1)	150

Chargers		£
QOP160100	Charger, Automatic, 220 V / 8 A, UK	0
QOP160103	Charger, Automatic, 220 V / 10 A, UK	214

(1) Requires QOP160103 Battery Charger 220 V / 10 A
(2) QOP150002 (Firm Suspension) recommended with QOP130009 Powered Lift/Tilt
(3) QOP150002 (Firm Suspension) recommended with QOP130004/QOP130012/QOP130007
Powered Tilt for clients over 75kg user weight

Lap Belt		£
QOP090128	Positioning Belt, Carseat Buckle, Short, 50 mm x 1700 mm	0
QOP090129	Positioning Belt, Carseat Buckle, Long, 50 mm x 2700 mm	0
Accessories		£
QOP090110	Rear View Mirror, Right (1)	53
QOP090118	Rear View Mirror, Left (1)	53
QOP090009	Crutch Holder (2)	62
QOP090119	Brackets For Bag On Backrest Frame (2)	32
QOP090120	Bag For Armrest (3)	100
QOP090122	Transfer Bar (4)	237
Tray Table (Mounted Opposite Joystick)		
QOP090040	Tray Table, Width-Adjustable, Swing Away	312
QOP090046	Tray Table, Half-Side, Swing Away (5)	249
Tray Table Accessories		
QOP090125	Partial Padding For Tray Table	19
QOP090126	Raised Edges, Acrylic For Tray Table	25
QOP090127	Elbow Pads, Padded, Pair For Tray Table	37
5-Switch Box		£
QOP110043	#1 Ctrl+5 Switch Box With Free Assignable Functions	809
#1 Ctrl+5 Switch Box Mounting Location		
QOP110180	#1 Ctrl+5 Switch Box Mounting Location, Behind Joystick (6)	0
QOP110181	#1 Ctrl+5 Switch Box Mounting Location, Opposite Joystick	0
#1 Ctrl+5 Switch Type		
QOP110175	#1 Ctrl+5 Switch Type With Push Buttons	0
QOP110176	#1 Ctrl+5 Switch Type With Toggle Switches	0
Omni2 Input Control Module		£
QOP110116	OMNI2 Input Control Module, R-Net, Bluetooth, Infrared (7)	1093
QOP110117	OMNI2 Input Control Module, R-Net, Scanning Mode, Bluetooth, Infrared	1093
Omni 2 Module Mount		£
QPO110119	OMNI2 Module, Fixed Mount, Right	0
QPO110120	OMNI2 Module, Fixed Mount, Left	0
QPO110123	OMNI2 Module, Swing Away Mount, Right	116
QPO110124	OMNI2 Module, Swing Away Mount, Left	116
Hand Control		£
QOP110301	Hand Control, Switch-It MicroGuide	2807
QOP110302	Hand Control, Switch-It MicroPilot	3119
QOP110304	Hand Control, Switch-It VersaGuide	998
QOP110305	Hand Control, Switch-It VersaGuide EZ (8)	1185
QOP110306	Hand Control, Allround Joystick	998
QOP110307	Hand Control, Allround Joystick Light (9)	1185
QOP110308	Hand Control, Allround Joystick Heavy Duty	1871
QOP110309	Hand Control, Allround Joystick, R-Net Bus Enabled (10)	1061
QOP110310	Hand Control, Allround Joystick Light, R-Net Bus Enabled (10)	1247
QOP110311	Hand Control, Allround Joystick Heavy Duty, R-Net Bus Enabled (10)	1996

(1) Not In Combination With Same Entry Side QOP090040 Tray Table, Swing Away, QOP090046 Tray Table, Half-Side, SA
(2) Not With QOP030086 Powered Recline, QOP030087 Powered Recline, 90° - 150° With Anti-Shear
(3) Not With QOP040123 Armrest, Reclining, Flip Back
(4) Not With QOP020342/QOP020341 Knee Support, Lateral, Swing Away, Right/Left & QOP050015-17 SEDEO PRO, Legrest, Swing Away & QOP050080 SEDEO PRO, Hanger Receivers Only
(5) Not With QOP090125-127 Tray Table Accessories

(6) Not Possible With QOP110056 Joystick Mounting, Parallel Swing Away, Capsulated (Here We Can't Swing Away)
(7) Display With Built In Infrared & Armrest Mount, Required For Specialty Input Devices
(8) Not With QOP090512 Joystick Knob, Stick, 80 x 15 mm, QOP090511 Joystick Knob, Mushroom, Ø 45 mm, QOP090506 Joystick Knob, T-Bar, 82 mm Wide
(9) Only Use With Default Joystick Knob Or Joystick Knob, Allround, Ball Ø 48 mm
(10) OMNI Module Not Required

Hand Control Mounting		£
QOP110342	Hand Control Mounting, Link-It To The Space Disk, Hand Tray (6)	624
Hand Control Mounting		
QOP110350	Hand Control Mounting, Link-It To The Armrest, Fixed	185
QOP110353	Hand Control Mounting, Link-It To The Armrest, Swing Away	196
Hand Control Mounting Location		
QOP110362	Hand Control Mounting Location, Right	0
QOP110361	Hand Control Mounting Location, Left	0
Head Control (1)		£
Proportional Input Devices		
QOP110370	Head Control, Switch-It Dual Pro, 3-Way Proportional, Facial Pads	5363
QOP110371	Head Control, Switch-It Dual Pro, 3-Way Proportional, Spot Pads	5363
Chin Control		£
QOP110470	Chin Control, Switch-It MicroGuide	2807
QOP110471	Chin Control, Switch-It MicroPilot	3119
QOP110472	Chin Control, Switch-It VersaGuide	998
QOP110473	Chin Control, Switch-It VersaGuide EZ	1185
QOP110474	Chin Control, Allround Joystick	998
QOP110475	Chin Control, Allround Joystick Light	1185
QOP110476	Chin Control, Allround Joystick, R-Net Bus Enabled (2)	1061
QOP110477	Chin Control, Allround Joystick Light, R-Net Bus Enabled (2)	1247
Chin Control Mounting		£
QOP110480	Chin Control Mounting, Link-It To Chin Harness, Small (3)	562
QOP110481	Chin Control Mounting, Link-It To Chin Harness, Large (3)	629
QOP110483	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Right	499
QOP110484	Chin Control Mounting, Link-It To Manual Chin Boom, Backrest Mount, Left	499
QOP110485	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Right	2495
QOP110486	Chin Control Mounting, Link-It To Power Chin Boom, Backrest Mount, Left	2495
QOP110487	Chin Control Mounting, Link-It To Headrest, Swing Away, Left	437
QOP110488	Chin Control Mounting, Link-It To Headrest, Swing Away, Right	437
Link-It Extension Kit		
QOP110590	Link-It Extension Kit, Base Module, 100 mm Extension Tube, 5 Cable Clips	150
QOP110591	Link-It Extension Kit, Base Module, 150 mm Extension Tube, 5 Cable Clips	187
Sip ,N Puff (4)		£
QOP110510	Sip 'N Puff With Link-It Hardware, Swing Away Mount To The Headrest (5)	624

(1) For switch selection see section H. Choose a Whitmyer headrest mount.
(2) Omni 2 module required
(3) If you need switches mounted on the harness, you can only use All Round joystick
(4) For switch selection see section H
(5) Headrest to be selected
(6) Only For Micro Series Joysticks

#1 Switches - Section E (1)		£
QOP110611	#1 Single Switch, Ø 64 mm, Black	108
QOP110612	#1 Single Switch, Ø 64 mm, Red	108
QOP110613	#1 Single Switch, Ø 64 mm, Green	125
QOP110523	#1 Egg Switch – Black	119
QOP110524	#1 MicroLite Switch – Black	174
QOP110651	#1 Single Switch, Ø 30 mm, Black	108
QOP110654	#1 Single Switch, Ø 30 mm, Blue	108
QOP110653	#1 Single Switch, Ø 30 mm, Green	108
QOP110652	#1 Single Switch, Ø 30 mm, Red	108
QOP110655	#1 Single Switch, Ø 30 mm, Yellow	108
QOP110534	#1 Twister Switch – Black	174
QOP110535	#1 Twister Switch – Green	174
QOP110536	#1 Twister Switch – Red	174
QOP110537	#1 Satellite Switch – Black (7)	187
QOP110538	#1 Satellite Switch – Green (7)	187
QOP110539	#1 Satellite Switch – Red (7)	187

#1 Switch Mounting		£
QOP110560	#1 Switch Mounting, Link-It To The Hand Control (2)	187
QOP110566	#1 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (3)	187
QOP110561	#1 Switch Mounting, Link-It To The Headrest	250
QOP110563	#1 Switch Mounting, Link-It To The Seat Rail (4)	187
QOP110565	#1 Switch Mounting, Link-It To Chin Boom (5)	187
QOP110569	#1 Omit Switch Mounting (6)	0

#1 Switch Mount Location		£
QOP110572	#1 Switch Mounting Location, Right	0
QOP110571	#1 Switch Mounting Location, Left	0

#2 Switches - Section F (1)		£
QOP110631	#2 Single Switch, Ø 64 mm, Black	108
QOP110632	#2 Single Switch, Ø 64 mm, Red	108
QOP110633	#2 Single Switch, Ø 64 mm, Green	108
QOP110088	#2 Egg Switch – Black	125
QOP110089	#2 MicroLite Switch – Black	174
QOP110671	#2 Single Switch, Ø 30 mm, Black	108
QOP110674	#2 Single Switch, Ø 30 mm, Blue	108
QOP110673	#2 Single Switch, Ø 30 mm, Green	108
QOP110672	#2 Single Switch, Ø 30 mm, Red	108
QOP110675	#2 Single Switch, Ø 30 mm, Yellow	108
QOP110095	#2 Twister Switch – Black	174
QOP110096	#2 Twister Switch – Green	174
QOP110097	#2 Twister Switch – Red	174
QOP110111	#2 Satellite Switch – Black (7)	187
QOP110112	#2 Satellite Switch – Green (7)	187
QOP110113	#2 Satellite Switch – Red (7)	187

(1) Select 1 Switch And Mounting Position Accordingly. For A Second Switch Please Use #2 Switch Selection Part Below. Single Function For Special Input Devices
(2) Not With HD Allround Joystick Or Space Disk Mount
(3) Armrest Rail Mount

(4) Knee, If 2 Switches Selected Only Left And Right Mount
(5) No Left / Right Switch Selection Required
(6) Switches Attached Lose
(7) Only For Allround Joystick, Maximum 2 Selectable

#2 Switch Mounting		£
QOP110098	#2 Switch Mounting, Link-It To The Hand Control (1)	187
QOP110099	#2 Switch Mounting, Link-It To The Armrest (Opposite Site As The Hand Control) (2)	187
QOP110100	#2 Switch Mounting, Link-It To The Headrest	250
QOP110104	#2 Switch Mounting, Link-It To The Seat Rail (3)	187
QOP110105	#2 Switch Mounting, Link-It To Chin Boom (4)	187
QOP110106	#2 Omit Switch Mounting (5)	0

#2 Switch Mount Location		£
QOP110108	#2 Switch Mounting Location, Right	0
QOP110107	#2 Switch Mounting Location, Left	0

Multi-Function Switches (7)		£
QOP110550	Mini 2 Button (Stereo Plug) (8)	343
QOP110551	Mini 2 Button (Dual Mono Plug) (8)	343
QOP110552	Mini 4 Button InLine (DB9 Connector) (8)	576
QOP110554	4-Way Toggle Switched Joystick	576

Accessories		£
QOP110580	Mono Extension Cable (3.5 mm Jack Socket) (2 m)	125
QOP110581	Stereo Extension Cable (3.5 mm Jack Socket) (2 m)	162
QOP110582	Extension Cable DB9 (1 m)	150
QOP110583	Y-Cable, 2 Female Mono To 1 Male Stereo	187

Safety Devices		£
QOP090147	Remote Stop Switch (30 m Activation Range, Fit Via OBC Port) (9)	916

Comments For Personalisations



Q500 M with SEDEO PRO		£
Base Price		
Total Of Additional Equipment		
Packaging And Transport		
Exclusive VAT		
VAT		
Total Price		

Date	Signature
------	-----------

Day Month Year

- (1) Not With HD Allround Joystick Or Space Disk Mount
 (2) Armrest Rail Mount
 (3) Knee, If 2 Switches Selected Only Left And Right Mount
 (4) No Left / Right Switch Selection Required
 (5) Switches Attached Lose

- (6) Only For Allround Joystick, Maximum 2 Selectable
 (7) Used In Combination With Assignable Button Programming
 (8) Attached Lose, When Ordered. No Link-It Fixation Possible
 (9) Not With VR2