

# Q6 edge<sup>®</sup> 2.0



The Q6 Edge<sup>®</sup> 2.0 features versatile powered solutions, which can be configured with the highly adaptable TRU-Balance<sup>®</sup> 3 Seating System. Meeting the needs of the most active user.

## UNIQUE TECHNOLOGY

to the industry.

12" iLevel<sup>™</sup> | TB3<sup>®</sup> Seating System  
Q-Logic 3 with Bluetooth | TB3 Flip-up armrests

**EXTRA**  
STABILITY TECHNOLOGY<sup>™</sup>



CRASH  
TESTED



@QUANTUMREHABUK



Tel: 01869 324600 | Email: [sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk) | Web: [www.quantumrehab.co.uk](http://www.quantumrehab.co.uk)

(Version 6: August 2018)



## ATX® adjustable suspension

Provides increased ride quality, obstacle climbing and stability  
Semi-independent articulating rear castor beam.

## 55 Ah Batteries

More miles per charge! 55 Ah batteries allow for increased range, for the most active user. Easy access through front shroud.

## Improved Aesthetics

Updated front & rear shroud design.  
New "snap on" plastic chassis shrouds

## Turning Radius

Access smaller spaces easier with a tight turning radius of 20.5" (52.07 cm)

## Motor Package

2-Pole advanced responsive motors  
4 Mph or 6 Mph motor options available

## Measurements

14" Drive wheel, 6" low impact solid castor wheels front and rear.  
Length: 35.5" (Without legrests) | Width: 24.4" |  
Ground Clearance: 2.7"

## Seating Solutions



### Synergy Seating

- 10" x 10" to 22" x 22" seat size options
- Visco padded or tension adjustable backrest options
- Durable and easily adjustable
- Easy centre of gravity adjustment
- Accepts a wide variety of after market seats and backs



### TRU-Balance® 3 Seating

- 12" x 12" to 22" x 22" seat size options
- Make fine adjustments while client is in the chair
- Accepts a wide variety of after market seats and backs
- Endomorph and mesomorph seating, allows for independent adjustment of seat width and backrest width
- Premium accessories available

## Power Options



### TRU-Balance® 3 Power Tilt

- 0-50° of tilt (Anterior tilt available)
- Lower standard seat-to-floor height
- Easy adjustment and durability
- Maintains centre of gravity shift



### TRU-Balance® 3 Power Recline

- 85° - 175° of recline
- Adjustable shear reduction
- Compatible with after-market backrest interface

### Powered Elevating leg rests

- Articulating foot platform
- Articulating swing-away leg rests
- Can be operated in conjunction with TB3 recline



@sophieslifeonilevel

### 12" iLevel™ Power Adjustable Seat Lift

- Lower standard seat-to-floor height
- iLevel™ allows 12" of seat elevation while at walking speed (3.5 mph.)
- iLevel is retrofittable on the Q6 Edge 2.0



## Additional Features



### TRU-Balance® 3 Cantilever Flip-up armrests

- Quick height adjustment
- Can be internally or externally rotated for optimal positioning
- Multiple armrest pad options



### TRU-Balance® 3 Centre mount foot platform

- Fixed or powered articulating
- Modular footplate extension kits
- Multi-adjustable calf pads
- Paediatric & adult footplate options



### Swing-away legrests

- 60° - 90° standard swing-away
- Paediatric & adult footplate options
- Powered Articulating
- Manual elevating



### Attendant/Dual control options

- Full proportional joystick
- Quick and easy to use control interface
- Control the chair at the attendants speed



### Body positioning

- A full range of multi adjustable swing-away laterals, offering advanced positioning of the trunk
- Designed to easily allow head to toe positioning using any of the Stealth positioning accessories



### After-market backrest interface

- Allows for an unprecedented level of compatibility between TRU-Balance 3 back rest and most after-market backrests
- Retains full functionality, including shear reduction



### Premium accessories

- Transfer bars
- Mobile phone holder
- USB Port (available with Q-Logic 2)
- Clothes guard
- Additional accessories available



### Full range of control devices available

- Switches
- Sensors
- Proportionals
- i-Drive™





## Quantum Q6 Edge 2.0 | Colour Options & Tech Spec

### Colour Options



Red Corvette



Back in Black



Pretty in Pink



Lime Light



Silver Bullet



Orange Crush



Electric Blue



White Satin



Deep Purple



Grey Street



Midnight Blue



Yellow Blaze

### Technical Specifications

<b>Drive Wheels (cm)</b>  35.56 cm pneumatic or solid	<b>Caster Wheels (cm)</b>  Front and Rear: 15.24 cm	<b>Max Speed (mph)</b>  4mph or 6 mph with lights <sup>1</sup>	<b>Ground Clearance (cm)</b>  Battery Tray: 6.9 cm	<b>Turning Radius (cm)</b>  52.1 cm	<b>Overall Length (cm)</b>  90.2 cm	<b>Overall Width (cm)</b>  62.2 cm	<b>Seating Options (inches)</b>  Synergy and TB3 Seating	<b>Batteries</b>  x 2 50 - 55 AH Deep Cycle	<b>Available Electronics</b>  VR2/RNet NE+/QLogic 2	<b>Weight Capacity (stone)</b>  21 Stone 6 lbs (137kg)
<b>Seat to ground height (")</b>  15.9"-19.5"	<b>Base Weight (kg)</b>  63 kg	<b>Battery Weight Each (kg)</b>  17.24 kg	<b>Warranty</b> <ol style="list-style-type: none"> <li>1. Lifetime limited warranty on frame.</li> <li>2. 1-year limited warranty on electronics.</li> <li>3. 18-month warranty on drive control system.</li> <li>4. 2-year limited warranty on drivetrain.</li> </ol>				<b>Notes</b> <ol style="list-style-type: none"> <li>1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.</li> <li>2. Battery weight may vary based on manufacturer.</li> <li>3. Not available with iLevel Seat lift</li> </ol>			

### Get Social with Quantum



FOLLOW THE LIVES OF OUR QUANTUM STARS, AND GET THEIR PERSPECTIVE OF THE Q6 EDGE 2.0

[www.quantumrehab.co.uk/quantum-stars](http://www.quantumrehab.co.uk/quantum-stars)



@QUANTUMREHABUK

Ensure you have the most up to date information about Quantum Rehab by subscribing to our social media pages

#ileveluk

#myquantumproduct

#socialmobility

Further information can be obtained through our web site, or prescription forms. All information is subject to change without notification. Please consult Quantum with any queries you may have.