

Pride[®]

REVO 2.0



Tough and built to last, ensuring a smooth and comfortable ride over varied terrain

Introducing the Revo 2.0 a new experience in enhanced outdoor performance. With 35 Ah batteries and a 26 Stone 7lbs weight capacity the Revo 2.0 is easily transportable and can be disassembled for portability. The Revo 2.0 features Comfort-Trac Suspension (CTS) allowing for a smooth, more enjoyable journey.



**EASILY
DISASSEMBLED**



@Pridemobilityuk



*Revo 2.0 in
True Blue*

Version 1.0 | April 2018



32 Wedgwood Road, Bicester, Oxon, OX26 4UL | Tel: 01869 324600 | Fax: 01869 323070 |
Web: pride-mobility.co.uk | Email: sales@pride-mobility.co.uk



Revo 2.0 | Standard Features

Revo 2.0 in
Grey Street



Comfort-Trac
front and rear
suspension



Front and Rear
Lights



Solid Alloy Style
Wheels



Underseat
Storage



Large Tiller
Basket



Easily
Disassembled

and much more..

Revo 2.0 | Colour Options




Grey Street



True Blue

Technical Specifications



Weight Capacity	Max Speed	Ground Clearance	Turning Radius	Overall Length	Overall Width	Front/Rear Tyres	Total Weight w/o batteries	Heaviest piece	Standard Seating	Batteries Required
 26 Stone 7lbs (170.09 kg)	 4 Mph (6.4 km/h)	 3.17 cm ⁴ (1.25 inch's)	 132.08 cm (52 inch's)	 119.38 cm (47 inch's)	 53.34 cm (21 inch's)	 Front: 22.86 cm Rear: 22.86 cm	 9.35 Stone 59.4 k.g.	 25.85 k.g. rear section	 Folded Seat 12 k.g ²	 x2 35 Ah deep cycle ⁶

Warranty

1. Life time warranty on frame
2. Two-year limited warranty on drive train, including: transaxle, motor and brake.
3. One-year limited warranty on all electrical parts
4. 12 month warranty for batteries

Notes

1. Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
3. Battery weight may vary based on manufacturer.
4. Measured at transaxle
5. Seat measurements are as follows: Width: 50.8 cm (20 in) - 43.2 cm (17 in.) / Depth: 45.7 cm (18 in.); 43.2 cm. (17 in)
6. Battery Pack Weight: U1 battery pack: 11.5 kg (25 lbs)
7. Information correct at time of printing
8. Images used for illustration purposes only.