
Rascal Rio

Compact powerchairs



Your new route to independence



Rascal Rio

Rascal Rio Specification:

'	
Type Class	B/2
Overall length	88cm / 35"
Overall width	55cm / 23"
Floor to seat height	50.8 - 55.8cm
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	2.3cm / 0.9"
Maximum kerb height	2.5cm / 1"
Range*	up to 13km / 8 miles
The maximum safe slope	6°
Turn-around width	1.06m/3′5″
Colour	White
Single Heaviest component	14 4kg / 31 75lhs

Single Heaviest component 14.4kg / 31./5lbs

*Range on full charge and flat ground

Main Features:

- Super lightweight & transportable
- Flip-up fully adjustable armrests
- Solid tyres all round
- Convenient joystick charging port
- Padded seat with folded back. adjustable arm height, angle &
- Flip-up footplate



Optional Extras:

• Optional Rascal arm rest bag



Dismantles easily into 3 pieces to fit in a car boot for transportation

Rascal Rio Profile:

The Rascal Rio is a modern & stylish powerchair ideal for indoor and outdoor use. The super lightweight portable powerchair easily dismantles making it the perfect choice for transportation and storage in small spaces. Getting the seat position right for you is simple with fully adjustable armrests, seat height and joystick.

The high strength base makes the Rio a sturdy option

for anyone looking for a durable,

compact powerchair.

Local Rascal Mobility Retailer:



Specifications subject to change without prior notice Images are for illustrative purposes only E. and O.E- v 2b Final-Sept 21





