

## RETAIL

Breezy Ibis (Transit)

Invoice to:		
Name:		
Street:		
Town:	Post code:	
E-Mail:		
Tel:		
Order date:		
Purchase order:		

 Deliver to:

 Name:

 Street:

 Town:
 Post code:

 E-Mail:

 Tel:
 Fax:

 Marked for:

**Base Price** £2,950

Fax number: +44 (0)845 605 66 89 Email: orders@sunmed.co.uk

₽ CE

Standard features - IBIS Pro Transit						
Y01117	Sedeo Pro Seat Frame	Y00135	12.5" solid rear wheels	Y01231	Side Supports PU	
Y01126	Height adjustable push handles	Y01122	Folding backrest, tool operated	Y01255	Lap strap	
Y00113	7" solid front wheels	Y01215	Removable armrests	Y00005	Anti-tip bracket, flip-up	
Crash te	sted and Passed	Max User	Max User weight 160Kg			

UNIT-ID	0255		£	
	BAS	E FRAME		
	SED	EO PRO SEAT FRAME		
Y01117	•	Sedeo pro seat frame		
Y01126	•	Height adjustable push handles (not with third party backrest)		
Y01124		Third Party Backrest - width adjustable	£490	
Y01127		Angle adjustable push bar under 120kg (not included with third party/seat shell backrest frame/PowerSupport)	£310	
	SEA	T HEIGHT (without cushion)		
Y00586	0	Low seat post for SH 32 - 42.5 / 12.5 - 16.5"		
Y00587	•	High seat post for SH 36.5 - 47cm / 14.5 - 18.5"		
Y00575	0	32cm measured to top of seat board Y00581 O 41cm measured to top of seat board		
Y00576	0	33.5cm measured to top of seat board Y00582 O 42.5cm measured to top of seat board		
Y00577	0	35cm measured to top of seat board Y00583 ○ 44cm measured to top of seat board		
Y00578	0	36.5cm measured to top of seat board Y00584 O 45.5cm measured to top of seat board		
Y00579	0	38cm measured to top of seat board Y00183 O 47cm measured to top of seat board		
Y00580	0	39.5cm measured to top of seat board		
Y00113	•	7" Solid front castors		
Y00085	•	Axle position middle 3		
Y00135	•	12.5" Solid rear wheels		
	BRA	KES & POWER SUPPORT		
Y00092	•	Footbrake		
Y00093	0	Footbrake for SW 54cm		
Y02047 (CSA)	0	Attendant brake for use 12.5" (not available with footbrake)		
Y00117		Power Support		
	ADJ	USTMENTS (tilt & recline, backrest options on next page)		
Y01512	•	Attendant operated, Gas spring tilt adjustable - max user weight 125kg. tilt 30°		
Y01513		Attendant operated, Gas spring tilt adjustable - Recommended for user weight 90kg - 160kg. tilt 30°		
Y00400		Adjustable tilt in space operated by user (60 - 90 kg) max. tilt 18° - RH	£505	
Y00401		Adjustable tilt in space operated by user (90 - 125kg) max. tilt 18° - RH		
Y02039 (CSA)		Negative Tilt. tilt 5°	£165	
Y01119		Attendant recline. recline 89-128°	£310	

				£
	SEA	T SIZE	(adjustability)	
Y01139	0	42 x 50cm (16" x 20")	16" - 21" Seat Width x 17" - 20" Depth	
Y02644	0	42 x 56 cm (16" x 22")	16" - 21" Seat Width x 19" - 22" Depth	
Y01142	0	46 x 50cm (18" x 20")	18" - 21" Seat Width x 17" - 20" Depth	
Y02645	0	46 x 56 cm (18" x 22")	18" - 21" Seat Width x 19" - 22" Depth	
Y01145	0	50 x 50cm (20" x 20")	19" - 21" Seat Width x 17" - 20" Depth	
Y02646	0	50 x 56 cm (20" x 22")	19" - 21" Seat Width x 19" - 22" Depth	
Y01148	0	54 x 50 cm (21" x 20")	21" Seat Width x 17" - 20" Depth	
Y02647	0	54 x 56 cm (21" x 22")	21" Seat Width x 19" - 22" Depth	
	SEA	T CUSHIONS		
Y01159	•	Board Seat Pan seat with fr	ont (50 mm) & depth adj. rear lip	
Y01156		Stabilizer seat, 90mm thick		£165
Y01166		Abduction support (pommel	(not included with central legrest)	£245
	BAC	KRESTS		
		KREST SIZE		
Y01173	0	42cm Back width - Standard	I (set at 50cm, adjustable to 58cm)	
Y01174	0	42cm Back width - Tall (set	at 56cm, adjustable to 64cm)	
Y01176	0	46cm Back width - Standard	(set at 50cm, adjustable to 58cm)	
Y01177	0		at 56cm, adjustable to 64cm)	
Y01179	0		(set at 50cm, adjustable to 58cm)	
Y01180	0		at 56cm, adjustable to 64cm)	
Y01182	0		I (set at 50cm, adjustable to 58cm)	
Y01183	0		at 56cm, adjustable to 64cm)	
		. ,		
	BAC	KREST OPTIONS		
Y01122	•	Folding backrest, tool opera	ted (Standard option if nothing selected)	
Y01123		Folding backrest, thumbwhe	pel	£100
		<u> </u>	<del>,</del>	
	BAC	KREST CUSHIONS		
Y01190	0	Flexi Deep Contour, Passive	e, Adj. Straps (not compatible with thoracic supports)	
Y01191	0	Flexi Deep Contour, Active,		
Y01520	0	Flexi Extra Deep Contour, A	ctive, Adj. Straps (not compatible with thoracic supports)	
Y01185	0	Profiled Backrest	( compatible with thoracic supports)	
Y01193	•	Removable incontinence ba	ckrest cushion cover (included with all backrests)	
Y01194			cover including incontinence inner liner)	£130
Y01195		-	ed with third party or seat shell backrest frame) (Use with padded backrest)	£475
	LEG	RESTS		
Y01198	0	Swingaway legrests		
Y01202		Calf Strap		
Y01204		□ Depth adj. calf suppo	ort, PU, set	£210
Y01209	0	Centre Mount Legrest		
Y01210	0	Centre Mount Articulating Lo	egrest (from seat height 35cm)	
Y01212		Central calf support		
Y01213		□ Central calf support \	vith side support	£100
Y01199		Manually elevating legrests	(Must select Y01204)	£280
Y01204		Depth adj. calf support		£210
Y01205		Knee pads (not available)	e with padded side supports, wedge-shape)	
	ABAE	HITEE CHIDDODTO	T	
Y01208	AMF	PUTEE SUPPORTS  Amputee support, right	(not available with centre mount articulating legrest)	£195
Y01207		Amputee support, left	(not available with centre mount articulating legrest)	£195
		cappert, lost	,	

			£
		ARMRESTS & ARMPADS	
Y01216	•	Removable armrests	
Y01223		RIGHT HAND - Armrest pad 30cm	
Y01222		LEFT HAND - Armrest pad 30cm	
Y01221		RIGHT HAND - Armrest pad 40cm	
Y01220		LEFT HAND - Armrest pad 40cm	
Y01548		O RIGHT HAND - Armrest pad wide 45cm (not i.c.w. SW 54 cm, not i.c.w. 24" rear wheels)	
Y01547		O LEFT HAND - Armrest pad wide 45cm (not i.c.w. SW 54 cm, not i.c.w. 24" rear wheels)	
Y01504		Raised Armrests +10 cm (only i.c.w. removable armrests with quick-lock)	
Y01227		RIGHT HAND - Hemiplegia armrest and pad	£260
Y01226		LEFT HAND - Hemiplegia armrest and pad	£260
		SIDE GUARDS	
Y01231	•	Side Supports PU (SW -2cm) (Standard Option if nothing selected)	
Y01232		Padded side supports, wedge-shape (SW -8cm)	
Y01234		Thigh supports (available with Standard Side support but not available with padded wedge side support)	£310
		HEADRESTS (please pick headrest & mounting)	
Y01235	0	Headrest with PU Cover	
Y01236		Headrest flexible	£115
		Headrest Mounting	
Y01239	0	Headrest socket on backrest frame (Standard Option if nothing selected)	
Y00333	0	Headrest socket on backrest frame - only with third party backrest	
Y01240		Width adjustable headrest bracket (allows greater range of lateral adjustment)	£115
		ACCESSORIES	
Y00005	•	Anti-tip bracket, flip up	
Y01255	•	Lap strap	
Y00012		Cane/crutch holder for 2 canes	£70
		WORKTOPS	
Y01249		Slide on Tray (Fixed acrylic with two insert tubes, abdominal recess 11cm)	£310

Base Price	£2,950	
Option price	£	
Total price	£	

