

Invoice to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Deliver to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Order date:
 Purchase order:

Marked for:

Base Price
£2,629

Fax number: +44 (0)845 605 66 89
 Email: orders@sunmed.co.uk
 Max User Weight 65Kg

Transit approved to ISO 7176-19
 (only when transit bracket option is selected)


Part No 062001-000 PQR01



CLIENT MEASUREMENT


		Seat width in cm																	
		22	24	26	28	30	32	34	36										
Seat depth in cm	24	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	Each seat width can be extended by 2 cm without extra parts									
	26	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	28	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	30	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	32	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	34	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	36	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	38	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
Only available with abduction frame																			
Each seat depth can be extended by 4 cm without extra parts																			
		Front seat height in cm																	
		Fork (length)	Frame:	36	37	38	39	40	41	42	43	44	45	46	47	48			
Front castor	3"	111 mm	Low	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>			
		111 mm	High					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
	4"	111 mm	Low		<input type="radio"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>									
		123 mm	Low				<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>									
		123 mm (1)	Low				<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>									
		111 mm	High					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>						
		123 mm	High						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
		123 mm (1)	High							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>				
	5"	111 mm	Low					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>									
		123 mm	Low						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
		123 mm (1)	Low						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
		111 mm	High							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
123 mm		High								<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
123 mm (1)		High									<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
6"	111 mm	Low							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
	123 mm	Low							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
	123 mm (1)	Low							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
	111 mm	High										<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	123 mm	High										<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	123 mm (1)	High										<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
		Rear seat height in cm																	
Rear wheel	Wheels	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46
	20" (2)	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	22"						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	24"							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

High frame



269 mm

Low frame



229 mm

Your rear seat height can only be between 0 cm and 13 cm below your front seat height:

i.e. rear seat height **CANNOT** be above front seat height

i.e. rear seat height **CANNOT** be more than 13 cm below front seat height

. More info, see first

(1) fork compatible with additional drive

(2) only with Proton wheels/Right Run tyres

										£	
SEAT UPHOLSTERY											
Seat sling										(with pad)	
SIB020002	<input type="radio"/>	1 Utility bag								(Std default option)	
SIB020004	<input type="checkbox"/>	Utility bag and catheter bag								(only from SW 34 cm)	
SIB020000	<input type="checkbox"/>	Without seat sling								92	
Seat cushion with cover, black											
a) Foam composite		Thickness:	<input type="checkbox"/>	3cm	<input type="checkbox"/>	5cm				(medium)	
b) Latex		Thickness:	<input type="checkbox"/>	5cm							(very soft)
Seat accessories											
SIB090018	<input type="checkbox"/>	Positioning belt with lock								43	
Backrest (range)											
Back height set at (cm)		20 - 25	25 - 30	30 - 35	35 - 40						
		<input type="radio"/> 20	<input type="radio"/> 25	<input type="radio"/> 30	<input type="radio"/> 35						
		<input type="radio"/> 22.5	<input type="radio"/> 27.5	<input type="radio"/> 32.5	<input type="radio"/> 37.5						
		<input type="radio"/> 25	<input type="radio"/> 30	<input type="radio"/> 35	<input type="radio"/> 40						
BACKREST											
Backrest frame (4 cm seat depth growth built in)											
SIB030015	<input type="radio"/>	Foldable, angle adjustable					(74° to 110° adjustable)				(-16° to +20°)
SIB030102	<input type="checkbox"/>	Double locking back	(6)	(not in combination with 1240)			(74° to 110° adjustable)				(-16° to +20°)
SIB030530	<input type="checkbox"/>	Relax position (+10°)								(not in combination with 4016)	
Push handles											
SIB030200	<input type="radio"/>	Without push handles (rounded tube end)									
SIB030201	<input type="radio"/>	Std. push handles (long)									
SIB030203	<input type="checkbox"/>	Fold-down push handles	(1)							147	
SIB030204	<input type="checkbox"/>	Height-adjustable push handles	(6)							197	
SIB030201	<input type="checkbox"/>	Prepared for JAY Zip back					(not with Distance 20 mm for stabilising bar)			218	
SIB030208	<input type="checkbox"/>	Push bar, height-adjustable, width and angle-adjustable					(not with SW 22)			380	
SIB030405	<input type="checkbox"/>	Fitting kit (spacer) for shells/backs	(6)							79	
Backrest upholstery											
SIB030316	<input type="radio"/>	EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black (Evo Exo Version)									
Colour, binding tape:		<input type="radio"/> Black	<input type="radio"/> Orange	<input type="radio"/> Silver	<input type="radio"/> Red	<input type="radio"/> Blue					
SIB030317	<input type="checkbox"/>	EXO PRO Backrest, synthetic fibre, padded, adj., breathable spacer fabric, 1 accessory bag, black (Evo Exo Version)								162	
Colour, binding tape:		<input type="radio"/> Black	<input type="radio"/> Orange	<input type="radio"/> Silver	<input type="radio"/> Red	<input type="radio"/> Blue					
SIB030300	<input type="radio"/>	Without backrest upholstery									
Headrest											
SIB030411	<input type="checkbox"/>	Headrest small	(2)	SIB030410	<input type="checkbox"/>	Headrest large	(2)				(not in combination with 4022, 4024)
SIDE GUARDS AND ARMREST											
SIB040103	<input type="checkbox"/>	Straight side guard aluminium with fender (Pair) - for use with Safari brakes								278	
SIB040055	<input type="checkbox"/>	Tool height adjustable armrest - short pad	(3)	(includes aluminium side guard)						376	
Accessories											
SIB090111	<input type="checkbox"/>	Height adjustable tray table, slide-on for SW 22 - 36 cm								145	
Therapy table, flip-down for SW 26 - 36 cm		(4)								(only with height adjustable armrests with pad)	
Attachment:		SIB090113	<input type="checkbox"/>	right hand	SIB090112	<input type="checkbox"/>	left hand				212
LEGREST & FOOTPLATES											
Footplates											
SIB050025	<input type="checkbox"/>	Platform footrest angle and depth adjustable (composite footboard)								(LLL between 22 - 36 cm)	
SIB050056	<input type="checkbox"/>	Aluminum, angle and depth-adjustable, calf strap, black								(LLL between 22 - 36 cm)	
SIB050032	<input type="radio"/>	Platform footrest swing back (composite footboard)								(LLL between 22 - 36 cm)	
SIB050048	<input type="checkbox"/>	Platform footplates, high mount (composite footboard)								(for lower leg length between 17 - 25 cm)	
SIB050127	<input type="checkbox"/>	Foot positioning plate								(not in combination with swing back footrest)	

(1): Only available from back height 30 cm

(2): Only available from back height 30 - 35cm

(3): Not available with the Safari wheel lock - wheel lock integrated into Side guard (6001)

(4): Only available with tool height adjustable armrest - short pad (4227)

(5): For lower leg length between 22 - 36 cm

(6): If the combination of 4016 and 4034 are chosen then 4046 must be added

										£					
Lower leg length (LLL) (Attention: LLL has to be 2.5 cm lower than seat height front)															
<input type="radio"/>	17	<input type="radio"/>	18	<input type="radio"/>	19	<input type="radio"/>	20	<input type="radio"/>	22	<input type="radio"/>	23	<input type="radio"/>	24	<input type="radio"/>	25
<input type="radio"/>	26	<input type="radio"/>	27	<input type="radio"/>	28	<input type="radio"/>	29	<input type="radio"/>	30	<input type="radio"/>	31	<input type="radio"/>	32	<input type="radio"/>	33
<input type="radio"/>	34	<input type="radio"/>	35	<input type="radio"/>	36	All measurements are in cm									
FRAME 95° knee angle															
SIB010018	<input type="radio"/>	Straight standard frame (1)													
SIB010019	<input type="radio"/>	Abduction standard frame (3cm)													
AXLE PLUG AND CENTRE OF GRAVITY (COG) (all measurements are in cm)															
Choose only one of the 9 choices											(measured from front of back post in cm)				
<input type="radio"/>	2.00 cm														
<input type="radio"/>	3.00 cm														
<input type="radio"/>	4.00 cm														
<input type="radio"/>	5.00 cm														
<input type="radio"/>	6.00 cm														
<input type="radio"/>	7.00 cm														
<input type="radio"/>	8.00 cm														
<input type="radio"/>	9.00 cm														
<input type="radio"/>	10.00 cm														
Frame accessories (choose either tip assist or anti- tip tubes, not both)															
SIB090100	<input type="radio"/>	Crash test brackets (tie downs)										(includes SIB030405)			
Tip assist (step tube)															
SIB090001	<input type="checkbox"/>	Left	SIB090002	<input type="checkbox"/>	Right						35				
SIB090009	<input type="checkbox"/>	Crutch holder with Velcro										(only in combination with 1230) 76			
Anti-tip tubes, swing-away															
SIB090004	<input type="checkbox"/>	Left	SIB090005	<input type="checkbox"/>	Right						96				
SIB090006	<input type="checkbox"/>	Pair										195			
CASTOR FORK															
SIB080013	<input checked="" type="radio"/>	Fork Carbotecture®													
SIB090093	<input type="checkbox"/>	Fork for power add ons, 123x45 mm										74			
STANDARD CASTORS (compare with your choice on the first page)															
SIB080301	<input type="radio"/>	Solid													
SIB080302	<input type="radio"/>	Soft solid										(only for 5" & 6")			
Castors with aluminium rim															
SIB080304	<input type="checkbox"/>	Solid Aluminium										(only for 4" and 5") 89			
SIB080306	<input type="checkbox"/>	Soft roll castor with aluminium rim wide version										(only for 4" and 5") 89			
REAR WHEELS															
SIB070001	<input type="radio"/>	Cross spoked										(pair) (only available with 22" & 24")			
SIB070002	<input type="radio"/>	Design wheels (straight spoked)										(pair) (not with 22")			
SIB070003	<input type="radio"/>	Lightweight wheels (straight spoked)										(pair) 22+24"			
SIB070004	<input type="checkbox"/>	Proton wheels										(pair) (only available with 20", 22" & 24") 579			
SIB070009	<input type="checkbox"/>	Mountain bike wheels										(only with 24", mountain bike tyres & aluminium rim, black coated) 369			
SIB070010	<input type="checkbox"/>	Additional mountain bike wheels (pair)										(only with 24" & mountain bike tyres) 498			
TYRES (choose only one tyre type)															
SIB070101	<input type="radio"/>	Right Run										(available in 20"; 22" and 24")			
SIB070102	<input type="radio"/>	Solid tyres										(only 22", not available with Safari wheel locks, lightweight or proton wheels)			
SIB070107	<input type="radio"/>	Schwalbe One (Slight puncture resistance)										(only 24")			
SIB070104	<input type="checkbox"/>	24" Marathon plus evolution										(not 20" or 22" - puncture resistant) 91			

(1): available from seat width 28cm

										£
HANDRIMS										
SIB070202	<input type="radio"/>	Aluminium silver anodised								
SIB070201	<input type="radio"/>	Aluminium silver powder coated								
SIB070204	<input type="radio"/>	Stainless steel 22" with thin diameter (15 mm)								
SIB070203	<input type="checkbox"/>	Coloured aluminium pushrim. colour code:								66
SIB070205	<input type="checkbox"/>	Titanium								205
SIB070206	<input type="checkbox"/>	Titanium clear powder coated								205
SIB070207	<input type="checkbox"/>	Handrims, Supergrip ®								244
SIB070208	<input type="checkbox"/>	Handrims, Maxgrepp ®								244
		Hand rim mounting	<input type="radio"/>	Standard = 17mm	<input type="radio"/>	Wide = 27mm	<input type="radio"/>	Close = 7mm		
SIB070250	<input type="checkbox"/>	Handrim cover								99
WHEEL CAMBER										
(choose only one option)										
		<input type="radio"/>	0°	<input type="radio"/>	3°	<input type="radio"/>	6°	<input type="radio"/>	9°	
Spoke guards Please refer to the brochure and provide number/description										
		<input type="checkbox"/>	Spoke guards	Number		Description		134		
Please consult the following page to choose your selection or to order a copy of our brochure please call 0845 605 66 88 or e-mail literaturerequests@sunmed.co.uk You can also view the brochure online by visiting www.sunrisemedical.com and navigating to Quicke Simba. You will find the brochure under Multimedia & Downloads										
WHEEL LOCKS										
(choose only one option)										
SIB060001	<input type="radio"/>	Standard wheel lock								
SIB060002	<input type="checkbox"/>	Knee lever wheel lock								63
SIB060003	<input type="checkbox"/>	Compact wheel lock								48
SIB060013	<input type="radio"/>	Safari wheel lock (wheel lock integrated into side guard - must pick side guard)								
SIB060014	<input type="checkbox"/>	Attendant wheel lock with Safari wheel lock (wheel lock integrated into side guard - must pick side guard)								152
Wheel lock accessories										
SIB060010	<input type="checkbox"/>	Extension for knee lever wheel lock, slip-on, swing-away								46
ACCESSORIES										
SIB090000	<input checked="" type="radio"/>	Tool kit								
SIB090024	<input type="checkbox"/>	High pressure pump								89
SIB090025	<input type="checkbox"/>	Air pump								17
SIB090030	<input type="checkbox"/>	Back pack								65
SIB090016	<input type="checkbox"/>	Impact guard								68
SIB090026	<input type="checkbox"/>	Mobile phone holder								48
COLOURS										
		Frame	Handrim (1)	Hub/rim (2)			Frame	Handrim (1)	Hub/rim (2)	
			£66	£188				£66	£188	
02	Brilliant Silver	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	21	Red Ruby	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
04	Silver	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	22	Brilliant Red	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
06	Sun Yellow	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	23	Matt Red	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
14	Blue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	27	Graphite metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
15	Brilliant Marine	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	29	Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
16	Matt Marine	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	31	Matt Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
19	Brilliant Purple	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	32	Signal White	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
40	Kiwi	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	33	Pink	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
41	Flamingo Glitter	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	34	Dark Chocolate Metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
42	Matt Anthracite	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	35	Skyblue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
46	Mamba Green	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	36	Cyan-Metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
					37	Mars Red (Orange)	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Upholstery colour										
(vented only in black)										
		Black	<input type="radio"/>	Purple	<input type="radio"/>	Dark blue	<input type="radio"/>			
		Red	<input type="radio"/>	Blue	<input type="radio"/>	Light Pink	<input type="radio"/>			
Style Package (coloured forks)										
(vented only in black)										
SIB100502	<input type="radio"/>	Black								
SIB100503	<input type="radio"/>	Silver								
SIB100501	<input type="radio"/>	Orange								
(1) Only available on Universal and Design wheels										
(2) Only available on Universal wheels										
Simba - Base Price							£2,629			
Option price							£			
Total Price							£			

Please follow measuring method:

The following method is recommended on the basis of extensive practical experience.

Deviations are possible in individual instances

Seat width:

a) For the user:

Allow space between side guard and between hip on each side of the chair (approx. 1 - 2 cm altogether).

b) On the wheelchair:

From inside side guard to inside side guard

Seat depth:

a) For the user:

When sitting upright: leave space for 2-3 fingers between front edge of seat upholstery and back of knee (approx. 3-6 cm). Care must be taken with adjustable upholstery.

b) On the wheelchair:

From back tube to front of seat upholstery.

Back height:

a) For the user:

Measured from the seat, including cushion, to be **no more** than one finger (1 - 2 cm) below the lower shoulder blade.

Important! The back height depends on the degree of disability.

b) On the wheelchair:

From seat tube to upper edge of backrest upholstery.

Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably.

Lower leg length:





a) For the user:

With legs angled at right angles, from back of knee to bottom of shoe sole.

Remember: to deduct the cushion height when providing the measurement

b) On the wheelchair:

From front of seat to top of footrest.



Overall width with:

	20"	22"	24"
For 0° camber:	SW+16cm	SW+16cm	SW+16cm
For 3° camber:	SW+20cm	SW+20cm	SW+21cm
For 6° camber:	SW+26cm	SW+26cm	SW+27cm
For 9° camber:	SW+31cm	SW+32cm	SW+33cm

each with handrim distance std

Please note that the matrix is based on a std configuration

Std chair 30 x 32 cm, 6° camber, 6 cm CoG

Seat angle of 6°

All changes in the configuration might change the seat heights.

Manufacturing tolerances of +/- 1 cm

NOTES:

Max seat height difference rear to front is 13 cm.

No negative seat angle is possible

Please note that the matrix is based on a standard configuration

Overall width:

Max overall height: 88 cm
Max overall length: 85 cm

Chair specification:

Long frame; seat depth: 42cm; COG: 2cm;
24" wheel and 6" front castor


Min overall height: 61 cm
Min overall length: 57 cm

Short frame; seat depth: 24cm; COG: 10cm;
20" wheel and 4" front castor



Weights:

Minimum weight:
approx from 8 kg

Maximum measurements





Minimum measurements




Frame abduction table


Seat width	22	24	26	28	30	32	34	36	38
Frame type	cm	cm	cm	cm	cm	cm	cm	cm	cm
Straight				23	25	27	29	31	33
Abduction	23	25	27	29	31	33	35	37	39



Seat width growth of 2 cm without new parts




Seat depth growth of 4 cm without new parts




By loosening those 2 screws

Backrest height: 20 – 25 cm range



Short. Backrest frame & short push handle tubes


Backrest height: 25 – 30cm range



Medium Backrest frame & short push handle tubes


2 new seat depths without using new parts. In Total additional 4cm

Backrest height: 30 – 35 cm range



Medium Backrest frame & medium push handle tubes

Backrest height: 35 – 40cm range



Long Backrest frame & long push handle tubes

SPOKE GUARDS - 22" & 24" ONLY



5000 - Princess Life



5007 - Monkeys



5002 - Jungle Fever



5004 - Happy Lions



5025 - Oriental



5008 - Pirates



5015 - Feather



5003 - Superhero



5018 - Graffiti



5017 - Chrome Wheel



5016 - Hot Wheels




5023 - Basketball



5024 - Basketball 2



5013 - Pool



5022 - Flying Water



5021 - Space



5020 - Silk



5010 - Zebra



5134 - Galactica



5009 - Tiger