

Invoice to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	

Deliver to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Marked for:	

E-mail: orders@sunmed.co.uk

Fax: +44 (0)845 605 66 89

Max user weight limit 140kg (22 stone)

<ul style="list-style-type: none"> ● Powered tilt (3° to 25°) ● Articulating drive wheel suspension ● 12" pneumatic or solid tyres ● 8" solid or pneumatic castors ● PG VR2 controller 	Standard Features <ul style="list-style-type: none"> ● Parallel swing-away control arm ● Width adjustable armrests ● Depth adjustable back (16.5" to 19.5") ● Fold down, angle adj. back (-3° to 12°) 	<ul style="list-style-type: none"> ● 70° swing-in/swing-out legrests ● 40Ah twin batteries ● 6kph (4mph) high torque motors ● Lap belt
---	--	--

Base Price
£4,646.00

Tango Seat Height Reference Matrix	
Seat Angle	Front Seat Height
3°	43cm (17")
6°	45cm (17.75")
9°	48cm (18.75")
12°	51cm (20")

Seat height matrix shows what front seat height can be achieved using power tilt (0° - 30°)

(!) All measurements are with 12" drive wheels - add 1cm to height when using 14" drive wheels

					£		
FRAME WIDTH:							
TANW043	● 43cm (17") frame width	(Adjustable outwards 43 - 50cm via armrests / inwards 43 - 38cm via side guard reduction pads)					
SEAT DEPTH							
Seat depth setting (42.5cm to 50cm seat depth adjustment via backrest)							
TAND0425	○ 42.5cm (16.5")	TAND0450	○ 45cm (17.5")	TAND0475	○ 47.5cm (18.5")	TAND0500	○ 50cm (19.5")
CUSHION							
TAN020011	○ Omit cushion						
TAN020001	○ 5cm (2") standard foam cushion 43cm wide						
	□ JAY Comfort seat (cushion & back)	(Seat cushion always supplied same width as frame: 43cm)			370		
TAN020002	○ 40cm width	TAN020003	○ 44cm width	TAN020004	○ 48cm width		
SEAT TILT							
TAN020007	○ Powered tilt (3° to 25°)						
TAN020005	○ Fixed tilt (0° - 9° adjustable, factory setting 3°)	(Factory setting 3°)					
TAN020008	□ Seat Height Kit (increases STF height by 30mm)					53	
BACKREST							
					(Fold down Backrest)		
TAN030000	○ Std. sling backrest (48cm (19") height)	(adjustable back height 18" to 21")					
TAN030001	○ Mechanical fixed recline back (6° factory setting)	(-3° to 12; in 3° steps)					
TAN030002	□ Manual adjustable recline back (-3° to 30°)	(locking gas cylinder)			296		
TAN030003	□ Powered recline back (-3° to 30°)	(max two powered options possible e.g. tilt & recline or tilt & powered legrests)			447		
TAN090024	□ Jay Shape Laterals (15cm x 12cm)				229		
HEADREST							
TAN030007	□ Comfort headrest	(Mounted via backrest when Comfort back selected, Mounted via push handles with std upholstery)			257		
ARMREST							
	● Width, depth, height adjustable armrest	(Rumba Modular type, moulded arm pad gives +7cm width adjustment)					
TAN040001	□ Side guard reduction pads	(Reduces seat width by 5cm/ 2")			68		

			£
LEGRESTS			
TAN050000	<input type="radio"/> 70° swing-away legrest		
TAN050001	<input type="checkbox"/> Manual elevating legrest		220
TAN050003	<input type="radio"/> 80° Uni style hangers and footplates	LLL 32-47cm (12.5"-18.5")	
	<input type="checkbox"/> Amputee support - RH	(UK market personalisation, process order via B-4-M: pt # 064280-001)	80
	<input type="checkbox"/> Amputee support - LH	(UK market personalisation, process order via B-4-M: pt # 064280-002)	80
Extension tubes for swing-away leg rests:			
	<input checked="" type="radio"/> Medium	34 - 48 cm (13" - 19")	(measured from seat pan to footplate)
	<input checked="" type="radio"/> Composite, angle and depth adjustable footplates		
VR2 60Ah CONTROLS			
			(including parallel swing-away controller mount)
TAN011011	<input type="checkbox"/> VR2 90Ah Controller Option	(Recommended with 10kph motor option)	174
TAN011008	<input type="radio"/> VR2 Joystick (inc actuator buttons)		
TAN011001	<input type="radio"/> VR2 Joystick for lights & Indicators		
	TAN011003 <input type="radio"/> RH	TAN011004 <input type="radio"/> LH	
TAN011002	<input type="checkbox"/> VR2 remote dual control (mounted RH)	(dual control module always mount right hand side)	538
TAN011009 / 10	<input type="checkbox"/> VR2 attendant only with actuator control	(not available with TA backrest)	(Lever drive / tiller style RH mount)
TAN011012 / 13	<input type="checkbox"/> VR2 Attendant only with actuator control		(Thumb drive / standard RH mount)
KERB CLIMBER			
TAN090002	<input type="checkbox"/> Kerb climber		196
DRIVE WHEELS & CASTORS			
TAN070000	<input type="radio"/> 12" Pneumatic drive wheels	8" Solid castors	
TAN070001	<input type="radio"/> 12" Solid drive wheels	8" Solid castors	
TAN070002	<input type="checkbox"/> 14" Pneumatic drive wheels V tread	9" Solid castors	144
TAN070003	<input type="checkbox"/> 14" Solid drive wheels V tread	9" Solid castors	144
TAN070007	<input type="checkbox"/> 14" Pneumatic drive wheels Knobby tread	9" Solid castors	144
TAN070008	<input type="checkbox"/> 14" Solid drive wheels Knobby tread	9" Solid castors	144
TAN070004	<input checked="" type="radio"/> Drive wheel mud guards		
MOTOR SPEED			
TAN010011	<input checked="" type="radio"/> 6kph (4mph)		
TAN010010	<input type="checkbox"/> 10kph (6mph)		322
FRAME COLOUR			
TAN090004	<input type="radio"/> Brilliant Silver		
TAN090005	<input type="radio"/> Brilliant Blue		
TAN090006	<input type="radio"/> Brilliant Red		
LIGHTS & INDICATORS			
TAN090000	<input type="checkbox"/> Lights & indicators		380
BATTERIES			
TAN016000	<input type="radio"/> 40Ah gel batteries		
TAN016001	<input type="checkbox"/> 50Ah gel batteries		220
CHARGER			
TAN016003	<input checked="" type="radio"/> 8A charger		
ACCESSORIES			
TAN090007	<input checked="" type="radio"/> Transportation tie down points		
TAN090008	<input type="checkbox"/> Crutch holder		73
TAN090009	<input type="checkbox"/> Swing-away tray table		206
TAN090010	<input checked="" type="radio"/> Adjustable lap strap		
TAN090013	<input type="checkbox"/> Rear view mirror		46

Quickie Tango - Base Price	£4,646
Option Price	£
Total Price	£