



▼ TIFIT. PERSONAL ERGONOMICS. ►

The difference between a TiFit® chair and all other chairs is that a TiFit chair is a wheeled prosthetic: Every frame is made to measure. A TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do. TiFit...not standard.

CUSTOM CONFIGURATION. ⊢

Custom Configuration allows you to add options and components to your chair to meet your needs and to optimize your riding experience. This includes choices of wheels, handrims, tires, wheel locks, color, tattoos, and more. TiLite offers more options than any other manufacturer!

TECHNICAL SPECIFICATIONS ⊢

| Frame Features | Titanium Mono-Tube Rigid |
|--------------------|----------------------------|
| | TiFit |
| | Limited Adjustability |
| Total Chair Width* | 23" |
| Seat Width | 10" – 20" |
| Seat Depth | 10" – 20" |
| Front Seat Height | 16" – 21" |
| Rear Seat Height | 12.5" – 20" |
| Front Frame Angle | 70°, 75°, 80°, 85°, or 90° |
| Seat Back Type | Folding Seat Back |
| | Adj. Height 8.5" – 21" |
| | Adj. Angle 80° – 101° |

| 0°, 2°, 4°, 6°, 8° or 12° |
|---|
| 6" |
| 0", 1.5", 2.5" |
| Transit Tie-Down |
| 10.7 lb. |
| 265 lb. Standard 350 lb. Heavy Duty |
| Five Years on Frame One Year on TiLite Components |
| |

^{*} Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.2 lb.

COLOR OPTIONS

HIGH GLOSS PAINT



TiLite Red Burgundy







Pink

Metallic









Deep

Blue Ice









Blue



Plum

NATURAL

Satin Polished Titanium Titanium

METALLIC PAINT



Metallic



Metallic







Metallic





ANODIZED COLORS









Black

PEARLESCENT PAINT





White Midnight Blue Pearl Pearl



Pepper

Orange

Sapphire Ultraviolet Fuchsia