



Deliver to:		Invoice to:	
Name:	(Name:	
Street:	`	Street:	
Town:	Postal code:	Town: Post	tal code:
E-Mail:		E-Mail:	
Tel:	Fax:	Tel: Fax:	
Order date:		Marked for:	
Purchase order:			
Number required:		Sunrise order no:	Continu
□ Order	□ Calculation □ Built 4 me	■ = Standard○ = No cost option	= Cost option
Part No: 475	500000 / 47500001 IC	Max. user weight: USER WEIGHT Mandatory information to ensure correct configuration	100 kg kg
	ATTITUDE - MOD	DEL SELECTION	
ATT010801 □	ATTITUDE Manual 30 speed gearing; Dual Drive with 3 speed integral gea casette and derailleur; trigger shifter; rear and front lam system; disk brake; parking brake; handle; v-crank.		£3095
ATT072102 ●	1 ATTITUDE Hybrid 10 speed casette and derailleur; trigger shifter; 17Ah Li motor support up to max 25kph; front and rear lamp; do wheels, handle; v-crank, disk brake; parking brake; spe Connect Display - coloured, includes smartphone conn	ocking system; stand with eedometer integrated in control display. ectivity via Bluetooth	£4815
FMG-MJ02-3187-02	B4Me option: Tetraplegia operation for Neodrives p	oower unit (1)	
ATT010803 □ ATT072111 ● ATT072112 ○	ATTITUDE Power Speed 6 kph (standard); handlebar; disk brake; 9Ah Lit front and rear lamp; docking system; speedometer integrational display; stand with wheels; parking brake. Hand throttle position right side Hand throttle position left side		£4540
	WHEELCHAIR MODEL - which will be us	sed in combination with the ATTITUDE	
	Suprice Wheelshaire	Other Wheelchairs	
	Sunrise Wheelchairs O Life T (fixed front frame only)		216
	Life T (fixed front frame only)Life RT (fixed front frame only)	☐ Küschall K-serie; KSL (2) ☐ Küschall Advance, Champion (2)	216
	O Life (fixed front frame only)	, , ,	216
	O Life R (fixed front frame only)	☐ Otto Bock Avantgarde CS/CLT (2/3) ☐ Otto Bock Motus CS (2/3)	216
	,	` '	216
	 Neon (fixed front frame only) Xenon/Xenon² (fixed front frame only) 	☐ Meyra Avanti Pro (2) ☐ Meyra Hurricane (2)	216
	O Argon² (fixed front frame only)	☐ Meyra Hurricane (2) ☐ Pro Activ Traveler (2)	216
	ÿ		216
	O RGK Tiga (2)	□ Pro Activ Speedy (2) □ Progeo Joker (2)	216
	O RGK Tiga (2)	□ Progeo Joker (2) □ Panthera S2, S3, U2, U3, (1) (4)	216
		□ Panthera 32, 33, 02, 03, (1) (4) □ Panthera U2 Light, U3 Light (1) (4)	216
	ORDER PROCESS	SINFORMATION	
ATT090213	O ATTITUDE only - without a related wheelchair order	er	
ATT090214	O ATTITUDE with wheelchair order - on a separate v		
FMG-MA03-3201-01	- B4Me option: Adaptation link -2 to + 6 cm (1)		178
(2) Please note: V-Fra(3) Only possible for fra	nt lead-times and discount structures imes are only suitable up to 5°; please check with your authorised deale ames without inset I hook bracket for short chairs to position away from the body	er if there are any further NoGo´s	

ATTITUDE = Standard O = No cost option = Cost option WHEELCHAIR FRONT FRAME WIDTH - measured from outside front frame tube to outside front frame tube 30 cm 26 cm 28 cm 32 cm 34 cm 36 cm 38 cm \circ \circ \circ \circ 0 0 \circ 40 cm 42 cm 44 cm 46 cm C 0 0 C Frame width MEASURED FROM OUTSIDE FRONT FRAME TUBE TO OUTSIDE FRONT FRAME TUBE (outside tubes) **Janual Hybrid** Power **FORK** ATT080035 Fork, angle and depth-adjustability of crank-position ATT080034 Fork, Standard (not angle and depth-adjustable on crank-position) **TYRES** ATT070104 Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution ATT070110 Tyre, pneumatic, tread, Mountainbike ATT070111 Tyre, pneumatic, tread, Big Apple **GEAR AND SHIFTERS** Dual drive (3 x 10 speed) for Manual Gear, dual drive, 3 speed integral hub, 10 speed casette / derailleur ATT070046 ATT070045 Shifter, trigger ATT070042 ☐ Shifter, grip shift, 10 speed 130 (1) ATT070043 Gear hub, 8 speed, coaster brake 10 speed casette with derailleur for Hybrid ATT070045 Shifter, trigger Shifter, grip shift, 10 speed 130 ATT070042 Shifter mount right ATT070048 O Shifter mount left ATT070047 ATT070044 Shifter, universal, 1 lever 270 (base to add levers) FMG-MJ02-3187-01 B4Me option: Tetraplegia operation for universal shifter (2) 119 CRANK AND HANDLE - special length, widths and shapes via B4Me ATT080022 V-crank, standard Length: 17 cm Width: 47 cm ATT080023 Handle, vertical, 10° angle, round, diameter 30 mm 130 ATT080027 ☐ Quad, hand support (Only available with ATT080023) ☐ Quad, arm support (Only available with ATT080023) ATT080028 130 Handle, vertical, 10° angle, round, diameter 40 mm ATT080024 ATT080025 Handle, vertical, 25° angle, oval, diameter 37 x 21 mm 130 ATT080026 Handle, vertical, 25° angle, oval, diameter 47 x 32 mm 130 HANDLEBAR & STEM - ATTITUDE Power only: ergonomic handlebar & height-adjustable stem standard ATT080050 Stem standard; height-adjustable ATT080051 Stem multi- adjustable; height- and depth adjustable DISC BRAKE - disc-diameter 180 mm, parking brake: V-brake (Power: 2nd caliper on the disc) ATT060050 Brake (manual) (not with gear hub, 8 speed, coaster brake) ATT060052 2 Driving brakes plus parking brake (3) ATT060053 Performance brake lever (offers improved adjustability of the brake) 54 ATT060102 Brake lever mounted to right handle ATT060101 Brake lever mounted to left handle 324 ATT060051 Brake, backpedal actuation (4)(not with gear hub, 8 speed, coaster brake) (1) only for 10 speed derailleur, trigger for 3 speed hub. (2) Built-4-me - different lead-times and discount structures; only in combination with ATT070044 (3) parking brake as rim brake and not as disk brake (4) not with ATT080035 Fork, angle and depth-adjustability of crank-position

	_						=				
ATTITUDE	Ξ						ına	rid	/er		
	= Standard	Ω	= No cost o	notion	□ = Cost option		Manual	Hybrid	Power	£	1
				<u> </u>	· .		2	<u> </u>	п.		
ATT000005		- Selection f	or Allilu		only			<u> </u>			1
ATT000025 •	Speed, ma			(1)			-	-	✓		1
ATT000027 O				(1)			-	-	1		
ATT000028 O	Speed, ma			(1)			-	-	1		
ATT000026 O	0 0 0 0 0, 1110			(1)			-	-	✓		
		R Y - chargeı		cluded							
ATT160020 ●	Battery, Li	thium-ion, 1	7Ah				-	✓	-		
ATT160021 ●	Battery, Li	thium-ion, 9	Ah				•	-	/		
ATT160022 □	Battery, Li	thium-ion, 1	3.4Ah				-	-	\	270	
ATT090088 □	Additional	battery, Lith	ium-ion, 1	7Ah			-	1	-	756	1
ATT090086 □	Additional	battery, Lith	ium-ion, 9	Ah			-	-	1	648	1
ATT090087 □		battery, Lith					-	-	1	702	1
ATT090091 🗆		charger, Hy					-	1	-	108	1
ATT090089 🗆		charger, Po					_	Ť.		108	1
7177000000	ACCESS		WCI						_	100	L
ATT090074 ●							1		1		l
	Front & rea				(introverted in the	diameter family desired and Danier		· /			-
ATT090080 •	Speedome				(integrated in the	display for Hybrid and Power)	1	✓	✓		
ATT090095 •	Mudguard						1	1	1		-
ATT090067 🗆							1	✓	✓	11	4
ATT090081 □	Bottle hold						1	\	\	16	1
ATT090082 □		carrier for sa					1	1	1	54]
ATT090024 □		sure pump (1	1	\	32]
ATT090084 O	Additional	weight, rem	ovable, 5k	g, to improv	e traction		1	1	\	0	
ATT090085 □	Bell						1	1	1	5	1
ATT090041 □	Multi tool						1	1	1	22	1
ATT090096 □	Protection	foil					1	1	1	43	1
			OUR AT	TITUDE (talanna fan ananka and akain e	guards are only available options for At	diducal a ma		ما الدياديا		
COLOUR CHOIC		FRAME CO			JARD COLOUR	CRANK COLOUR	uluue III	anuai an	u nybnu	,	
01 White RAL 90		ATT100101	0	ATT100801	O	ATT100901 O					
02 Brilliant Silver	10	ATT100101	0	ATT100801	0	ATT100901 O					
03 Silver Matt			0		0						
		ATT100103		ATT100803		ATT100903 O					-
04 Silver		ATT100104	0	ATT100804	0	ATT100904 O					
06 Sun Yellow		ATT100105	0	ATT100805	0	ATT100905 O					
07 Brilliant Lime		ATT100106	0	ATT100806	0	ATT100906 O					
08 Matt Lime		ATT100107	0	ATT100807	0	ATT100907 O					
14 Blue		ATT100110	0	ATT100810	0	ATT100910 O					
15 Brilliant Marine	e	ATT100111	0	ATT100811	0	ATT100911 O					
16 Matt Marine		ATT100112	0	ATT100812	0	ATT100912 O					
19 Brilliant Purple)	ATT100113	0	ATT100813	0	ATT100913 O					1
20 Matt Purple		ATT100114	0	ATT100814	0	ATT100914 O					1
21 Red Ruby		ATT100115	0	ATT100815	0	ATT100915 O					1
22 Brilliant Red		ATT100116	0	ATT100816	0	ATT100916 O					1
23 Matt Red		ATT100110	0	ATT100817	0	ATT100917 O					†
24 Tutti Frutti Ora	nge	ATT100117	0	ATT100817	0	ATT100917 O					1
			0	1	0	ATT100918 O					1
27 Graphite meta		ATT100120		ATT100820							1
28 Platinum meta	IIIC	ATT100121	0	ATT100821	0	ATT100921 O					-
29 Black		ATT100122	0	ATT100822	0	ATT100922 O					4
31 Matt Black		ATT100123	0	ATT100823	0	ATT100923 O					
32 Signal White		ATT100124	0	ATT100824	0	ATT100924 O					1
33 Pink		ATT100125	0	ATT100825	0	ATT100925 🔘]
34 Dark Chocolat	e Metallic	ATT100126	0	ATT100826	0	ATT100926 ○					
35 Skyblue		ATT100127	0	ATT100827	0	ATT100927 🔘]
36 Cyan-Metallic		ATT100128	0	ATT100828	0	ATT100928 O					1
37 Marsred (Oran	nge)	ATT100129	0	ATT100829	0	ATT100929 O					1
39 Cappuccino	<u> </u>	ATT100131	0	ATT100831	0	ATT100931 O					1
40 Kiwi		ATT100131	0	ATT100831	0	ATT100931 O					†
41 Flamingo-Glitte	or	ATT100132 ATT100133	0	ATT100832	0	ATT100932 O					1
				1							1
42 Matt Anthracite	Ħ	ATT100134	0	ATT100834	0	ATT100934 O					1
43 Azure Blue		ATT400405								•	i
		ATT100135	0	ATT100835		ATT100935 O					1
You can choose sep		for the frame,	the cranks	and the chain	guard. (co	olours for cranks and chain guard	only fo	or ATTI	TUDE M	lanual & Hybrid))
You can choose sep	which allow a h	for the frame, higher speed or	the cranks a	and the chain private ground		olours for cranks and chain guard	only fo	or ATTI	TUDE M	lanual & Hybrid))


 $if the \ confirmation \ on \ the \ back \ page \ is \ signed \ by \ the \ endorser!$

ATTITUI	DE	
		= Standard O = No cost option □ = Cost option
	-	
		STICKER COLOURS - for "ATTITUDE" and "Manual / Hybrid / Power" on the fork
ATT100550	•	Black
	0	Silver
	Ō	Orange
		COMMENTS FOR SPECIAL BUILD - please apply for a quote; different lead time and discount structure.
		built4me
Quickie A	TT	E £ £ £
Option price	!	
Total Price		Į£

ATTITUDE


ATTITUDE MANUAL & HYBRID ADJUSTMENT AND SET-UP

ATTITUDE Manual / Hybrid with Standard fork: crank position range (adjustment via spacers in the upper docking clamp)




ATTITUDE MANUAL & HYBRID fork: Angle and depth adjustability

ATTITUDE Manual / Hybrid - crank position range with option ATT080035 "Fork, angle and depth-adjustability of crank-position "



Depth adjustability: 100mm (+/- 50mm compared to standard fork) continuously adjustable

Angle adjustment: from 0° (parallel to the ground) to 28° continuously adjustable



B4ME-OPTIONS

FMG-MJ02-3187-02 - Tetra operation for ATTITUDE Hybrid



FMG-MA03-3201-01 Adaptation link -2 to + 6 cm



To move the pedal position 2 cm closer or up to 6 cm further away

ATTITUDE



Fork for Attitude Manual & Hybrid, angle and depth-adjustability of crank-position

TYRES



Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution



Tyre, pneumatic, tread, Mountainbike



Tyre, pneumatic, tread, Big Apple



Stem, multi adjustable (ATTITUDE Power only)

SHIFTER AND HANDLE



Shifter, trigger



Shifter, grip shift, 10 speed



Shifter, universal, 1 lever



FMG-MJ02-3187-01 Tetra operation



Handle, vertical, 10° angle, round, diameter 30 mm



Quad, hand support



Quad, arm support



Handle, vertical, 25° angle, oval,

DISC BRAKE



Performance brake lever

ATTITUDE

ACCESSORIES



Brake, backpedal actuation



Mudguard, bolt on, black



Multi tool



Rear view mirror



Additional weight, removable, 5 kg, to increase traction



Bell



Connect Display
Connects with your smartphone (download of neodrives app necessary) Records trip data such as average speed,maximum speed, trip distance, and total distance.



High pressure pump (0 - 11 bar)



Low rider carrier



Bottle holder





LETTER OF ACKNOWLEDGEMENT

We would like to point out that using the 10/15/21 Kph power version of the attitude in road traffic is only permitted in certain countries, otherwise it may only be used on private land.

In the absence of approval, in accordance with the road traffic laws, the *ATTITUDE Power* version 10/15 and 21kph version may not be driven on public roads or even cycle ways. You need to check your country specific legislation about the use and insurance duty of the attitude 10/15/21kph power version.

When placing an order, please sign below to confirm that you have understood the information and that we should deliver the 10/15/21kph version of Attitude Power.

If you want to use the *ATTITUDE Power*10/15/21kph versions in another country, you must first obtain information about the regulations and legislation in force there.

I hereby confirm, that I am aware of these circumstances and that I will strictly the country specific legislation.
Customer signature and date:

