

Drive
DeVilbiss
Powered
Mobility



drive

DeVilbiss
HEALTHCARE

www.drivedevilbiss.co.uk



Our Powered Mobility Range:

SCOOTER RANGE

04 - 18

ST1	04
ST5D	05
ST6	06
SCOUT RANGE	07
PRISM	08
PRISM SPORT	09
ENVOY RANGE	10 - 11
VIPER	12
COBRA	14
KING COBRA	15
ROYALE 3 & 4	16
ROYALE 4 SPORT	17
SPORT RIDER	18
EASY RIDER	19

POWERCHAIR RANGE

21 - 32

ENERGI RANGE	21 - 22
COBALT	23
CIRRUS	24
TITAN	25
SUNFIRE PLUS GT	26
IMAGE EC	27
SUNFIRE GENERAL	28
SEREN RANGE	29 - 30
MULTEGO	31
VOLT	32

POWERSTROLL RANGE

34 - 36

S-DRIVE	34
LIGHTWEIGHT POWERSTROLL	35
U-DRIVE	36



The ST1 Scooter is the lightest in the ST range and has an ultra compact design which is ideal for nipping to the shops, or for putting in the car boot to take on day trips and holidays. The streamlined ST1's compact and stylish frame uniquely splits into two or five pieces for effortless lifting. Despite being the lightest scooter in the range, the ST1 has plenty of power to easily take you up slopes and hills.

- › Top speed of 4mph*
- › Maximum range of 6 miles on a full battery charge*
- › Height adjustable, hard wearing and fully cushioned profiled seat
- › Plug and wire free connections for quick and simple assembling and dismantling
- › Splits easily for storage and transportation
- › Width adjustable armrests and adjustable handlebars to suit comfort of user
- › Comfortable lightweight and rotating seat which swivels 360° for easy transfers on and off the scooter
- › Delta bars as standard
- › Removable seat, easy tiller adjustment and an integral lifting handle for easy transportation
- › Lightweight scooter weighing only 40kg (88lb)
- › Horn to alert others of your presence for added safety
- › Fitted with reflectors and illuminated controls to improve outdoor visibility at dusk or when dark
- › Front removable storage basket as standard
- › Grey puncture proof non marking tyres
- › Rear anti tip wheels as standard
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Easy to remove battery pack fitted with 2 x 12 Ah airline friendly batteries
- › Includes an on or off board charging facility as standard
- › Available in Ruby Red or Sapphire Blue
- › Optional extras - Rear basket and additional storage options
- › Weight capacity 115kg (18st)
- * Subject to conditions of use

Colour Options

- Ruby Red
- Sapphire Blue

“Tight turning circle for ease of use indoors”

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	WEIGHT	HEAVIEST COMPONENT	MAXIMUM WEIGHT CAPACITY
ST1BE	ST1 Scooter	100cm (39")	89cm (35")	51cm (20")	6°	40kg (88lb)	14kg (31lb)	115kg (18st)



The ST5D Scooter has a modern, robust design that gives you the latest in scooter styling and performance. It features bespoke alloy wheels, sleek paintwork and touch screen controls, and despite being one of the larger scooters in the range, it is compact and easy to manoeuvre with a tight turning circle.

- › Top speed of 8mph*
- › Maximum range of 26 miles on a full battery charge*
- › Luxurious fully adjustable Captains seat with padded adjustable headrest
- › High / low speed selector to switch between 4mph (pavement only) and 8mph (road legal)
- › Excellent ground clearance 21cm (8") for a smooth ride over rough terrain
- › LED touch screen control panel
- › Fully functioning indicators and twin adjustable rear view mirrors for extra visibility
- › Advanced all round suspension ensures a smooth ride and optimal handling
- › Fitted with innovative Rhino 2 control system for a smooth, powerful ride
- › Comfortable lightweight and rotating seat which swivels 360° for easy transfers on and off the scooter
- › Removable seat and pneumatic tiller adjustment for easy transportation
- › Sliding and reclining seat with adjustable armrests
- › Delta bars as standard
- › Tie down points for non-occupant transportation
- › Horn to alert others of your presence for added safety
- › Fitted with reflectors, illuminated controls and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark
- › Power / Self-Diagnostic and Hazard warning lights for added safety
- › Front removable storage basket as standard
- › Bespoke alloy wheels with black pneumatic tyres for a smoother ride
- › Rear anti tip wheels as standard
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Easy to remove 2 x 79 Ah airline friendly batteries
- › Includes an on board charging facility as standard
- › Emergency brake for added safety
- › Available in Anthracite Grey and Spirit Red
- › Optional extras - Rear basket and additional storage options
- › Weight capacity 180kg (28st)
- * Subject to conditions of use

Colour Options

- Anthracite Grey
- Spirit Red

“Gives you an incredibly smooth ride, even on rougher terrain”

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
ST5DBE	ST5D Scooter	144cm (57")	136cm (53.5")	68cm (27")	10°	180kg (28st)

ST6



The ST6 is a perfectly formed scooter with its sweeping curves and enviable design, and a high performance which allows the scooter to be used both on and off road. The scooter features extra ground clearance and oil damped suspension, both of which make your journey smoother and more comfortable. On top of all these features, this scooter has an excellent grade climb to ensure you can effortlessly go up hills.

- > Top speed of 8mph*
- > Maximum range of 31 miles on a full battery charge*
- > Offers superb comfort with a luxurious, reclining, angle adjustable and sliding seat with a padded adjustable headrest and lumbar support
- > Luxurious Captains seat with two tone cover
- > High / low speed selector to switch between 4mph (pavement only) and 8mph (road legal)
- > High ground clearance of 10cm (5") for a smooth ride over rough terrain
- > Fully functioning indicators and adjustable twin rear view mirrors for extra visibility
- > Oil damped suspension ensures a smooth and comfortable ride
- > Backlit LCD illuminated controls
- > Fitted with innovative Rhino 2 control system for a smooth, powerful ride
- > Comfortable lightweight and rotating seat which swivels 360° for easy transfers on and off the scooter
- > Removable seat and pneumatic tiller adjustment for easy transportation
- > Comfortable height adjustable seat with adjustable armrests
- > Delta bars as standard
- > Four tie down points for safe transportation (without occupant)
- > Horn to alert others of your presence for added safety
- > Fitted with reflectors and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark
- > Power / Self-Diagnostic and Hazard warning lights for added safety
- > Alloy wheels with black pneumatic tyres for a smoother ride
- > Rear anti tip wheels as standard
- > Anti roll back system safety device prevents scooter from rolling backwards on hills
- > Freewheel facility allows movement of scooter without switching the motor on
- > Easy to remove 2 x 79 Ah airline friendly batteries
- > Includes an on board charging facility as standard.
- > Emergency brake for added safety
- > Available in Silver and Spirit Red
- > Additional storage in the modesty panel and seat pocket
- > Wide range of optional extras including rear basket and additional storage options
- > Weight capacity 197kg (31st)
- * Subject to conditions of use

Colour Options

- Silver
- Spirit Red

"Excellent on and off road!"



Bariatric Model

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
ST6BE	ST6 Scooter	160cm (63")	128cm (50")	72cm (28")	9°	197kg (31st)

Scout Range



The Scout has many features which are only normally found on larger scooters. It's the first Drive mini scooter to include Delta Bars as standard, it also features the next generation Drive splitting mechanism.

- > Top speed of 4mph*
- > Maximum range of 10 miles (MS003 / 008), 13 miles (MS009) on full battery charge*
- > Comfortable lightweight seat which swivels 360° for easy transfers on and off the scooter
- > Removable battery pack with easy to carry handle and locking mechanism
- > On board and in house charging as standard
- > Automatic electromagnetic brake system
- > Height adjustable seat to suit user comfort requirements
- > Delta bars as standard
- > Easy to change body panels
- > Easy tiller adjustment to suit user requirements
- > Removable seat, fold down tiller and splits easily for storage and transportation
- > Stylish silver wheel hubs with black puncture proof non marking tyres
- > Width adjustable flip up armrests for added flexibility
- > Front removable storage basket as standard
- > Rear anti tip wheels as standard
- > Anti roll back system safety device prevents scooter from rolling backwards on hills
- > Freewheel facility allows movement of scooter without switching the motor on
- > Easy to remove battery pack fitted with 2 x 12 Ah or 2 x 20 Ah airline friendly batteries
- > Available in Spirit Red and Electric Blue
- > Weight capacity 130kg (20st)
- * Subject to conditions of use

Colour Options

- Spirit Red
- Electric Blue

"Compact in size, big on features"

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	WEIGHT	HEAVIEST COMPONENT	MAXIMUM WEIGHT CAPACITY
MS008	Scout 4 Wheel - 12 Ah	108cm (42.5")	91cm (36")	48cm (19")	8°	44kg (97lb)	16kg (35lb)	130kg (20st)
MS009	Scout 4 Wheel - 20 Ah	108cm (42.5")	91cm (36")	48cm (19")	8°	49kg (107lb)	16kg (35lb)	130kg (20st)
MS003	Scout 3 Wheel - 12 Ah	108cm (42.5")	91cm (36")	48cm (19")	8°	44kg (97lb)	16kg (35lb)	130kg (20st)



The Prism range of scooters are complete with exciting features which include a stylish and modern body design, a huge choice of standard colours and the facility to order any colour of your choice.

- › Top speed of 4mph*
- › Maximum range of 10 miles on a full battery charge*
- › Second battery pack extends the range up to 20 miles*
- › Wireless connector means the scooter can be split apart for transportation in seconds
- › A huge choice of colours
- › Stylish and appealing body design with easy to change body panels
- › Comfortable lightweight seat which swivels 360° for easy transfers on and off the scooter
- › Includes an on or off board charging facility as standard
- › Intelligent electromagnetic brake system
- › Comfortable height adjustable seat with width adjustable armrests
- › Removable seat, fold down tiller and an integral lifting handle for easy transportation
- › Stylish silver wheel hubs with black puncture proof non marking tyres
- › Rear anti tip wheels as standard
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Front removable storage basket as standard
- › Delta bars available as an optional extra
- › Easy to remove battery pack fitted with 2 x 12 Ah or 2 x 20 Ah airline friendly batteries
- › Weight capacity 115kg (18st)
- Subject to conditions of use

Colour Options



*Only available on the Prism 3 wheel model

“Optional battery pack extends the maximum range up to 20 miles”



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	WEIGHT	HEAVIEST COMPONENT	MAXIMUM WEIGHT CAPACITY
MS017	Prism 3 Wheel	95cm (37")	87cm (34")	47cm (18")	8°	42kg (93lb)	13kg (29lb)	115kg (18st)
MS018	Prism 4 Wheel	104cm (41")	88cm (35")	48cm (19")	8°	47kg (103lb)	17kg (38lb)	115kg (18st)



The Prism Sport has all the performance and reliability of the Prism plus additional features. These include larger wheels, a maximum range of 12 miles, a greater weight capacity and mud flaps.

- › Top speed of 4mph*
- › Maximum range of 12 miles on full battery charge*
- › Excellent wireless connector means it can be split apart in seconds for easy transportation
- › Improved styling and appealing body design
- › Comfortable Bosun seat which swivels 360° for easy transfers on and off the scooter
- › On board and in house charging as standard
- › Intelligent electromagnetic brake system
- › Removable seat, fold down tiller and integral lifting handles for easy transportation
- › Height adjustable seat
- › Width and angle adjustable flip up armrests for easy transfers
- › Tiller head with headlight and LED battery gauge for ease of use
- › Excellent mud flaps
- › Larger 22cm x 10cm (9" x 4") grey solid tyres provide greater ground clearance
- › Rear anti tip wheels as standard
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Front removable storage basket as standard
- › Delta bars available as optional extra
- › Easy to remove battery pack fitted with 2 x 20 Ah airline friendly batteries
- › Weight capacity 135kg (21st)
- * Subject to conditions of use

“The high performance version of the popular Prism 4”



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	WEIGHT	HEAVIEST COMPONENT	MAXIMUM WEIGHT CAPACITY
MS019	Prism Sport	104cm (41")	90cm (36")	55cm (22")	8°	52kg (113lb)	29kg (66lb)	135kg (21st)



Envoy 8+



Envoy 6



Envoy 4

The Envoy range of Scooters are excellent additions to the Drive range, with the performance usually associated with much larger models. These scooters enjoy a whole host of excellent features including all round suspension, LED headlights and longer range batteries.

Envoy 8+

- › Top speed of 8mph*
- › Maximum range of 30 miles*
- › Fully functioning indicators and mirror
- › Powerful 500W motor
- › Secondary cable brake
- › Sport style seat
- › Modern dashboard design
- › Black pneumatic tyres with large alloy wheels
- › One touch hydraulic tiller adjustment
- › Delta bars as standard
- › Recommended MK batteries: 2 x 12V 50Ah

Envoy 6

- › Top speed of 6mph*
- › Maximum range of 30 miles*
- › Fully functioning indicators and mirror
- › All round suspension
- › Powerful 470W motor
- › Stylish silver wheel hubs
- › Delta bars as standard
- › Black pneumatic tyres with large alloy wheels
- › LED front headlight
- › Recommended MK batteries: 2 x 12V 50Ah

Envoy 4

- › Top speed of 4mph*
- › Maximum range of 30 miles*
- › Fully functioning mirror
- › Powerful 350W motor
- › All round suspension
- › Stylish silver wheel hubs
- › Delta bars as standard
- › Black pneumatic tyres with large alloy wheels
- › LED front headlight
- › Recommended MK batteries: 2 x 12V 50Ah
- * Subject to conditions of use

Colour Options

- Spirit Red
- Electric Blue
- White (4 & 8 Only)
- Silver (4 & 8 Only)

“The Envoy brings a new stylish design and offers a huge 30 mile range”



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	WEIGHT	MAXIMUM WEIGHT CAPACITY
MS050	Envoy 4	121cm (48")	102cm (40")	60cm (24")	8°	94kg (207lb)	160kg (25st)
MS051	Envoy 6	121cm (48")	102cm (40")	60cm (24")	8°	94kg (207lb)	160kg (25st)
MS058	Envoy 8+	121cm (48")	122cm (48")	60cm (24")	8°	100kg (220lb)	160kg (25st)



Viper



The Viper Scooter offers modern styling and a whole host of excellent features. These include modern state of the art controls, one touch tiller adjustment and delta bars as standard. Fully independent front and rear suspension, large 14" stylish alloy wheels providing a 10cm (4") high ground clearance, and an ergonomic adjustable seat. The result is a smooth and comfortable driving experience.



- › Top Speed of 8mph*
- › Maximum range of 25 miles (50ah) or 30 miles (79ah)*
- › All round suspension ensures a smooth ride and optimal handling
- › 14" stylish alloy wheels with black pneumatic tyres for a smoother ride
- › State of the art controls
- › Comfortable rotating seat with sliding and reclining facility
- › One touch hydraulic tiller adjustment to suit user requirements and assist easy transportation
- › Delta bars as standard
- › Height adjustable headrest to suit user requirements
- › Armrests are angle adjustable and flip up for easy transfers
- › Twin rear view mirrors for extra visibility
- › Rear anti-tip wheels as standard
- › Front removable storage basket as standard
- › Freewheel facility allows movement of scooter without switching the motor on
- › Large seat back storage pocket for storing essentials
- › Additional hand operated cable brake for added safety
- › Anti-roll back system safety device prevents scooter from rolling backwards on hills
- › Recommended MK batteries: 2 x 12v 50ah or 2 x 12v 79ah**
- › Weight capacity 160kg (25st)

* Subject to conditions of use
** Batteries not included

Colour Options



Spirit Red

"Large wheels and all round suspension for a smoother ride"



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	WEIGHT	MAXIMUM WEIGHT CAPACITY
VIPERRD	VIPER	143cm (56")	132cm (52")	56cm (26")	8°	120kg (264lb)	160kg (25st)



The Cobra is a revolutionary scooter in terms of design, aesthetics and styling. It has been aerodynamically designed to give an outstandingly stable drive and maximise battery consumption. It's a remarkable scooter that makes a huge statement and sets a new benchmark with its striking dynamic design.

- › Top speed of 8mph*
- › Maximum range of 30 miles on full battery charge*
- › Comfortable rotating seat with sliding and reclining facility
- › Electromagnetic brake system and emergency brakes
- › Height adjustable headrest to suit user requirements
- › Armrests are angle and width adjustable and flip up for easy side transfers
- › All round suspension ensures a smooth ride and optimal handling
- › One touch hydraulic tiller adjustment to suit user requirements and assist easy transportation
- › Excellent turning circle for use in small areas and busy environments
- › LED lighting and high level indicators and front and rear lights for maximum visibility at dusk or when dark
- › Twin rear view mirrors for extra visibility
- › Stylish 'Blade' alloy wheels with black pneumatic tyres for a smoother ride
- › Delta bars as standard
- › Rear anti tip wheels as standard
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Available in Graphite Grey, Spirit Red and White
- › Recommended MK batteries: 2 x 12V 79Ah
- › Weight capacity 182kg (28st)
- Subject to conditions of use

Colour Options



"This revolutionary scooter will take your adventures to the next level!"



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
COBRA	Cobra	160cm (63")	135cm (53")	73cm (29")	10°	182kg (28st)





King Cobra

King
COBRA

Following the success of our Cobra scooter we've developed the new King Cobra. This offers an evolved dynamic design, generous spaciousness, a 31 mile range and a maximum weight limit of 32 stone.



- › Top speed of 8mph*
- › Maximum range of 32 miles on full battery charge*
- › Stylish & distinctive design
- › Comfortable rotating seat with sliding and reclining facility
- › Height adjustable headrest and angle adjustable armrests to suit user requirements
- › Electromagnetic brake system and emergency cable brakes
- › All round suspension ensures a smooth ride and optimal handling
- › Generous tiller adjustment to suit user requirements and assist easy transportation
- › User friendly controls with a modern digital LCD display
- › Automatic speed reduction when steering for added safety
- › Twin USB power socket
- › Fully functioning indicators and twin adjustable rear view mirrors for extra visibility
- › High level indicators for added safety
- › Fitted with reflectors and LED front, rear and brake lights to improve outdoor visibility at dusk or when dark
- › Stylish 'Blade' alloy wheels with black pneumatic tyres for a smoother ride
- › Delta bars as standard
- › Bottle/cup holder included
- › Rear anti tip wheels as standard
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Available in Graphite Grey
- › Optional rear pannier box
- › Recommended MK Gel batteries: 2 x 12V 96Ah
- › Weight capacity 200kg (32st)
- * Subject to conditions of use

Colour Options

- Graphite Grey



"Powerful performance and dynamic design defines this luxury scooter"

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
KCOBRAGG	King Cobra	170cm (67")	140cm (54")	73cm (29")	10°	200kg (32st)

Royale 3 & 4

ROYALE



The Drive Royale 3 and 4 are two of the finest mobility scooters available today, both models have an impressive performance; powered by 70Ah batteries they have an excellent maximum range of 32 miles. These stylish scooters also feature a modern illuminated LCD display, providing Speedo, Battery Gauge, Clock, Odometer, Temperature, Reverse Light and Intelligent Diagnostics.

- › Top speed of 8mph*
 - › Maximum range of 32 miles on full battery charge*
 - › Stylish modern design
 - › Comfortable rotating seat with sliding and reclining facility
 - › Height adjustable headrest to suit user requirements
 - › Armrests are angle and width adjustable and flip up for easy side transfers
 - › Automatic speed reduction when steering for added safety
 - › All round suspension ensures a smooth ride and optimal handling
 - › Tiller adjustment by lever to suit user requirements and assist easy transportation
 - › Digital dash display with user friendly controls
 - › High level indicators and lighting to improve outdoor visibility
- › at dusk or when dark
 - › Twin rear view mirrors for extra visibility
 - › Large 40cm (16") alloy wheels with black pneumatic tyres for a smoother ride
 - › T-Bar style handlebars
 - › Anti roll back system safety device prevents scooter from rolling backwards on hills
 - › Freewheel facility allows movement of scooter without switching the motor on
 - › Rear anti tip wheels as standard
 - › Lap belt included for added safety
 - › Optional rear pannier box
 - › Available in Red, White and Black
 - › Recommended MK batteries: 2 x 12V 79Ah
 - › Weight capacity 185kg (29st)
 - * Subject to conditions of use

Colour Options

- Red
- White
- Black

"The result is an exceptional driving experience"

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
HW001	Royale 3 Wheel	158cm (62")	135cm (53")	73cm (29")	10°	185kg (29st)
HW002	Royale 4 Wheel	158cm (62")	135cm (53")	73cm (29")	10°	185kg (29st)





Royale 4 Sport

ROYALE

The Royale range of scooters continues to evolve with the introduction of this state of the art Royale 4 Sport. This stunning scooter features a stylish canopy, powerful 1500W motor and a maximum range up to an impressive 37 miles, providing the perfect balance of luxury comfort and performance.

- › Top speed of 8mph*
 - › Maximum range of 28 miles (70Ah batts) or 37 miles (100Ah batts) on full battery charge*
 - › Stylish modern design with canopy
 - › Large 51cm (20") seat width accommodates larger users providing added comfort
 - › Comfortable rotating seat with sliding and reclining facility
 - › Height adjustable headrest to suit user requirements
 - › Armrests are width adjustable and flip up for easy side transfers
 - › Automatic speed reduction when steering for added safety
 - › All round suspension ensures a smooth ride and optimal handling
 - › Tiller adjustment by lever to suit user requirements
 - › Digital dash display with user friendly controls
 - › High level indicators and lighting to improve outdoor visibility at dusk or when dark
 - › Twin rear view mirrors for extra visibility
 - › Large 36cm (14") front and 48cm (19") rear alloy wheels with black pneumatic tyres for a smoother ride and excellent ground clearance
 - › Anti roll back system safety device prevents scooter from rolling backwards on hills
 - › Freewheel facility allows movement of scooter without switching the motor on
 - › Lap belt included for added safety
 - › Highly powerful 1500W motor
 - › Available in Bronze
 - › Recommended MK batteries: 2 x 12V 79Ah or 2 x 12V 100Ah
 - › Weight capacity 160kg (25st)
- * Subject to conditions of use

Colour Options

- Bronze

"The stylish canopy adds to this bold and dynamic design"



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
HW006	Royale 4 Sport	169cm (67")	173cm (68")	75cm (30")	10°	160kg (25st)



The Sport Rider is a modern design 8mph scooter that includes all the tried and tested benefits of the Drive range. Using technology borrowed from the motorcycle industry the all round suspension provides excellent ride quality and handling, making a cruise on the Sport Rider both exhilarating and comfortable.

- › Top speed of 8mph*
- › Maximum range of 31 miles on full battery charge*
- › Comfortable rotating sport style seat with sliding and reclining facility
- › Height adjustable headrest and width adjustable armrests to suit user requirements
- › Angle adjustable handle bars allow comfortable positioning for the user
- › Throttle lever with automatic safe-braking system
- › Emergency cable brake for added security
- › Telescopic front motorbike suspension for optimal handling
- › Excellent manoeuvrability and turning circle for use in small areas and busy environments
- › Illuminated controls feature a twin dial digital dash
- › Trip computer with odometer for measuring distance travelled
- › Easy to use LED controls
- › All round lighting, inc large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark
- › Twin adjustable rear view mirrors for extra visibility
- › Large accommodating footwell with anti-slip rubber mat
- › Motorbike style wheels with modern chrome mud guards
- › Rear bumper provides additional protection
- › Rear anti tip wheels as standard
- › Lap belt included for added safety
- › Large seat back storage pocket for keeping essentials safe when out and about
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Optional lockable and removable rear pannier box
- › Recommended MK batteries: 2 x 12V 79Ah or 2 x 12V 100Ah
- › Weight capacity 180kg (28st)

* Subject to conditions of use

Colour Options

- Black

“The stylish Sport Rider was inspired by some of the leading designs in the motorcycle industry”

CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
SR001BLK	Sport Rider	165cm (65")	112cm (44")	78cm (31")	10°	180kg (28st)





Easy Rider

Inspired by the success of the Sport Rider, the design team got to work on a new model that would continue to redefine the typical luxury scooter. The Easy Rider has a striking new design, stunning instrumentation and offers precise handling that must be experienced to be appreciated.



- › Top speed of 8mph*
- › Maximum range of 31 miles on full battery charge*
- › Comfortable rotating sport style seat with sliding and reclining facility
- › Height adjustable headrest and angle adjustable armrests to suit user requirements
- › Angle adjustable handle bars allows comfortable positioning for the user (requires tools)
- › Throttle lever with automatic safe-braking system
- › Emergency cable brake for added security
- › Telescopic front motorbike suspension for optimal handling
- › Illuminated controls feature a twin dial digital dash
- › Trip computer with odometer for measuring distance travelled
- › Easy to use controls
- › All round lighting, including large headlight, front/rear indicators and rear brake lights to improve outdoor visibility at dusk or when dark
- › Twin adjustable rear view mirrors for extra visibility
- › Large accommodating footwell with anti-slip rubber mat
- › Motorbike style wheels with modern chrome mud guards
- › Rear bumper provides additional protection
- › Rear anti tip wheels as standard
- › Lap belt included for added safety
- › Large seat back storage pocket for keeping essentials safe when out and about
- › Anti roll back system safety device prevents scooter from rolling backwards on hills
- › Freewheel facility allows movement of scooter without switching the motor on
- › Optional lockable and removable rear pannier box
- › Recommended MK batteries: 2 x 12V 79Ah or 2 x 12V 96Ah
- › Weight capacity 180kg (28st)

Colour Options

- Black

“State of the art performance and design”



CODE	DESCRIPTION	LENGTH	HEIGHT	WIDTH	MAXIMUM GRADIENT	MAXIMUM WEIGHT CAPACITY
SR003BLK	Easy Rider	166cm (65")	115cm (45")	76cm (30")	10°	180kg (28st)

Energi Powerchair



The Enigma Energi Lightweight Powerchairs have many excellent features normally found on more expensive chairs yet are surprisingly affordable. This Powerchair range provides great comfort with extra features as standard.

- › Maximum speed of 4mph*
- › Maximum range of 15 miles on a full battery charge*
- › Available in 41cm (16"), 46cm (18"), 51cm (20") and 56cm (22") seat widths
- › Designed for indoor use and limited outdoor use such as shopping centres
- › Rear wheel drive system for optimum manoeuvrability
- › Detachable height and width adjustable armrests with removable PU arm pads to suit comfort of user
- › Double cross brace for responsive manoeuvrability and added support
- › Half folding and removable backrest for easy storage and transportation
- › Tension and height adjustable backrest and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences
- › De-clutch operation to enable you to safely freewheel the chair where and when required
- › Padded upholstery with choice of four seat widths
- › 20cm (8") front castors and 32cm (12.5") rear wheels with pneumatic tyres provide durability and low maintenance
- › Comfortable seat cushion included as standard
- › Swing in / swing out detachable footrests for maximum comfort and convenience when transferring
- › PG Drives VR2 controller provides an easy and enjoyable driving experience
- › Powerful 320W motors
- › Tried and tested industry standard electronics
- › Detachable batteries and basket for easy storage
- › Built in storage bag for battery charger
- › On board charging (2 x 34Ah batteries)
- › Rear anti tip wheels as standard
- › Available in three colour options – Metallic Blue, Red or Black
- › This powerchair has been crash tested providing added safety when being transported in a vehicle
- › Various optional extras are available such as attendant controls, dual controls, kerb climbers and elevating footrests
- › Maximum weight capacity 135 – 180kg (21 - 28st)
- * Subject to conditions of use

Colour Options

- Metallic Blue
- Red (16", 18" & 20")
- Black



Crash Tested

CODE	DESCRIPTION	CODE	DESCRIPTION
PCA16	Energi 41cm (16") Powerchair	PCKC006	Kerb Climber 51cm (20")
PCA18	Energi 46cm (18") Powerchair	PCKC007	Kerb Climber 56cm (22")
PCA20	Energi 51cm (20") Powerchair	Z15982-01	Elevating Footrest - Right
PCA22	Energi 56cm (22") Powerchair	Z15983-01	Elevating Footrest - Left
PCKC004	Kerb Climber 41cm (16")	Z21076PG	Attendant Control
PCKC005	Kerb Climber 46cm (18")	Z21065 & Z21071	Dual Control (Bracket & Joystick)

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
PCA16	61cm (24")	118cm (46")	94cm (37")	41cm (16")	40cm (16")	53cm (21")	56kg (123lb)	135kg (21st)
PCA18	66cm (26")	118cm (46")	94cm (37")	46cm (18")	40cm (16")	53cm (21")	56kg (123lb)	135kg (21st)
PCA20	71cm (28")	118cm (46")	94cm (37")	51cm (20")	40cm (16")	53cm (21")	56kg (123lb)	135kg (21st)
PCA22	76cm (30")	118cm (46")	94cm (37")	56cm (22")	40cm (16")	53cm (21")	58kg (128lb)	180kg (28st)

Energi+ Powerchair



The Enigma Energi+ features stylish oval tubing and has an electric blue finish with vivid artwork. It also features a handy slide out battery tray.

- › Maximum speed of 4mph*
- › Maximum range of 15 miles on a full battery charge*
- › Available in 46cm (18") seat width
- › Stylish oval tubing with an Electric Blue finish
- › Excellent slide out battery tray for easy charging
- › Designed for indoor use and limited outdoor use such as shopping centres
- › Rear wheel drive system for optimum manoeuvrability
- › Detachable height and width adjustable armrests with removable PU arm pads to suit the comfort of the user
- › Double cross brace for responsive manoeuvrability and added support
- › Half folding and removable backrest for easy storage and transportation
- › Tension and height adjustable backrest and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences
- › De-clutch operation to enable you to safely freewheel the chair where and when required
- › Padded upholstery for additional comfort and support
- › 20cm (8") front castors and 32cm (12.5") rear wheels with solid puncture proof tyres provide durability and low maintenance
- › Comfortable seat cushion included as standard
- › Swing in / swing out detachable footrests for maximum comfort and convenience when transferring
- › PG Drives VR2 controller provides an easy and enjoyable driving experience
- › Powerful 320W motors
- › Tried and tested industry standard electronics
- › Detachable batteries and basket for easy storage
- › Built in storage bag for battery charger
- › On board charging (2 x 34Ah batteries)
- › Rear anti tip wheels as standard
- › Various optional extras are available such as attendant controls, dual controls, kerb climbers and elevating footrests
- › Maximum weight capacity 135kg (21st)
- * Subject to conditions of use



Slide out battery tray

CODE	DESCRIPTION
PCP18	Energi+ 46cm (18") Powerchair
PCKC005	Kerb Climber 46cm (18")
ELR002R	Elevating Footrest - Right
ELR002L	Elevating Footrest - Left
Z21076PG	Attendant Control
Z21065 & Z21071	Dual Control (Bracket & Joystick)

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT LENGTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
PCP18	66cm (26")	114cm (45")	96cm (38")	46cm (18")	40cm (16")	47cm (19")	60kg (132lb)	135kg (21st)

Cobalt



The Cobalt Powerchair is a compact, lightweight and highly transportable chair which is simple and easy to use due to its user friendly programmable controller. This powerchair splits into 3 sections for easy handling and storage when not in use.

- › Maximum speed of 4mph*
- › Maximum range of 8 miles on a full battery charge*
- › Available in 46cm (18") seat width
- › Designed for indoor and outdoor use
- › Compact and manoeuvrable for use in tight spaces and busy environments
- › Splits easily for storage and transportation
- › Rear wheel drive system for optimum manoeuvrability
- › Detachable width adjustable armrests with removable arm pads to suit comfort of user
- › Folding backrest for easy storage and transportation
- › Height adjustable seat and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences
- › Swivel seat facilitates easy transfers on and off the chair
- › De-clutch operation to enable you to safely freewheel the chair where and when required
- › Padded upholstery for additional comfort and support
- › 20cm (8") front castors and 23cm (9") rear wheels with solid puncture proof tyres provide durability and low maintenance
- › Height adjustable, flip up footplate to suit users leg length and assist with easy transfers
- › Dynamic Shark programmable controller provides an easy and enjoyable driving experience
- › Rear anti tip wheels as standard
- › Seat belt included for additional safety and security
- › Available in two colour options – Metallic Blue or Metallic Red
- › Maximum weight capacity 113kg (18st)
- * Subject to conditions of use

Colour Options

- Metallic Blue
- Metallic Red

"Highly compact and ultra-manoeuverable"



CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	WEIGHT	WEIGHT CAPACITY
COBALTX16	61cm (24")	98cm (38.5")	84 - 91cm (33 - 36")	46cm (18")	38cm (15")	51 - 58cm (20 - 23")	55kg (121lb)	113kg (18st)

Cirrus



The Cirrus Powerchair is foldable and lightweight making it ideal for transporting. It comes in an attractive, lightweight, carbon steel frame with a silver finish which is easy to maintain.

- › Maximum speed of 4mph*
- › Maximum range of 15 miles on a full battery charge*
- › Available in 46cm (18") seat width
- › Attractive steel cross-brace frame
- › Designed for indoor and limited outdoor use such as shopping centres
- › Rear wheel drive system
- › Flip up, height adjustable and removable padded armrests to suit comfort of user
- › Comes with 1.5" seat cushion and front pocket for storing essentials
- › Comes with back cushion which includes a larger pocket for additional storage
- › Tension adjustable backrest for added comfort
- › Adjustable length controller mount
- › Detachable and height adjustable back post and push handles
- › Foldable lightweight frame for easy storage and transportation
- › Electromagnetic brakes and emergency brakes
- › Programmable controller
- › De-clutch operation to enable you to safely freewheel the chair where and when required
- › Padded nylon upholstery for additional comfort and support
- › Solid front 8" castors with 12.5" pneumatic rear tyres provide durability and low maintenance
- › Swing away, removable footrests to assist with easy transfers
- › Height adjustable legrests suit users leg length
- › Detachable anti tip wheels as standard
- › One touch lap belt and calf strap included for additional safety and security
- › Off board charging with folding battery tray (2 x 34Ah Batteries)
- › Available in Silver
- › Maximum weight capacity 136kg (21st)
- * Subject to conditions of use

"Robust yet lightweight steel frame"



Adjustable length controller mount

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	WEIGHT	WEIGHT CAPACITY
PCC18SIL	64cm (25")	108cm (42.5")	92 - 100cm (36 - 39")	45cm (18")	40cm (16")	65kg (143lb)	136kg (21st)

Titan



The Titan Powerchair is perfect for users on the go as it can be disassembled, without tools, and loaded easily into a car for easy transportation. Its ultra-manoeuvrable and compact size makes it ideal for around the home and in busy shopping centres.

- › Maximum speed of 4mph*
 - › Maximum range of 15 miles on a full battery charge*
 - › Available in 46cm (18") seat width
 - › Designed for indoor and outdoor use
 - › Front wheel drive system for optimum manoeuvrability
 - › Detachable padded and adjustable armrests to suit comfort of user
 - › Adjustable length controller mount
 - › Folding backrest for easy storage and transportation
 - › Height adjustable and reclining seat to suit comfort preferences
 - › Swivel seat facilitates easy transfers on and off the chair
 - › Height adjustable and removable headrest to suit user requirements
 - › De-clutch operation to enable you to safely freewheel the chair where and when required
 - › Padded upholstery for additional comfort and support
 - › Solid non-marking tyres provide durability and low maintenance
 - › Height adjustable, flip up footplate to suit users leg length and assist with comfort
 - › Dynamic Shark programmable controller provides an easy and enjoyable driving experience
 - › Front anti tip wheels as standard
 - › Seat belt included for additional safety and security
 - › Available in Spirit Red
 - › Maximum weight capacity 136kg (21st)
- * Subject to conditions of use

"Featuring the new front wheel drive system for optimum manoeuvrability"



Easy length adjustable controller

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	WEIGHT	WEIGHT CAPACITY
TITAN18LC	63.5cm (25")	102cm (40")	112cm (44")	46cm (18")	46cm (18")	72kg (159lb)	136kg (21st)

Sunfire Plus GT



The Sunfire Plus GT Powerchair offers a traditional rear wheel drive configuration with modern styling and an automotive style Captains seat for additional comfort.

- › Maximum speed of 4mph*
 - › Maximum range of 20 miles on a full battery charge*
 - › Available in 46cm (18") seat width
 - › Mid-size rear wheel drive for optimum manoeuvrability
 - › Reclining padded seat provides comfort for user
 - › Articulating front beam enabling easy use over gradients and uneven surfaces
 - › Swivel seat facilitates easy transfers on and off the chair
 - › Height adjustable and removable headrest to suit user requirements
 - › Flip up armrests allow improved access to tables, desks, etc and assist to facilitate easy transfers
 - › Designed for indoor and outdoor use for maximum flexibility
 - › Moulded PU armrests that are width, height and angle adjustable to suit comfort preferences
 - › De-clutch operation to enable you to safely freewheel the chair where and when required
 - › 20cm (8") front and 25cm (10") rear castors with grey solid puncture proof tyres provide durability and low maintenance
 - › Angle and height adjustable flip up footplate to facilitate transfers
 - › PG Drives VR2 controller provides an easy and enjoyable driving experience
 - › Powerful 350W motors
 - › On or off board charging (2 x 34Ah batteries)
 - › Rear anti tip wheels as standard
 - › Seat belt included for additional safety and security
 - › Maximum weight capacity 136kg (21st)
- * Subject to conditions of use

"Sleek and robust design"



CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	BACKREST HEIGHT (EXC HEADREST)	WEIGHT	WEIGHT CAPACITY
PCMM07	66cm (26")	104cm (41")	114cm (45")	46cm (18")	42cm (17")	43cm (17")	96kg (201lb)	136kg (21st)

Image EC



The Image EC Powerchair is ideal for use indoors due to its highly manoeuvrable mid wheel drive configuration and small footprint.

- › Maximum speed of 4mph*
- › Maximum range of 20 miles on a full battery charge*
- › Available in 46cm (18") seat width
- › Mid wheel drive for optimum manoeuvrability
- › Height adjustable reclining padded seat provides comfort for user
- › Six wheels provide added stability when using the powerchair
- › Swivel seat facilitates easy transfers on and off the chair
- › Height adjustable and removable headrest to suit user requirements
- › Flip up armrests allow improved access to tables, desks, etc and assist to facilitate easy transfers
- › Designed for indoor and outdoor use for maximum flexibility
- › Moulded PU armrests that are width, height and angle adjustable to suit comfort preferences
- › De-clutch operation to enable you to safely freewheel the chair where and when required
- › Padded upholstery to provide optimum comfort for the user
- › 25cm (10") middle wheels and 15cm (6") front and rear castors with solid puncture proof tyres provide durability and are low maintenance
- › Angle and height adjustable flip up footplate to facilitate transfers
- › PG Drives VR2 controller provides an easy and enjoyable driving experience
- › Powerful 200W motors
- › On or off board charging (2 x 34Ah batteries)
- › Rear anti tip wheels as standard
- › Seat belt included for additional safety and security
- › Maximum weight capacity 136kg (21st)
- * Subject to conditions of use

"Highly manoeuvrability and agile"



Sunfire General



The Sunfire General features a reclining padded seat and headrest for ultimate comfort, and is upholstered in an attractive two tone fabric. The innovative one touch armrests, swivelling seat and flip up footplate all assist in providing easy transfers on and off the powerchair.

- › Maximum speed of 4mph*
- › Maximum range of 20 miles on a full battery charge*
- › Available in 46cm (18") seat width
- › Larger rear wheel drive for optimum manoeuvrability
- › Large height adjustable reclining padded seat provides comfort for user
- › Innovative one touch height adjustable armrests
- › Articulating front beam enabling easy use over gradients and uneven surfaces
- › Swivel seat facilitates easy transfers on and off the chair
- › Height adjustable and removable headrest to suit user requirements
- › Flip up armrests allow improved access to tables, desks, etc and assist to facilitate easy transfers
- › Designed for indoor and outdoor use for maximum flexibility
- › Moulded PU armrests that are width, height and angle adjustable to suit comfort preferences
- › Single lever de-clutch operation to enable you to safely freewheel the chair where and when required
- › 23cm (9") front and 33cm (13") rear castors with solid puncture proof tyres provide durability and low maintenance
- › Angle and height adjustable flip up footplate for comfort
- › PG Drives VR2 controller provides an easy and enjoyable driving experience
- › Powerful 450W motors
- › On or off board charging (2 x 34Ah batteries)
- › Rear anti tip wheels as standard
- › Seat belt included for additional safety and security
- › Maximum weight capacity 180kg (28st)
- * Subject to conditions of use

"Comfortable & accommodating large seat"



CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	BACKREST HEIGHT (EXC HEADREST)	WEIGHT	WEIGHT CAPACITY
PCMM08	63cm (25")	106cm (42")	110cm (43")	46cm (18")	42cm (16")	43cm (17")	87kg (192lb)	136kg (21st)

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	BACKREST HEIGHT (EXC HEADREST)	WEIGHT	WEIGHT CAPACITY
PCMM512	64cm (25")	103cm (41")	112cm (44")	51cm (20")	43cm (17")	50cm (20")	100kg (220lb)	180kg (28st)



Bariatric Model

Seren with Captain's Seat



The nimble Seren is an ultra-manoeuvrable powerchair which can be used both indoors and outdoors. Compact and agile, the Seren moves effortlessly around tight corners, so it's ideal in areas where space is limited. The powerchair uses a dynamic control system to give a superb non-stop driving performance.

- › Maximum speed of 4mph*
 - › Maximum range of 12 miles on a full battery charge*
 - › Available in 46cm (18") seat width
 - › Padded Captain's seat and adjustable headrest provide additional comfort and support
 - › Designed for indoor use and limited outdoor use such as shopping centres
 - › Compact, agile and ultra manoeuvrable for use in tight spaces and busy environments
 - › Front-wheel drive system
 - › Height and width adjustable armrests with PU arm pads to suit comfort of user
 - › Flip up armrests facilitate side transfers
 - › Front and rear suspension to ensure a smooth and safe ride
 - › Swivel seat facilitates transfers on and off the chair
 - › Aluminium 25cm (10") middle wheels and 15cm (6") front and back castors with solid puncture proof tyres
 - › Single centre mount footplate
 - › Dynamic Shark controller
 - › Powerful 420W motors
 - › Off board charging (2 x 36Ah batteries)
 - › Fitted with a seat belt for extra security
 - › Maximum weight capacity 135kg (21st)
- * Subject to conditions of use

"Superb non-stop driving performance"



Front and rear suspension to ensure a smooth and safe ride

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	BACK HEIGHT	WEIGHT	WEIGHT CAPACITY
SERENBE	61cm (24")	100cm (39")	105cm (41")	46cm (18")	45cm/47cm/50cm (18"/19"/20")	43cm/46cm/48cm/51cm (17"/18"/19"/20")	53cm (21")	80kg (176lb)	135kg (21st)

Seren with Clinical / Adjustable Seat



The nimble Seren is an ultra-manoeuvrable powerchair which can be used both indoors and outdoors. Compact and agile, the Seren moves effortlessly around tight corners, so it's ideal in areas where space is limited. The powerchair uses a dynamic control system to give a superb non-stop driving performance.

- › Maximum speed of 4mph*
 - › Maximum range of 12 miles on a full battery charge*
 - › Available in 41cm (16") and 46cm (18") seat widths
 - › Fold down back mechanism for easy transportation and storage
 - › Designed for indoor use and limited outdoor use such as shopping centres
 - › Compact, agile and ultra manoeuvrable for use in tight spaces and busy environments
 - › Front-wheel drive system
 - › Height and width adjustable armrests with PU arm pads to suit comfort of user
 - › Flip up armrests facilitate side transfers
 - › Tension adjustable backrest to suit comfort preferences
 - › Width, height and depth adjustable seat with cushion and adjustable back width to suit various user sizes
 - › Front and rear suspension to ensure a smooth and safe ride
 - › Swivel seat facilitates transfers on and off the chair
 - › Aluminium 25cm (10") middle wheels and 15cm (6") front and back castors with solid puncture proof tyres
 - › Single centre mount footplate
 - › Dynamic Shark controller
 - › Powerful 420W motors
 - › Off board charging (2 x 36Ah batteries)
 - › Fitted with a seat belt for extra security
 - › Maximum weight capacity 135kg (21st)
- * Subject to conditions of use

Additional options available –

- › Attendant Control
- › Dual Control
- › Swing Away Joystick Holder
- › Elevating Legrest
- › Standard Legrest

"Compact and agile"



Dynamic Shark Controller

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	BACK HEIGHT	WEIGHT	WEIGHT CAPACITY
SERENBE	61cm (24")	100cm (39")	105cm (41")	41cm/46cm (16"/18")	45cm/47cm/50cm (18"/19"/20")	43cm/46cm/48cm/51cm (17"/18"/19"/20")	43cm (17")	80kg (176lb)	135kg (21st)

Multego



The Multego is a dynamic mid-range powerchair designed for simple yet powerful mobility. This powerchair features a height and tension adjustable backrest, kerb climber and seat belt.

- › Maximum speed of 4mph*
- › Maximum range of 16 miles on a full battery charge*
- › Available in 41cm (16") and 46cm (18") seat widths
- › Side panels provide a snug fit, increased safety and protection from draughts when using the chair outdoors
- › Rear-wheel drive system
- › Fold down back mechanism for easy storage and transportation
- › Width, depth and angle adjustable seat with cushion and adjustable back width to suit various user sizes
- › Legrests, armrests and backrest can also be adjusted independently of each other for maximum comfort and support
- › Height and width adjustable armrests with PU arm pads to suit comfort of user
- › Swing in / swing out detachable and height adjustable legrest for maximum comfort and convenience when transferring
- › Choice of three seat depths to suit user preferences
- › Aluminium rear wheels with 32cm (13") solid puncture proof tyres
- › Front castors with 20cm (8") solid puncture proof tyres
- › Dynamic Shark Series control for an easy and enjoyable driving experience
- › Powerful 420W motors
- › Off board charging (2 x 36Ah batteries)
- › Fitted with a seat belt for extra security
- › Fitted with a kerb climber as standard
- › Maximum weight capacity 136kg (21st)
- * Subject to conditions of use

Additional options available –

- › Attendant Control
- › Dual Control
- › Swing Away Joystick Holder

“Fully adjustable for maximum comfort and support”



Swing in / swing out armrests



Crash Tested

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	BACK HEIGHT	WEIGHT	WEIGHT CAPACITY
MULTEGOBE	58cm (23")	104cm (41")	92cm (36")	41cm/46cm (16"/18")	41cm/43cm/46cm (16"/17"/18")	41cm - 51cm (16" - 20")	74kg (163lb)	136kg (21st)

Volt



The Volt Powerchair has a fixed aluminium frame which is fitted with solid tyres to ensure you'll never have a puncture, and a manual seat tilt and back recliner as standard to help prevent discomfort and pressure sores.

- › Maximum speed of 4mph*
- › Maximum range of 16 miles on a full battery charge*
- › Available in 41cm (16"), 46cm (18") and 50cm (20") seat widths
- › Calf rest provides added comfort and support for the user
- › Side panels provide a snug fit, increased safety and protection from draughts when using the chair outdoors
- › Attendant handlebar allows carers to manoeuvre the chair when transporting users
- › Small carry pouch is positioned just below the front of the seat for handy access to essentials such as keys
- › Tension adjustable backrest and length adjustment on armrests allow chair to be easily adapted to suit comfort preferences
- › Detachable height and width adjustable armrests with PU arm pads to suit comfort of user
- › Single lever de-clutch operation to enable you to safely freewheel the chair where and when required
- › Width and depth adjustable seat (requires tools)
- › Swing in / swing out detachable and angle adjustable legrests
- › Choice of three seat depths to suit user preferences
- › Aluminium rear wheels with 35cm (14") solid puncture proof tyres
- › Front castors with 22cm (9") solid puncture proof tyres
- › Dynamic Shark Series Controller for an easy and enjoyable driving experience
- › Powerful 450W motors
- › On or off board charging (2 x 50Ah batteries)
- › Seat cushion included for added comfort
- › Fitted with a kerb climber as standard.
- › Fitted with a seat belt for extra security
- › Maximum weight capacity 135kg (21st)
- * Subject to conditions of use

Additional options available –

- › Attendant Control
- › Dual Control
- › Elevating Legrest

“Comfortable & accommodating large seat”



Swing in / swing out legrests and kerb climber as standard



Crash Tested

CODE	WIDTH	LENGTH	HEIGHT	SEAT WIDTH	SEAT DEPTH	SEAT HEIGHT	BACK HEIGHT	WEIGHT	WEIGHT CAPACITY
VOLTBE	62cm (24")	103cm (40.5")	102cm (40")	41cm/46cm/50cm (16"/18"/20")	41cm/43cm/46cm (16"/17"/18")	48.5cm (19")	53cm (21")	71kg (157lb)	135kg (21st)

Powerstroll S Drive



This excellent Powerstroll will turn most manual wheelchairs into carer controlled powered wheelchairs. It takes the strain out of pushing a heavy wheelchair and makes easy work of hills, ramps and uneven surfaces. It can be fitted or removed easily in seconds.

- › Allows speed of up to 4mph*
- › Travels up to 10 miles (16 km) on a full battery charge*
- › The powerstroll can easily be fitted or removed within seconds
- › The excellent P G Drives S Drive controller features a fully programmable and comprehensive fault management system
- › New push and slide battery removal system
- › It is very portable, making transportation simple
- › Fits most standard wheelchairs up to 50cm (20") seat width
- › Additional parts available for wider wheelchairs
- › Hand control with variable speed dial, forward/reverse function, illuminated battery gauge and fingertip trigger control
- › Comes complete with a carry holdall with shoulder strap
- › Simple to manoeuvre and reverse
- › The powerstroll has a powerful 150 watt (200 watt-HD), 24V motor, capable of coping with most terrains
- › Maintenance free airline friendly 24v 12Ah valve regulated sealed lead acid battery pack
- › 2A charger as standard for off board charging
- › Weight capacity 135kg (21st), this includes the wheelchair and user combined
- › Heavy Duty model - maximum weight capacity 180kg (28st)
- * Subject to conditions of use

"Features our unique patented one click docking station"



Bariatric Model

CODE	TOTAL WEIGHT	WEIGHT WITHOUT BATTERY	WEIGHT BATTERY PACK ONLY	MAXIMUM GRADIENT	MAXIMUM RANGE	WEIGHT CAPACITY
PWCPP009	19kg (42lb)	10kg (22lb)	9kg (20lb)	8°	10 Miles (16km)	135kg (21st)
PWCPP009HD	19kg (42lb)	10kg (22lb)	9kg (20lb)	8°	10 Miles (16km)	180kg (28st)

Lightweight Powerstroll



This dual wheel Lightweight Powerstroll can easily be fitted and removed within seconds, and adapts to fit most wheelchair sizes. This powerstroll will turn most manual wheelchairs into carer controlled powered wheelchairs and take the strain out of pushing a heavy wheelchair.

- › Allows speed of up to 3mph*
- › Travels up to 10 miles (16km) on a full battery charge*
- › Fits most wheelchairs with seat widths from 35cm (14") to 50cm (20") (including narrow track wheelchairs such as the Lightweight Aluminium)
- › Fits 2.2cm - 2.5cm (7/8" - 1") diameter tubing
- › Uses existing wheelchair handgrips
- › Lightweight at 14kg (31lb) and less than 8kg (17lb) with battery pack removed for easy transportation
- › Trigger handset with variable speed dial and direction button to suit the attendant and user requirements
- › The powerstroll can be easily fitted and removed
- › Twin 'Drive' wheels with 20cm (8") black non marking PU tyres
- › Removable fabric battery pack for easy storage and transportation
- › Uses 12v 20Ah battery (airline friendly)
- › Lightweight 100w 12v DC motor making the unit very portable
- › Easy to use lift up strap to increase manoeuvrability

"Lightweight, portable, and is available at a great price"



CODE	TOTAL WEIGHT	WEIGHT WITHOUT BATTERY	WEIGHT BATTERY PACK ONLY	MAXIMUM GRADIENT	MAXIMUM RANGE	WEIGHT CAPACITY
PWCPP010	14kg (31lb)	8kg (17lb)	7kg (14lb)	6°	10 Miles (16km)	115kg (18st)



The unique U-Drive system converts the XS aluminium self propelled wheelchair, into a user controlled electric powerchair. This revolutionary new product is an excellent addition to the Powerstroll range from Drive Medical.

- > Allows speed of up to 4mph*
- > Travels up to 6 miles (10km) on a full battery charge*
- > The U-Drive can easily be fitted or removed
- > The excellent P G Drives S-Drive Joystick provides true user control
- > New push and slide battery removal system to facilitate easy charging
- > Very portable, making transportation simple
- > When detached the U-Drive and wheelchair will fit into most car boots for easy transportation
- > Simple to manoeuvre and reverse making the U-Drive very user friendly
- > Maintenance free airline friendly 24v 12Ah valve regulated sealed lead acid battery pack
- > On or off board charging
- > 2A charger provided as standard for off board charging
- > Dual control option is available to allow both the attendant and the user to control the wheelchair
- * Subject to conditions of use

"The U-Drive offers the best of both worlds"



CODE	TOTAL WEIGHT	WEIGHT WITHOUT BATTERY	WEIGHT BATTERY PACK ONLY	MAXIMUM GRADIENT	MAXIMUM RANGE	WEIGHT CAPACITY
U-DRIVE	31kg (67lb)	22kg (49lb)	9kg (19lb)	6°	6 Miles (10km)	115kg (18st)



Leaflet Code: LL218

Comfort, Style, Performance & Reliability

We understand the true value of innovation within this market place; and know that our powered mobility development will have a direct influence on improving an individual's life.

www.drivedevilbiss.co.uk

All sizes and weights stated are nominal. All errors and omissions excepted.
Due to the limitations of photography and the printing process colours shown may not be 100% accurate.