



eFOLDi[®]

eFOLDi Power Chair

The eFOLDi Power Chair is the lightest, most compact power wheelchair in its class. Weighing only 15kg (without batteries) and fully airline compliant, it also has an easy to use manual set-up and fold operation.

Making life simple with both user and attendant control, the eFOLDi Power Chair is the perfect travel companion.

Enjoy freedom, friendship and fun!



eFOLDi[®]

Features

Speed	4mph (6kmh)
Range	8miles / 13 Km
Carry capacity	135Kg / 297 lbs / 21 Stone
Climbing ability	1:10 Gradient (6 degrees)
Folded dimensions	58 x 34 x 92 cm / 22 x 13 x 36 inches
Brake	Electro Magnetic
Tyre type	Front solid / Rear solid
Battery	5.2AH Li-ion x 2
Motor	200W x2 24V
Controller	User joystick with additional attendant control