



Invacare® | **Matrix® Flo-tech™ Lite Visco**

Low profile visco elastic memory foam cushion offering exceptional comfort


matrix®
SEATING SERIES

INVACARE
Yes, you can.

The *Invacare **Matrx Flo-tech Lite Visco*** cushion is made of high-density visco elastic foam, allowing it to conform to the individual's shape. This affords exceptional levels of comfort and increases and redistributes pressure. This low profile cushion format is ideal for both chair and wheelchair use, facilitating side transfers.

Low memory foam designed to maximise body contact area

The visco-elastic foam conforms to the individual's body, maximising the contact area, optimising pressure distribution and delivering exceptional levels of comfort.



Features and Options



Ingress resistant cover

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base with carry handle for easy transportation.



Available in 64 sizes

The *Invacare **Matrx Flo-tech Lite Visco*** is available in a range of 64 sizes, from 10-20" widths and 10-21" depths. (See prescription chart for details).



Three foam densities available

Available in three densities to accommodate a range of client weights – Soft (up to 32 kg), Medium (up to 51 kg) and Firm (up to 100 kg).

Contoured seating surface withommel and side support

Offers mild lateral support and pelvic stability to enhance the individual's posture.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk

Various sizes	Height	Total product weight	Max. user weight ¹	Max. cleaning temperature	Warranty	Manufactured to European Standards: Interior Foam: BS EN 1021-1, BS EN 1021-2, BS 5852: 1982 Crib 5
Up to 64 sizes	78 mm	1.26 kg	Soft: 32kg (5 stone) Medium: 51kg (8 stone) Firm: 100kg (16 stone)	80°	3 years	

1. Guide maximum user weight limits are set assuming that the cushion selected is of appropriate width and depth for the intended user/s.

Colours


Black Smooth



Black wipedown



Black spacer fabric



Baltic Blue

