



Colours Available:



# Just Fold And Go!

## The Adventure of a Lifetime!

This compact power chair has plenty of storage and folds in just a few steps to make travel effortless. When your plane touches down at your destination, you can just unfold and go!

## Features



Compact  
Lightweight Design



Suspension



Micon Control System

Specifications	
<b>Chair</b>	
Weight Capacity	17 st 12 lbs (113 kg)
Overall Length <sup>3</sup>	1016 mm (40")
Overall Width <sup>3</sup>	597 mm (23.5")
Ground Clearance <sup>3</sup>	44mm
Suspension	Front
Colour(s)	Silver
<b>Power</b>	
Turning Radius <sup>3</sup>	616 mm (24.25")
Max Speed <sup>1</sup>	3.6 mph (5.8 km/h)
Drivetrain	2 x Motors
Dual Breaking System	Electronic Regenrative, Disc Park Brake
<b>Battery</b>	
Battery Requirements <sup>3</sup>	18.2 Ah Lithium-Ion Battery
Battery Weight <sup>4</sup>	3.4 kg
Battery Charger	2 - Amp, Off-Board
Total Weight Without Batteries <sup>2, 4</sup>	24 kg

### NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.