

INDIANA

Rise & Recliner /
Standard & Petite

**Restwell.**

*“Standard
and petite
options”*



CE

The Indiana Rise Recliner is both stylish and comfortable, providing excellent value for money. Available in standard or petite size options.

www.restwellchairs.co.uk

driveTM
www.drivemedical.co.uk

The Indiana Riser Recliner is upholstered in a soft feel Chenille fabric which provides added luxury and comfort. The riser has a modern design which will suit any décor, and is available in four attractive colours.

Superb Comfort, Style & Performance...



- » Soft feel Chenille fabric for luxurious comfort.
 - » Attractive and modern design to suit any decor.
 - » Single motor activates a distinct lift, rest and recline action.
 - » It will rise and tilt users to their feet.
 - » Excellent snooze position for supreme comfort and relaxation.
 - » Fully padded chaise for increased comfort.
 - » Quiet motors provide a smooth and discreet operation.
 - » Ergonomically contoured seat and back and padded upholstery for added comfort and support.
 - » Convenient side pocket for easy storage.
 - » Quality assured, meeting all relevant fire retardancy standards.
 - » Comes with a simple to operate two button handset.
 - » Ideal for individuals who have difficulty standing up and sitting down due to restricted mobility e.g. those with MS, Arthritis, Rheumatism, circulation problems.
 - » Available in Petite size for proportional fit to people 1.7m (5'6") tall or less.
 - » Available in 4 colour options -Berry, Biscuit, Charcoal and Mushroom.
 - » 2 year warranty*.
 - » Maximum weight capacity 150kg (23st).
- * Warranty conditions apply, for more information please see user manual.
** Batteries not included.

INDIANA

CODE: CLR19C (Standard)
CLR29C (Petite)

Max Weight Capacity:
150kg (23st)




 **Height:**
108cm (42½") - ST
103cm (40½") - PE

 **Width:**
81cm (32") - ST
77cm (30¼") - PE

 **Depth:**
96cm (37¾") - ST
88cm (34½") - PE

 **Seat Depth:**
53cm (21") - ST
48cm (19") - PE


 **Seat Height:**
48cm (19") - ST
44cm (17½") - PE

 **Seat Width:**
51cm (20") - ST
47cm (18½") - PE

Colour Options:



Berry



Biscuit



Charcoal



Mushroom

Leaflet Code: LL163 If you would like to view this product please contact your local stockist, details below.



Two Button Handset