



Invacare[®] | **Propad[™] Profile Cushion**

Profiled cut-foam cushion designed for comfort with
castellations and support



The *Invacare Propad Profile cushion* provides individuals with an option of having a lower 3" high cushion, rather than a standard height (4") original cushion. Some individuals prefer a lower profile cushion, finding it more compatible with their daily living activities.

Castellated Foam Cushion

Each castellation delivers enhanced support and pressure area care.

Reduces shear and friction forces

The castellations in the *Invacare Propad Profile cushion* reduce shear and friction forces on tissue.



Features and Options



Cover helps to prevent a build up of shear and friction forces

Cover is two-way stretch, water-resistant and vapour permeable. Cover features an anti-slip base.



Available in a range of sizes

The *Invacare Propad Profile Cushion* is available in 96 sizes, ranging from 10-24" widths and 10-22" depths.



Full length flap conceals the zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk

	Size range	Max. user weight	Max. cleaning temperature	Warranty	Manufactured to European Standards:
Propad Profile Cushion	 A range of 102 sizes	 95 kg (15 stone)	 80°	 1 year	Interior Foam: BS EN 597-1 Cig, BS EN597-2 Match, BS 5852 (Crib 5) Cover: BS 7175 (Crib 5)

Colours

Fuscia

Black Smooth

