

Permobil

# Progeo Joker Energy

Made-to-measure performance



Aluminum



Rigid



Minimized weight



Made-to-measure

permobil

# A made-to-measure agile ride

Progeo Joker Energy is a made-to-measure rigid wheelchair. The one-piece aluminum frame with optimized low weight ensures high performance while harnesses energy for propulsion. The oval tubing provides a grip friendly surface for easy transfers. The Joker Energy is designed for the agile users.



## Progeo Joker Energy

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Made-to-measure
Footplate distance	Made-to-measure (standard) Adjustable (optional)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	Made-to-measure
Total length	85 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Backrest angle	From 74° to 94° to the seat (adjustable)
Camber	0° - 2° - 4° - 8° - 12°
Weight without rear wheels	Starting from 4,7 kg approx.
Wheelchair weight	Starting from 7,3 kg approx.
Maximum load	125 kg

OPAQUE  
BLACK

PEARL  
WHITE

OPAQUE  
ORANGE

METALLIC  
MATTE BLUE

BRIGHT  
SILVER

Personalise your ride\*

Permobil reserves the right to change specifications and design without prior notice. \*More color options available.

Actual product colors may vary.

permobil

permobil.com | customer\_support\_uk@permobil.com | T. +44 (0)1484 722 888



PROGEOJOKERENERGYDATASHEETUK210909