

Permobil

# Progeo Ego Custom

Designed for autonomy and independence



Aluminium



Foldable



Minimized  
weight



Made-to-  
measure

permobil

# Light and foldable with impressive drivability

Progeo Ego Custom is an extremely light, foldable, active wheelchair with the rigidity of a fixed frame offering impressive drivability. The innovative easy-to-use horizontal folding system combined with the foldable backrest makes it truly compact and facilitates storage, transfers and traveling. The angle adjustable backrest supports positioning and trunk stability.



## Progeo Ego Custom

Seat width	33 - 36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 cm
Seat height	Front from 40 to 58 cm Rear from 38 to 45 cm
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	85° - 90° - 100°
Total length	88 cm approx. (standard configuration)
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	From 78° to 94° to the seat (adjustable)
Camber	0° - 2°
Weight without rear wheels	Starting from 6,4 kg approx.
Wheelchair weight	Starting from 9,0 kg approx.
Maximum load	125 kg

OPAQUE  
BLACK

PEARL  
WHITE

OPAQUE  
ORANGE

METALLIC  
MATTE BLUE

BRIGHT  
SILVER

Personalise your ride\*

Permobil reserves the right to change specifications and design without prior notice. \*More color options available

Actual product colors may vary.

**permobil**

permobil.com | customer\_support\_uk@permobil.com | T. +44 (0)1484 722 888



PROGEOEGOCUSTOMDATASHEETUK210909