

Your new route to independence



Rascal Vista & Vista DX Profile:

Introducing the exclusive Rascal Vista and Vista DX mid-size pavement scooters. Benefiting from the quick release feature it dismantles into a few easy-to-lift sections for transporting in your car. The compact & slim-line style offers plenty of legroom yet it is still extremely robust with generous turning manoeuvrability, making it ideal for everyday use.

Each available in two eye catching colours with appealing features including easy access tiller charging & delta handle bar plus the all-round suspension offers added comfort. When taking apart, the split battery pack is also a real advantage. It comes with front & rear LED lights maximising visibility, giving you both confidence & support wherever you go. Available with solid puncture proof or pneumatic tyres - the choice is yours!



Enjoy day trips without the worry!

The Vista DX's increased 50Ah battery size offers an impressive range of up to 34.5 miles & the luxury Captain's seat makes it an excellent solution for active people who require a multi-purpose scooter & want to stay out for longer with friends & family.

Vista DX shown in Sunset Orange



Main Features Vista:

- All-round suspension
- Easy take apart action car transportable
- Dual compound solid puncture proof tyres
- Comfort seat with folding backrest, adjustable arms & integrated swivel function
- Delta tiller with useful charging port
- Convenient front basket
- Class leading weight capacity
- Off-Board Battery Charging capability

Main Features Vista DX:

- All of the above features PLUS
- Captain's seat with folding backrest, headrest, adjustable arms with swivel & slide functions
- Range of up to 34.5 miles with 50Ah Battery capacity



Travel in comfort & style with the all-round suspension and slide & swivel seating



OPTIONAL

Options:

- Rascal Armrest bag
- Rear Saddle bag
- Metal Cane holder kit



Also available with Pneumatic



Fully dismantles into 6 easy to handle parts for storage & transportation





Tyre options & colour choices available on both models

Vista & Vista DX

Specification:	Vista	Vista DX
Type Class	B/2	\checkmark
Overall length	118cm / 46.5"	\checkmark
Overall width	58.5cm / 23"	✓
Floor to seat height	57cm - 62cm	✓
Seat width	44.5cm / 17.5"	45.7cm / 18"
Seat depth	43.1cm / 17"	43.1cm / 17"
Maximum carrying capacity	150kgs / 23.6 stone	140kgs/22 stone
Battery voltage and capacity	12V / 35Ah x 2	12V / 50Ah x 2
Maximum speed	6.4km/h / 4mph	\checkmark
Ground clearance	7.5cm / 3"	✓
Maximum kerb height	6cm / 2.3"	✓
Range*	up to 38km / 24 miles	up to 55.5km / 34.5 miles
The maximum safe slope	6°	✓
Turn-around width	1.7m / 5'8"	✓
Colour	Oxford Blue / Sunset Orange	Oxford Blue / Sunset Orange



*Range on full charge and flat ground

(front Chassis)

Single Heaviest component 25kgs / 55.12 lbs



Local Rascal Mobility Retailer:







