

The Scout has many new features which are only normally found on larger scooters. Its the first Drive mini scooter to include Delta Bars as standard, it also features the next generation Drive splitting mechanism.



The Scout is ideal for use indoors and outdoors where the terrain is forgiving (eg. smooth paths, shopping centres, etc). The Scout has been developed to ensure a focus on Comfort, Style, Performance and Reliability.




Superb Comfort, Style & Performance...

- » Top speed of 4mph*
- Maximum range of 10 miles (MS003 / 008), 13 miles (MS009) on full battery charge*
- » Comfortable lightweight seat which swivels 360° for easy transfers on and off the scooter
- » Removable battery pack with easy to carry handle and locking mechanism
- On board and in house charging as standard
- » Automatic electromagnetic brake system
- » Height adjustable seat to suit user comfort requirements
- » Delta bars as standard
- » Easy to change body panels
- » Easy tiller adjustment to suit user requirements
- » Removable seat, fold down tiller and splits easily for storage and transportation

- Stylish silver wheel hubs with black puncture proof non marking tyres
- » Width adjustable flip up armrests for added flexibility
- Front removable storage basket as standard
- » Rear anti tip wheels as standard
- Anti roll back system safety device prevents scooter from rolling backwards on hills
- » Freewheel facility allows movement of scooter without switching the motor on
- » Easy to remove battery pack fitted with 2 x 12 Ah or 2 x 20 Ah airline friendly batteries**
- » Weight capacity 130kg (20st)
- * Subject to conditions of use
- ** Batteries included as standard



Splits easily for storage and transportation.



DESCRIPTION: Scout



Length: 108cm (42.5")



Height: 91cm (36")



Width: 48cm (19")



Max Capacity: 130kg (20st)



Weight: 43kg (95lbs)



Max Gradient:




Delta Bars as standard.

360° swivelling seat.



Splits easily in seconds.



Removable battery.

Colour Options





ctric Blue Spi

Leaflet Code: LL170

If you would like to view this product please contact your local stockist, details below.