

Permobil

Progeo Tekna Advance Swing Away

Designed for ease of use



Aluminium



Foldable



Adjustable

permobil

Balanced movement and ease of use

Tekna Advance is a lightweight, easy to use, folding wheelchair designed for the active user. Open design, shortened folding frame and the swing-in/out footrests enhance easy transfers and transportation. A performance oriented variety of options to support seating and positioning needs is available. The user will experience a balanced and energy efficient ride.



Progeo Tekna Advance Swing Away

Seat width	36 - 39 - 42 - 45 - 48 cm
Seat depth	35 - 37,5 - 40 - 42,5 - 45 - 47,5 - 50 cm
Seat height	Front from 40 to 58 cm (adjustable) Rear from 38 to 46 cm (adjustable)
Footplate distance	From 35 to 54 cm (adjustable)
Backrest height	From 24,5 to 47 cm (adjustable)
Front frame angle	95° - 110° - 120°
Total length	85 cm approx.
Total width	Seat width + 17 cm (with 0° camber)
Encumbrance	28 cm approx. (with 0° camber)
Backrest angle	84° - 90° - 96° (adjustable)
Camber	0° - 2° - 4°
Weight without rear wheels	Starting from 8,6 kg approx.
Wheelchair weight	Starting from 11,2 kg approx.
Maximum load	125 kg

Personalise your ride*

PEARL
WHITE

OPAQUE
BLACK

BRIGHT
SILVER

Permobil reserves the right to change specifications and design without prior notice. *More color options available.

Actual product colors may vary.

permobil

permobil.com | customer_support_uk@permobil.com | T. +44 (0)1484 722 888



PROGEOTEKNAADVANCESWINGAWAYDATASHEETUK210909