

R-SERIES "HP 16"

Public price starting from £4950.00

MDR: 254494

Basic UDI-DI: **805378246HPSO69F3**

Order

Date

Quote

Dealer

Order reference

MODEL FEATURES

- Triride Software
- Customisable electronics
- Display with 5 levels of power and speed
- Electronic brake
- Reverse
- Cruise Control
- Assisted load
- Weight: 12.8 kg
- Range: up to 50 km
- Torque: 50Nm
- Lithium-ion battery 48V 10.5A
- Tyre City or Cross
- ECODRIVE-Energy saving mode
- Two driving modes: Normal/Sport
- Charging of the battery during downhill and braking

Top of the line performance and a larger diameter wheel for great driveability over longer distances, a real "mileage grinder" capable of satisfying even the most demanding users.

Great for tackling the most demanding situations, such as steeper inclines and uneven routes.



TO CHALLENGE DISTANT HORIZONS

Notes

WHEELCHAIR FEATURES

Brand
 Model
 Front tube's diameter
 Seat width
 Seat depth (E)
 Toe to front tube distance (A)
 Ground/knee height (D)
 Length of stretched arms (C)
 Front tubes angle (degree)
 Frame length (B)
 Front tubes width(F)*
 User's weight

Footrest

Fixed ☐
 Removable ☐

Internal ☐
 External ☐

USER'S PERSONAL DATA

Name
 Tel.:
 E-mail:
 Age Pathology:

DELIVERY ADDRESS


Delivery to:

 Postcode
 Address

 Town
 County

* The measurements for the width of the front tubes must be taken from the outside of the tubes at a height between 30 and 35cm

All fields are obligatory



N.B.: For wheelchairs with non-standard measurements or for users with special needs and dimensions the frame and driving controls can be customized with a quote (eg: dwarfism, bariatrics, etc)

Assembly, installation and instruction for use carried out by a Triride specialist technician

(on request) ☐

MTW: TREKKING POWER

Adaptation for MTW Trekking Power Kit	purchase (complete and attach dedicated form)	no cost
T0013	Adaptation for MTW Trekking Power (future purchase)	£185.00

COMMANDS FOR TETRAPLEGIC USERS

Configuration for Tetraplegic Commands N. _____
(For the desired configuration fill and attach the dedicated form)

Personalised configuration (attach picture or drawing)

"For more information, see the attached page with various configurations"

TRIRIDE CONFIGURATION

Battery

BS-R-368-S	Battery 48V 6.2 USB (297 Wh)	£0.00
BS-R-395-S	Battery 48V 10.5 USB (504 Wh)	default

Tyre

BS-R-334	Tyre City long life 16" x 2.125	default
BS-R-355	Tyre City long life 16" x 2.5	£0.00
BS-R-298	Tyre Cross long life 16" x 2.125	£0.00

Electronics

BS-O-096	PerformancePack : IBS+ICC+IBS button	£171.00
BS-O-101	IBS - Intelligent Braking System	£68.00
BS-O-097	ICC - Intelligent Cruise Control	£68.00
BS-O-086	ECODRIVE - Increased battery range	£68.00
BS-R-366	Electronic brake button / IBS	£44.00
BS-O-076	Analogue brake with thumb control	£95.00
BS-O-077	Automatic brake on accelerator release	£41.00
BS-O-079	Emergency brake with remote control	£95.00

Display LCD

BS-C-077	Display LCD Special with USB	default
BS-C-078	Display LCD Colour with USB	£53.00
BS-C-079	Display Colour with USB with customisable Triride programming	£126.00

customisation: _____ (21 character)

Throttle

BS-R-299	Thumb throttle ambidextrous	 LHS  RHS	£0.00	
L: BS-R-349 R: BS-R-348	Palm throttle with long lever	 LHS 	£20.00	
L: BS-R-333 R: BS-R-335	Half twist throttle	 LHS 	£38.00	
L: T0062 R: BS-R-324	Full twist throttle	 LHS 	£38.00	

Brakes

L: BS-R-319 R: BS-R-318	Single mechanical brake	 LHS 	£0.00	
L: BS-O-203 R: BS-O-273	Single mechanical parking brake (not available with hydraulic brake)	 LHS 	£14.00	
L: BS-O-274 R: BS-O-275	Single hydraulic brake	 LHS 	£122.00	
BS-O-218	Mechanical brake - double caliper		£149.00	
BS-O-194	Hydraulic brake - double caliper		£224.00	

N.B: The electronic brake (including IBS when required) can be actioned also by the brake levers of the device

Handlebar

BS-O-246	Straight (indicate width between 46cm and 50cm)	 cm	£0.00	
BS-O-208	4 Bends riser		£0.00	
BS-O-107	4 Bends high riser		£0.00	

Handlebar stem

BS-O-222	Adjustable 110 mm		£0.00	
BS-O-220	Adjustable 130 mm		£0.00	
BS-O-185	Adjustable 160 mm		£81.00	

Custom frame

BS-O-189	Extended frame (to bring the handlebar closer)		£81.00	
----------	--	--	--------	---

Frame colour

BS-O-078	Titanium		default	
BS-O-078-P	Polished		£252.00	
BS-O-217-N	Black Ral 9005		£135.00	
BS-O-217-R	Red Ral 3020		£135.00	
BS-O-217-W	White Ral 9003		£135.00	
BS-O-217-G	Yellow Ral 1003		£135.00	
BS-O-217-A	Orange Ral 2004		£135.00	
BS-O-217-B	Blue Ral 5017		£135.00	
BS-O-217	RAL on request		£207.00	



Personalized frame wraps

BS-O-087-M	Camouflage	£504.00	☐
BS-O-087-C	Carbon Effect	£504.00	☐
BS-O-087-F	Graffiti wrap Black/White	£504.00	☐
	Personalisation: Two-tone, gradients, etc.	on request	☐

Mudguard

BS-O-245	Mudguard Black	default	☐
BS-O-167	Mudguard Blue	£0.00	☐
BS-O-169	Mudguard Orange	£0.00	☐
BS-O-166	Mudguard Green	£0.00	☐
BS-O-164	Mudguard Red	£0.00	☐
BS-O-215	Mudguard Carbon	£96.00	☐

Handlescrews

BS-O-301	Handlescrews - Black	£0.00	☐
BS-O-128	Handlescrews - Blue	£0.00	☐
BS-O-127	Handlescrews - Orange	£0.00	☐
BS-O-130	Handlescrews - Green	£0.00	☐
BS-O-129	Handlescrews - Red	£0.00	☐
BS-O-139	Lever handlescrews for Tetraplegic user	£5.00	☐

SPEED SETTING

We herewith advise, that using the Triride with max. speed settings higher than 4mph is only permitted on public roads in certain countries and otherwise may only be used on private property. In the absence of approval, in accordance with the road traffic laws, the Triride with speed-settings >4mph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Triride with the speed-settings greater than 4mph.

When placing an order for a max. speed-setting higher than 4mph please confirm with your signature below, that you have understood the above information, you will follow strictly the country specific legislation and you will intend the usage of your Triride on private property.

Up to 4mph (6 Kph) ☐

Up to 8mph (13 Kph) ☐

>8mph (13kph) ☐

Date

Customer Name

Customer signature

OPTIONAL

Attachment Kit for SECOND WHEELCHAIR - please provide details below (prices vary as per specifics)

Tilting frame ☐

Connecting arms ☐

Clamps ☐

Brand

Model

Seat width cm

Seat depth cm

Front tubes diameter mm

ACCESSORIES / OPTIONAL

BS-R-368	Battery 48V 6.2 USB (297 Wh)	£651.00	
BS-R-395	Battery 48V 10.5 USB (504 Wh)	£753.00	
BS-O-012	Neoprene 20/19 battery cover	£24.00	
BS-R-294	Front spare battery holder with switch	£224.00	
BS-O-083	Central arch stand	£149.00	
BS-O-278	Lateral stand with wheels	£113.00	
BS-O-229	Rearview mirror	£29.00	LHS RHS
BS-O-251	Carry Trailer for fixed frame	£279.00	Rear axle diameter mm
BS-O-019	Handlebar bag with B-Twin adaptor 10 L	£144.00	
BS-O-020	Handlebar bag with B-Twin adaptor 6 L	£144.00	
BS-O-011	Handlebar basket with fabric cover 24 L	£144.00	
BS-O-214	Leg cover	£19.00	
BS-O-292	Calf support with adaptation to tilting frame	£142.00	
BS-O-213	Triride bag	£153.00	
BS-O-114	Arms stopper rings, pair	£16.00	
BS-O-174	Laser customised engraving on frame (on inox and black frame only)	£68.00	
BS-O-156	Digital engraving	£68.00	
BS-O-157	Personalised sticker	£59.00	
BS-O-081	Integrated front light, activated from display	default	
BS-O-022	Rear light with handlebar switch (combined with the front light)	£95.00	
BS-O-232	Light - Front 1200 Lumen	£100.00	
BS-O-198	Light Front USB	£29.00	
BS-O-180	Light - Rear USB	£33.00	
BS-O-015	Rear axle extension for wheelchair (please complete dedicated form)	on request	

Prices are intended for public and before tax

For details and photos of accessories and spare parts, see the relevant "spare parts and accessories" price list or visit our website [web www.triride.uk.com](http://www.triride.uk.com)

N.B.: The Company reserves the exclusive right to bring changes, without notice, if it deems it useful for the improvement of the Product